

Making a Difference in the Journey of Life



Product Catalog

Updated September 2018



At Hollister Incorporated, our Mission is to **make life more rewarding and dignified** for people who use our products and services. This Mission is based on our unwavering commitment to our customers and the communities we serve. We work in the healthcare industry, developing innovative medical products, services and education materials. Our promise and the overriding commitment of every Hollister Associate is to **make a positive difference in the lives of others.**

General Information

HCPCS Codes Disclaimer

The reimbursement information provided herein is intended to provide general information concerning coding of Hollister products only. Hollister does not guarantee coverage or payment for any products. The ultimate responsibility for proper coding, satisfying reimbursement requirements, and obtaining reimbursement remains with the provider. Coding and coverage policies and guidelines are complex, can vary from one carrier or region to another, and are updated frequently. Providers should check with their local carriers or intermediaries often and should consult with counsel, a reimbursement specialist, and/or the DMECS website for any coding, coverage, reimbursement, or billing questions.

HCPCS codes shown in this Product Catalog are current as of January 2018

Latex Statement

The specifications for the Hollister Continence Care products listed below include natural rubber latex as a component; our suppliers add natural rubber latex as part of the production process.

If a Hollister product(s) found in this catalog is not listed below, the specification for that product(s) and packaging does not include natural rubber latex as a component, and our component suppliers do not add natural rubber as part of their production process.

Current Stock Numbers:

1051	9107	9207	B14R	B12RB
1067	9108	9208	B16R	B14RB
1084	9109	9209		B16RB
9106	9206			

General Information

How to Use The Hollister Product Catalog

Hollister has created this simplified product catalog to help our customers find what they need in a quick and efficient manner.

If you...

know the stock number (or stock series number)

Go to the **Index** at the back of the catalog to find the page where it appears.

If you...

know the product brand or series name

Go to the **Table of Contents** to find where it appears or reference the Product Name Index.

If you...

know the HCPCS code

Go to the **HCPCS Code Index** at the back of the catalog to find the page where it appears.

Prior to use, be sure to read the Instructions for Use for information regarding Intended Use, Contraindications, Warnings, Precautions, and Instructions.

Table of Contents

Ostomy Care Products

Ceramide Infused Products

New Image Two-Piece Pouching System

Premier One-Piece Pouching System

Pouchkins Pediatric Pouching System

Ostomy Accessories

Irrigation Products

CenterPointLock Two-Piece Pouching System

Karaya 5 One-Piece Pouching System

Continence Care Products

Intermittent Catheters

Male External Catheters

Leg Bag Systems

Leg Bag Accessories

Collection Systems

Bowel Care Products

Fecal Collectors

Tube Fastener Products

Tube Fasteners

Wound Care Products

FastCast OLS Total Contact Cast System

Foam Dressings

Contact Layers with **TRIACT** Technology

CalciCare and **Restore** Calcium Alginates

Restore Hydrogel

Restore Wound Cleanser and Skin Care

Restore Hydrocolloids

Education • Glossary • Resources

Stock Number Index

Product Name Index

HCPCS Code Index

Choosing the Right Skin Barrier

Choosing the Right Skin Barrier Makes All the Difference

- No single skin barrier works for everyone
- Matching individual needs is at the core of skin barrier selection
 - Skin type: healthy, normal, thin, fragile, sensitive, dry, flaky, oily, sweaty
 - Type of stomal output: semi-formed/formed stool, liquid stool, urine
 - Change frequency

Types of Skin Barriers and Their Performance

CeraPlus *extended wear* skin barrier

- Contains all skin barrier characteristics (adhesion, absorption and erosion resistance) and is infused with ceramide – a naturally occurring lipid in the skin which protects against dryness
- Its formulation is designed to help protect the skin's natural moisture barrier and help maintain good peristomal skin health

Flexlend *extended wear* skin barrier

- Durable with high resistance to erosion
- Swells when it comes in contact with liquid to help form a seal between the skin and the stoma
- High initial tack and suitable for those who perspire heavily
- Strong adhesion and flexibility

FormaFlex *extended wear* skin barrier

- Helps to provide a secure peristomal skin seal as it will maintain its shape after being stretched or shaped
- Provides an easy-to-customize opening due to barrier material
- Easy to shape to individual stomas without the use of scissors

Choosing the Right Skin Barrier

FlexWear *standard wear* skin barrier

- Good balance of durability and gentleness providing moderate adhesion, and can handle frequent or infrequent pouch changes
- Resistant to fluid erosion and tends to hold its shape in the presence of fluids
- Suitable for all skin types, depending on change frequency
- Designed to be worn more than one day

SoftFlex *standard wear* skin barrier

- Appropriate choice for:
 - Those who have thin, fragile, or sensitive skin (such as those who have undergone chemo/radiation therapy)
 - People who desire frequent removal
- Offers low erosion resistance
- Designed to limit skin damage as its gentle to the skin, which allows for frequent pouch changes
- Thin and flexible without compromising security

Karaya 5 *standard wear* skin barrier

- Well-suited for sensitive skin or in some cases where allergic reactions to synthetic barriers has occurred
- Mildly adhesive – needs a belt or tape to secure pouching system

Ostomy Care

Ostomy Care

Continence
Care

Bowel Care

Tube
Fasteners

Wound Care

Education

Hollister Secure Start Services

Personalized Patient Support

Ostomy Care

Continence Care

Bowel Care

Tube Fasteners

Wound Care

Education



Support Throughout the Continuum of Care

We are proud to offer dedicated support for each and every patient we serve along the continuum of care. Hollister Secure Start services provide a lifetime of personalized support for our members.

Our team is focused on helping patients with:

- Finding the right product fit
- Navigating insurance coverage
- Identifying supplier options
- Providing condition specific education
- And more!

To learn more, call us at **1.888.808.7456**.



The Hollister logo and Secure Start are trademarks of Hollister Incorporated.
©2017 Hollister Incorporated.

Medicare Allowables Quick Reference Guide

Reimbursement for Fistula Management

Medicare Utilization Chart

Ceramide Infused Products

New Image Two-Piece Pouching System

CeraPlus Flat Skin Barriers with Integrated Floating Flange

CeraPlus Convex Skin Barriers with Integrated Floating Flange

Premier One-Piece Pouching System

Drainable Pouches with **CeraPlus** Flat Barriers

Drainable Pouches with **CeraPlus** Soft Convex Barriers

Drainable Pouches with **CeraPlus** Convex Barriers

Urostomy Pouches with **CeraPlus** Soft Convex Barriers

Non Sterile Kits — Drainable

Non Sterile Kits — Urostomy

Ostomy Accessories

Adapt CeraRing Flat Barrier Rings

Adapt CeraRing Convex Barrier Rings

Adapt CeraRing Oval Convex Barrier Rings

New Image Two-Piece Pouching System

Flat Skin Barriers with Integrated Floating Flange

Convex Skin Barriers with Integrated Floating Flange

Drainable Pouches

High Output Drainable Pouches

Closed Pouches

Urostomy Pouches

Non Sterile Kits — Drainable

Non Sterile Kits — Urostomy

Sterile Kits — Drainable

Sterile Kits — Urostomy

Premier One-Piece Pouching System

Drainable Pouches with Flat Barriers

Drainable Pouches with Soft Convex Barriers

Drainable Pouches with Convex Barriers

High Output Drainable Pouches

Closed Pouches with Flat Barriers

Stoma Cap

Urostomy Pouches with Flat Barriers

Urostomy Pouches with Convex Barriers

Urostomy Pouches with **CeraPlus** Soft Convex Barriers

Non Sterile Kits — Drainable

Non Sterile Kits — Urostomy

Pouchkins Pediatric Pouching System

Premie Pouch

Newborn Pouch

One-Piece Pouches

Two-Piece Skin Barrier

Two-Piece Pouches

Pediatric Ostomy Belt

Ostomy Accessories

- Adapt** Skin Barrier Rings
- Adapt** Skin Barrier Strips
- Adapt** Barrier Extenders
- Adapt** Skin Barrier Paste
- Adapt** Stoma Powder
- Adapt** Odor Eliminator and Lubricant
- Adapt** Ostomy Belts
- Adapt** Skin Prepping Product
- Adapt** Adhesives and Removers
- Flextend** (extended wear) Skin Barrier
Odor Eliminators
Cleaner/Decrystallizer
Drainable Pouch Clamp
Urostomy Drain Tube Adapter
- Hollister** Universal Catheter Access Port
Loop Ostomy Bridge
Bedside Drainage Collection System with Anti-Reflux Valve
Skin Barrier Pastes and Powders
Skin Barriers
Replacement Filters

Irrigation Products

CenterPointLock Two-Piece Pouching System

- Flat Skin Barriers with Floating Flange
- Convex Skin Barriers with Floating Flange
- Drainable Pouches
- Closed Pouches
- Urostomy Pouches
- Stoma Cap
- Irrigator Sleeve

Karaya 5 One-Piece Pouching System

- Drainable Pouches
- Closed Pouches
- Urostomy Pouches

General Medicare Product Maximum Allowables of Selected Ostomy Products

Two-Piece Pouching Systems	Maximum Quantity per Month
Flat, Standard Wear Skin Barriers with Floating Flanges	20
Drainable Pouches	20
Closed Pouches	60
Urostomy Pouches	20
Irrigator Sleeves	4

One-Piece Pouching Systems	Maximum Quantity per Month
Drainable Pouches with Flat, Standard Wear Skin Barriers	20
Closed Pouches with Flat, Standard Wear Skin Barriers	60
Stoma Caps	31

Pediatric Pouching Systems	Maximum Quantity per Month
One-Piece Drainable and Urostomy Pouches	20
Two-Piece Skin Barriers	20
Two-Piece Drainable and Urostomy Pouches	20
Two-Piece Mini Closed Pouches	60
Pediatric Ostomy Belt	1

Accessories	Maximum Quantity per Month
Skin Barrier Rings	*
Skin Barrier Paste	4 oz
Stoma Powder	10 oz/6 months
Lubricating Deodorant	*
Ostomy Belt	1
Skin Protective Wipes	150/6 months

Note: Medicare allows a usual maximum quantity of supplies each month. The table above provides the maximum allowable amounts. The amount of supplies you need may vary from the maximum quantities. That will be determined by your physician and must be on your order and in your medical record. If you require more ostomy supplies than what is listed in the table, your physician will need to provide a medical justification for the additional supplies. The ostomy supplier will need to have that medical justification on file when billing the Medicare program for the additional supplies.

* No usual maximum amount established.

Reimbursement for Fistula Management

Reimbursement For Fistula Management by Medicare

Reimbursement for fistula management under the Medicare program is described in the Surgical Dressings Policy Article. Specifically the article states:

“Surgical dressings are covered when either of the following criteria are met:

1. They are required for the treatment of a wound caused by, or treated by, a surgical procedure; **or**
2. They are required after debridement of a wound.”

The article goes on to explain:

“Examples of situations in which dressings are noncovered under the Surgical Dressings benefit are:

- a. Drainage from a cutaneous fistula which has not been caused by or treated by a surgical procedure”

For information see The Fistula Management Care Tip on www.hollister.com

Ostomy Care

Continence Care

Bowel Care

Tube Fasteners

Wound Care

Education

Sources: Local Coverage Article: Surgical Dressings - Policy Article - Effective 7-24-2017 (A52491)

Medicare's Usual Maximum Quantity for Ostomy Supplies

The quantity of ostomy supplies needed by a patient is determined primarily by the type of ostomy, its location, its construction, and the condition of the skin surface surrounding the stoma. There will be variation according to individual patient need and their needs may vary over time. The table below lists the maximum number of items/units of service that are usually medically necessary. The actual quantity needed for a particular patient may be more or less than the amount listed depending on the factors that affect the frequency of barrier and pouch change. The medical necessity for use of a greater quantity of supplies than the amounts listed must be clearly documented in the patient's medical record, and must be available upon request. If adequate documentation is not provided when requested, the excess quantities will be denied as not medically necessary. The HCPCS and reimbursement information provided herein is intended to provide general information concerning coding of Hollister products only. Please refer to the **HCPCS Code Disclaimer** for more information.

Per Medicare LCD for Ostomy Supplies (L33828), effective 1-1-2017.

HCPCS	Long Description	Maximum Quantity
A4331	Extension drainage tubing, any type, any length, with connector/adaptor, for use with urinary leg bag or urostomy pouch, each	*
A4357	Bedside drainage bag, day or night, with or without anti-reflux device, with or without tube, each	2 ea /mo
A4361	Ostomy faceplate, each	3 ea/6 mo
A4362	Skin barrier, solid, 4 x 4 or equivalent, each	20 ea/mo
A4363	Ostomy clamp, any type, replacement only, each	*
A4364	Adhesive, liquid or equal, any type, per ounce	4 oz/mo
A4366	Ostomy vent, any type, each	*
A4367	Ostomy belt, each	1 ea/mo
A4368	Ostomy filter, any type, each	*
A4369	Ostomy skin barrier, liquid (spray, brush, etc.), per ounce	2 oz/mo
A4371	Ostomy skin barrier, powder, per ounce	10 oz/6 mo
A4372	Ostomy skin barrier, solid 4 x 4 or equivalent, standard wear with built-in convexity, each	*
A4373	Ostomy skin barrier, with flange (solid, flexible, or accordion), with built-in convexity, any size, each	*
A4375	Ostomy pouch, drainable, with faceplate attached, plastic, each	*
A4376	Ostomy pouch, drainable, with faceplate attached, rubber, each	*
A4377	Ostomy pouch drainable for use on faceplate, plastic, each	10 ea/mo
A4378	Ostomy pouch, drainable, for use on faceplate, rubber, each	*
A4379	Ostomy pouch, urinary, with faceplate attached, plastic, each	*
A4380	Ostomy pouch, urinary, with faceplate attached, rubber, each	*
A4381	Ostomy pouch, urinary, for use with faceplate, plastic, each	10 ea/mo
A4382	Ostomy pouch, urinary, for use on faceplate, heavy plastic, each	*
A4383	Ostomy pouch, urinary, for use on faceplate, rubber, each	*
A4384	Ostomy faceplate equivalent, silicone ring, each	*

* No established usual maximum, refer to Medicare Utilization Guidelines.

Medicare Utilization Chart

Ostomy Care

HCPCS	Long Description	Quantity
A4385	Ostomy skin barrier, solid 4 x 4 or equivalent, extended wear, without built-in convexity, each	*
A4387	Ostomy pouch, closed, with barrier attached, with built-in convexity (one-piece), each	*
A4388	Ostomy pouch, drainable, with extended wear barrier attached, (one-piece), each	*
A4389	Ostomy pouch, drainable, with barrier attached, with built-in convexity (one-piece), each	*
A4390	Ostomy pouch, drainable, with extended wear barrier attached, with built-in convexity (one-piece), each	*
A4391	Ostomy pouch, urinary, with extended wear barrier attached (one-piece), each	*
A4392	Ostomy pouch, urinary, with standard wear barrier attached, with built-in convexity (one-piece), each	*
A4393	Ostomy pouch, urinary, with extended wear barrier attached, with built-in convexity (one-piece), each	*
A4394	Ostomy deodorant, with or without lubricant, for use in ostomy pouch, liquid, per fluid ounce	*
A4395	Ostomy deodorant for use in ostomy pouch, solid, per tablet	*
A4396	Ostomy belt with peristomal hernia support	*
A4397	Irrigation supply, sleeve, each	4 ea/mo
A4398	Ostomy irrigation supply, bag, each	2 ea/6 mo
A4399	Ostomy irrigation supply, cone/catheter, with or without brush	2 ea/6 mo
A4402	Lubricant, per ounce	4 oz/mo
A4404	Ostomy ring, each	10 ea/mo
A4405	Ostomy skin barrier, non-pectin based, paste, per ounce	4 oz/mo
A4406	Ostomy skin barrier, pectin-based, paste, per ounce	4 oz/mo
A4407	Ostomy skin barrier, with flange (solid, flexible, or accordion), extended wear, with built-in convexity, 4 x 4 inches or smaller, each	*
A4408	Ostomy skin barrier, with flange (solid, flexible, or accordion), extended wear, with built-in convexity, larger than 4 x 4 inches, each	*
A4409	Ostomy skin barrier, with flange (solid, flexible, or accordion), extended wear, without built-in convexity, 4 x 4 inches or smaller, each	*
A4410	Ostomy skin barrier, with flange (solid, flexible, or accordion), extended wear, without built-in convexity, larger than 4 x 4 inches, each	*
A4411	Ostomy skin barrier, solid 4 x 4 or equivalent, extended wear, with built-in convexity, each	*
A4412	Ostomy pouch, drainable, high output, for use on a barrier with flange (two-piece system), without filter, each	*
A4413	Ostomy pouch, drainable, high output, for use on a barrier with flange (two-piece system), with filter, each	*
A4414	Ostomy skin barrier, with flange (solid, flexible, or accordion), without built-in convexity, 4 x 4 inches or smaller, each	20 ea/mo
A4415	Ostomy skin barrier, with flange (solid, flexible, or accordion), without built-in convexity, larger than 4 x 4 inches, each	20 ea/mo

* No established usual maximum, refer to Medicare Utilization Guidelines.

Continence Care

Bowel Care

Tube Fasteners

Wound Care

Education

Medicare Utilization Chart

HCPCS	Long Description	Quantity
A4416	Ostomy pouch, closed, with barrier attached, with filter (one-piece), each	60 ea/mo
A4417	Ostomy pouch, closed, with barrier attached, with built-in convexity, with filter (one-piece), each	60 ea/mo
A4418	Ostomy pouch, closed, without barrier attached, with filter (one-piece), each	60 ea/mo
A4419	Ostomy pouch, closed, for use on barrier with non-locking flange, with filter (two-piece), each	60 ea/mo
A4420	Ostomy pouch, closed, for use on barrier with locking flange (two-piece), each	60 ea/mo
A4421	Ostomy supply; miscellaneous	*
A4422	Ostomy absorbent material (sheet/pad/crystal packet) for use in ostomy pouch to thicken liquid stomal output, each	*
A4423	Ostomy pouch, closed, for use on barrier with locking flange, with filter (two-piece), each	60 ea/mo
A4424	Ostomy pouch, drainable, with barrier attached, with filter (one-piece), each	20 ea/mo
A4425	Ostomy pouch, drainable, for use on barrier with non-locking flange, with filter (two-piece system), each	20 ea/mo
A4426	Ostomy pouch, drainable, for use on barrier with locking flange (two-piece system), each	20 ea/mo
A4427	Ostomy pouch, drainable, for use on barrier with locking flange, with filter (two-piece system) each	20 ea/mo
A4428	Ostomy pouch, urinary, with extended wear barrier attached, with faucet-type tap with valve (one-piece), each	*
A4429	Ostomy pouch, urinary, with barrier attached, with built-in convexity, with faucet-type tap with valve (one-piece), each	20 ea/mo
A4430	Ostomy pouch, urinary, with extended wear barrier attached, with built-in convexity, with faucet-type tap with valve (one-piece), each	*
A4431	Ostomy pouch, urinary, with barrier attached, with faucet-type tap, with valve (one-piece), each	20 ea/mo
A4432	Ostomy pouch, urinary, for use on barrier with non-locking flange, with faucet-type tap with valve (two-piece), each	20 ea/mo
A4433	Ostomy pouch, urinary, for use on barrier with locking flange (two-piece), each	20 ea/mo
A4434	Ostomy pouch, urinary, for use on barrier with locking flange, with faucet-type tap with valve (two-piece), each	20 ea/mo
A4435	Ostomy pouch, drainable, high output, with extended wear barrier (one-piece system), with or without filter, each	*
A4450	Tape, non-waterproof, per 18 square inches	40 ea (18 sq in)/mo
A4452	Tape, waterproof, per 18 square inches	40 ea (18 sq in)/mo
A4455	Adhesive remover or solvent (for tape, cement, or other adhesive), per ounce	16 oz/6 mo
A4456	Adhesive remover, wipes, any type, each	*
A5051	Ostomy pouch, closed, with barrier attached (one-piece), each	60 ea/mo
A5052	Ostomy pouch, closed, without barrier attached, (one-piece), each	60 ea/mo
A5053	Ostomy pouch, closed, for use on faceplate, each	60 ea/mo

* No established usual maximum, refer to Medicare Utilization Guidelines.

Medicare Utilization Chart

Ostomy Care

HCPCS	Long Description	Quantity
A5054	Ostomy pouch, closed, for use on barrier with flange (two-piece), each	60 ea/mo
A5055	Stoma Cap	31 ea/mo
A5056	Ostomy pouch, drainable, with extended wear barrier attached, with filter, (one-piece), each	40 ea/mo
A5057	Ostomy pouch, drainable, with extended wear barrier attached, with built-in convexity, with filter, (one-piece), each	40 ea/mo
A5061	Ostomy pouch, drainable, with barrier attached (one-piece), each	20 ea/mo
A5062	Ostomy pouch, drainable, without barrier attached (one-piece), each	20 ea/mo
A5063	Ostomy pouch, drainable, for use on barrier with flange (two-piece system), each	20 ea/mo
A5071	Ostomy pouch, urinary, with barrier attached (one-piece), each	20 ea/mo
A5072	Ostomy pouch, urinary, without barrier attached (one-piece), each	20 ea/mo
A5073	Ostomy pouch, urinary, for use on barrier with flange (two-piece), each	20 ea/mo
A5081	Continent device, plug for continent stoma	31 ea/mo
A5082	Continent device, catheter for continent stoma	1 ea/mo
A5083	Continent device, stoma absorptive cover for continent stoma	150 ea/mo
A5093	Ostomy accessory, convex insert	10 ea/mo
A5102	Bedside drainage bottle with or without tubing, rigid or expandable, each	2 ea/6 mo
A5120	Skin barrier, wipes or swabs, each	150 ea/6 mo
A5121	Skin barrier, solid, 6 x 6 or equivalent, each	20 ea/mo
A5122	Skin barrier, solid, 8 x 8 or equivalent, each	20 ea/mo
A5126	Adhesive or non-adhesive, disk or foam pad	20 ea/mo
A5131	Appliance cleaner, incontinence and ostomy appliances, per 16 ounce	1 ea (16 oz btl)/mo
A6154	Wound pouch, each	3/wk
A6209	Foam dressing, wound cover, sterile, pad size, 16 sq in or less, without adhesive boarder, each dressing	3/wk
A6216	Gauze, non-impregnated, non sterile, pad size 16 sq in or less, without adhesive border, each dressing	60 ea/mo
A9270	Non-covered item or service	*

* No established usual maximum, refer to Medicare Utilization Guidelines.



DID YOU KNOW...

How to Size **Hollister** Ostomy Pouches

A properly-fitted pouching system is essential to minimize stomal and skin irritation, and to obtain optimum pouch wear time. A stoma measuring guide allows correct measurement of the stoma. Measure the stoma before every barrier application. Select the barrier opening which fits where the skin and stoma meet. Verify that no skin is showing between the skin barrier and the stoma. Too large of an opening will leave skin exposed and could lead to skin irritation or erosion.

Continence Care

Bowel Care

Tube Fasteners

Wound Care

Education

Ostomy Care
Ceramide-Infused
Products

The science of skincare

is infused in our ostomy

care products



Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Ostomy Accessories

Continence
Care

Bowel Care

Tube
Fasteners

Wound Care

Education

Ceramide Products

[Education](#)

[Wound Care](#)

[Tube
Fasteners](#)

[Bowel Care](#)

[Continence
Care](#)

[Ostomy Care
Ostomy Accessories](#)

[Ostomy Care
Pouchkins Pediatric
Pouching System](#)

[Ostomy Care
Premier One-Piece
Pouching System](#)

[Ostomy Care
New Image Two-Piece
Pouching System](#)

[Ostomy Care
Ceramide Infused
Products](#)

Maintaining the Right Environment

Makes all the Difference



**Available on Premier One-Piece and
New Image Two-Piece Pouching Systems**

The CeraPlus Skin Barrier with Remois Technology*

At Hollister, we realize the importance of healthy peristomal skin. Our newest skin barrier is infused with ceramide, the skin's naturally occurring protection against dryness. The CeraPlus skin barrier is designed to help protect the skin's own moisture barrier, and is available on Premier one-piece and New Image two-piece pouching systems. This means you can enjoy the same pouch and skin barrier features from Hollister that you've come to rely on—and more!

The CeraPlus skin barrier from Hollister—because peristomal skin deserves better.

www.ceraplus.us



*Remois is a technology of Alcare Co., Ltd.

Ceramide Infused Products

Ceramide is a natural component of human skin, consisting of waxy lipid molecules. Ceramide links cells of the outer skin (epidermis) together to form a waterproof, protective barrier. Ceramides help prevent water loss that can lead to skin damage and dryness.

All stock numbers shaded in gray throughout this catalog denote ceramide infused products.

New Image Two-Piece Pouching System

Color Matching System

To help you select compatible New Image pouching system components, every box label is color-coded. First, select the proper skin barrier; then select a matching color-coded pouch ("red" skin barrier box label is compatible to "red" pouch box label).



CeraPlus Flat Skin Barriers with Integrated Floating Flange



1120x (CUT-TO-FIT)

Box of 5

- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ With Tape Border

HCPCS: A4409

Stock No	Flange Size	For Stomas Up To
11202	1-3/4" (44 mm)	1-1/4" (32 mm)
11203	2-1/4" (57 mm)	1-3/4" (44 mm)
11204	2-3/4" (70 mm)	2-1/4" (57 mm)



1510x (CUT-TO-FIT)

Box of 5

- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ Without Tape Border

HCPCS for 15102, 15103: A4409

HCPCS for 15104: A4410

Stock No	Flange Size	For Stomas Up To
15102	1-3/4" (44 mm)	1-1/4" (32 mm)
15103	2-1/4" (57 mm)	1-3/4" (44 mm)
15104	2-3/4" (70 mm)	2-1/4" (57 mm)

CeraPlus Convex Skin Barriers with Integrated Floating Flange



1140x (CUT-TO-FIT)

Box of 5

- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ With Tape Border



HCPCS: A4407

Stock No	Flange Size	For Stomas Up To
11402	1-3/4" (44 mm)	1" (25 mm)
11403	2-1/4" (57 mm)	1-1/2" (38 mm)
11404	2-3/4" (70 mm)	2" (51 mm)

CeraPlus Convex Skin Barriers with Integrated Floating Flange (continued)



1530x (CUT-TO-FIT)

- ▶ CeraPlus (extended wear) Skin Barrier
- ▶ Without Tape Border

Box of 5



HCPCS for 15302, 15303: A4407

HCPCS for 15304: A4408

Stock No	Flange Size	For Stomas Up To
15302	1-3/4" (44 mm)	1" (25 mm)
15303	2-1/4" (57 mm)	1-1/2" (38 mm)
15304	2-3/4" (70 mm)	2" (51 mm)



1150x (PRE-SIZED)

- ▶ CeraPlus (extended wear) Skin Barrier
- ▶ With Tape Border

Box of 5



HCPCS: A4407

Stock No	Flange Size	Barrier Opening
11503	1-3/4" (44 mm)	7/8" (22 mm)
11504	1-3/4" (44 mm)	1" (25 mm)
11505	2-1/4" (57 mm)	1-1/8" (29 mm)
11506	2-1/4" (57 mm)	1-1/4" (32 mm)

Premier One-Piece Pouching System

Drainable Pouches with CeraPlus Flat Barriers



8931 (CUT-TO-FIT)

- ▶ Drainable Pouch
- ▶ CeraPlus (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ Without Integrated AF300 Filter
- ▶ Lock 'n Roll Microseal Closure
- ▶ With ComfortWear Panels — Body Side Only

Box of 10

HCPCS: A4388

Ultra-Clear Stock No	For Stomas Up To	Pouch Length
8931	2-1/2" (64 mm)	12" (30 cm)



8901 (CUT-TO-FIT)

- ▶ Drainable Pouch
- ▶ CeraPlus (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Integrated AF300 Filter
- ▶ Lock 'n Roll Microseal Closure
- ▶ With ComfortWear Panels — Both Sides

Box of 10

HCPCS: A5056

Beige w/Viewing Option Stock No	For Stomas Up To	Pouch Length
8901	2-1/2" (64 mm)	12" (30 cm)

Ceramide Infused Products

Drainable Pouches with **CeraPlus** Flat Barriers (continued)



89xx (PRE-SIZED)

Box of 10

- ▶ Drainable Pouch
- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With *ComfortWear* Panels — Both Sides

HCPCS: A5056

Beige w/Viewing Option

Stock No	Barrier Opening	Pouch Length
8925	1" (25 mm)	12" (30 cm)
8930	1-3/16" (30 mm)	12" (30 cm)
8935	1-3/8" (35 mm)	12" (30 cm)
8940	1-9/16" (40 mm)	12" (30 cm)

Drainable Pouches with **CeraPlus** Soft Convex Barriers



895xx (CUT-TO-FIT)

Box of 5

- ▶ Drainable Pouch
- ▶ Soft Convex
- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels — Both Sides



HCPCS Code: A5057

Beige w/Viewing Option

Stock No	For Stomas Up To	Pouch Length
8954	1" (25 mm)	12" (30 cm)
8958	1-1/2" (38 mm)	12" (30 cm)
89511	2-1/8" (55 mm)	12" (30 cm)



896x (PRE-SIZED)

Box of 5

- ▶ Drainable Pouch
- ▶ Soft Convex
- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels — Both Sides



HCPCS: A5057

Beige w/Viewing Option

Stock No	Barrier Opening	Pouch Length
8960	3/4" (19 mm)	12" (30 cm)
8961	7/8" (22 mm)	12" (30 cm)
8962	1" (25 mm)	12" (30 cm)
8963	1-1/8" (29 mm)	12" (30 cm)

Drainable Pouches with CeraPlus Convex Barriers



891x (CUT-TO-FIT)

- ▶ Drainable Pouch
- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels — Both Sides

Box of 5



HCPCS Code: A5057

Beige w/Viewing Option

Stock No	For Stomas Up To	Pouch Length
8914	1" (25 mm)	12" (30 cm)
8918	1-1/2" (38 mm)	12" (30 cm)
89111	2" (51 mm)	12" (30 cm)

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Premier One-Piece
Pouching System



898x (CUT-TO-FIT)

- ▶ Drainable Pouch
- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ Without Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels — Body Side Only

Box of 5



HCPCS Code: A4390

Ultra-Clear

Stock No	For Stomas Up To	Pouch Length
8984	1" (25 mm)	12" (30 cm)
8988	1-1/2" (38 mm)	12" (30 cm)
89811	2" (51 mm)	12" (30 cm)

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Ostomy Accessories



899x (PRE-SIZED)

- ▶ Drainable Pouch
- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels — Both Sides

Box of 5



HCPCS: A5057

Beige w/Viewing Option

Stock No	Barrier Opening	Pouch Length
8990	3/4" (19 mm)	12" (30 cm)
8991	7/8" (22 mm)	12" (30 cm)
8992	1" (25 mm)	12" (30 cm)
8993	1-1/8" (29 mm)	12" (30 cm)

Continence
Care

Bowel Care

Tube
Fasteners

Wound Care

Education

Ceramide Infused Products

Urostomy Pouches with CeraPlus Soft Convex Barriers

Upgraded Design



8413xx (CUT-TO-FIT)

- ▶ Urostomy Pouch
- ▶ Anti-Reflux Valve
- ▶ CeraPlus (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Belt Tabs
- ▶ With ComfortWear Panels — Body Side Only

**Box of 5
includes 1 adapter**



NEW

Upgraded Design Features:

- ▶ Multi-Chamber Design
- ▶ Softer, Streamlined Tap
- ▶ Drain Valve Indicator

HCPCS: A4430

Ultra Clear Stock No	For Stomas Up To	Pouch Length
84134	1" (25 mm)	9" (23 cm)
84138	1-1/2" (38 mm)	9" (23 cm)
841311	2 1/8" (55 mm)	9" (23 cm)

Upgraded Design



841x (PRE-SIZED)

- ▶ Urostomy Pouch
- ▶ Anti-Reflux Valve
- ▶ CeraPlus (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Belt Tabs
- ▶ With ComfortWear Panels — Body Side Only

**Box of 5
includes 1 adapter**



NEW

Upgraded Design Features:

- ▶ Multi-Chamber Design
- ▶ Softer, Streamlined Tap
- ▶ Drain Valve Indicator

HCPCS: A4430

Beige w/Viewing Option Stock No	Barrier Opening	Pouch Length
8412	3/4" (19 mm)	9" (23 cm)
8413	1/8" (22 mm)	9" (23 cm)
8414	1" (25 mm)	9" (23 cm)
8415	1-1/8" (29 mm)	9" (23 cm)

Ceramide Infused Products

Non Sterile Kits — Drainable



1960x

Box of 5

Colostomy/Ileostomy Kit includes:

- ▶ (1) 1120x *CeraPlus* (extended wear) Skin Barrier with Tape Border
- ▶ (1) 1800x (*Lock 'n Roll* Microseal Closure) Ultra-Clear
- ▶ Single-Use Package

HCPCS: A5063 + A4414

Ultra-Clear Stock No	Flange Size	For Stomas Up To
19602	 1-3/4" (44 mm)	1-1/4" (32 mm)
19603	 2-1/4" (57 mm)	1-3/4" (44 mm)
19604	 2-3/4" (70 mm)	2-1/4" (57 mm)



1980x

Box of 5

Colostomy/Ileostomy Kit includes:

- ▶ (1) 1120x *CeraPlus* (extended wear) Skin Barrier with Tape Border
- ▶ (1) 1817x (Clamp Closure) Ultra-Clear
- ▶ Single-Use Package

HCPCS: A5063 + A4414

Ultra-Clear Stock No	Flange Size	For Stomas Up To
19802	 1-3/4" (44 mm)	1-1/4" (32 mm)
19803	 2-1/4" (57 mm)	1-3/4" (44 mm)
19804	 2-3/4" (70 mm)	2-1/4" (57 mm)

Non Sterile Kits — Urostomy



1990x

Box of 5

Urostomy Kit includes:

- ▶ (1) 1120x *CeraPlus* (extended wear) Skin Barrier with Tape Border
- ▶ (1) 1892x with *Comfortwear* panels (Ultra-Clear Pouch Film)
- ▶ (1) 7331 Urostomy Drain Tube Adapter
- ▶ Single-Use Package

HCPCS: A4432 + A4409

Ultra-Clear Stock No	Flange Size	For Stomas Up To
19902	 1-3/4" (44 mm)	1-1/4" (32 mm)
19903	 2-1/4" (57 mm)	1-3/4" (44 mm)
19904	 2-3/4" (70 mm)	2-1/4" (57 mm)

[Education](#)

[Wound Care](#)

[Tube
Fasteners](#)

[Bowel Care](#)

[Continence
Care](#)

[Ostomy Care
Accessories](#)

[Ostomy Care
Pouchkins Pediatric
Pouching System](#)

[Ostomy Care
Premier One-Piece
Pouching System](#)

[Ostomy Care
New Image Two-Piece
Pouching System](#)

[Ostomy Care
Ceramide Infused
Products](#)

Infused with the Science

of Skin Health



Adapt CeraRing barrier ring

Sealed with the confidence of ceramide

Adapt barrier ring technology was introduced more than 10 years ago. Its ability to be stretched and shaped to help minimize the risk of leakage and protect the peristomal skin makes it a favorite among WOC nurses.



Now, innovation has reached the next level, with the introduction of Adapt CeraRing.

Infused with ceramide—a natural component of the skin—Adapt CeraRing is designed to support the healthy peristomal skin your patients deserve. Ceramides help protect the skin's natural moisture barrier and help prevent water loss that may lead to damage and dryness.

Experience the confidence of ceramide—try Adapt CeraRing today.



Ceramide Infused Products

Ostomy Accessories

Adapt CeraRing Flat Barrier Rings



88xx

- ▶ Barrier Rings (Flat)
- ▶ *Flexlend M* (extended wear) with Ceramide

Box of 10

HCPCS: A4385

Stock No	Outer Diameter	Thickness
8805	2" (48 mm)	4.5 mm
8815	2" (48 mm)	2.3 mm SLIM


Adapt CeraRing Convex Barrier Rings



895xx

- ▶ Convex Barrier Rings
- ▶ *Flexlend* (extended wear) with Ceramide

Box of 10



HCPCS: A4411

Stock No	Inner Diameter
89520	13/16" (20 mm) — can be stretched to 1" (25 mm)
89530	1-3/16" (30 mm) — can be stretched to 1-3/8" (35 mm)
89540	1-9/16" (40 mm) — can be stretched to 1-3/4" (45 mm)

Adapt CeraRing Oval Convex Barrier Rings



896xx

- ▶ Oval Convex Barrier Rings
- ▶ *Flexlend* (extended wear) with Ceramide

Box of 10



HCPCS: A4411

Stock No	Inner Diameter
89601	7/8" x 1-1/2" (22 mm x 38 mm) — can be stretched to 1-1/8" x 1-3/4" (27 mm x 43 mm)
89602	1-3/16" x 1-7/8" (30 mm x 48 mm) — can be stretched to 1-3/8" x 2-1/8" (35 mm x 53 mm)
89603	1-1/2" x 2-3/16" (38 mm x 56 mm) — can be stretched to 1-3/4" x 2-3/8" (43 mm x 61 mm)

Ostomy Care
Ceramide-Infused
Products

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Ostomy Accessories

Continence
Care

Bowel Care

Tube
Fasteners

Wound Care

Education

Keeping peristomal skin

healthy is the inspiration

for all we do



New Image
Two-Piece Pouching System

Why choose a two-piece pouching system?

With a two-piece pouching system, the skin barrier and pouch are separate and connect together with a plastic ring, or flange. Here are a few reasons why you might want to use a two-piece pouching system:

- You may want to change your pouch more often than your skin barrier
- You may want to alternate between different pouch styles. For example, you might want to wear a smaller pouch during exercise or in the shower, and a larger pouch at night
- You may find it easier to center your skin barrier over your stoma using a two-piece system

New Image Two-Piece Pouching System

Color Matching System

To help you select compatible New Image pouching system components, every box label is color-coded. First, select the proper skin barrier; then select a matching color-coded pouch ("red" skin barrier box label is compatible to "red" pouch box label).



Flat Skin Barriers with Integrated Floating Flange



1120x (CUT-TO-FIT)

Box of 5

- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ With Tape Border

HCPCS: A4409

Stock No	Flange Size	For Stomas Up To
11202	1-3/4" (44 mm)	1-1/4" (32 mm)
11203	2-1/4" (57 mm)	1-3/4" (44 mm)
11204	2-3/4" (70 mm)	2-1/4" (57 mm)



1510x (CUT-TO-FIT)

Box of 5

- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ Without Tape Border

HCPCS for 15102, 15103: A4409

HCPCS for 15104: A4410

Stock No	Flange Size	For Stomas Up To
15102	1-3/4" (44 mm)	1-1/4" (32 mm)
15103	2-1/4" (57 mm)	1-3/4" (44 mm)
15104	2-3/4" (70 mm)	2-1/4" (57 mm)



1460x (CUT-TO-FIT)

Box of 5

- ▶ *Flextend* (extended wear) Skin Barrier
- ▶ With Tape Border

HCPCS for 14602, 14603 & 14604: A4409

HCPCS for 14606: A4410

Stock No	Flange Size	For Stomas Up To
14602	1-3/4" (44 mm)	1-1/4" (32 mm)
14603	2-1/4" (57 mm)	1-3/4" (44 mm)
14604	2-3/4" (70 mm)	2-1/4" (57 mm)
14606	4" (102 mm)	3-1/2" (89 mm)



1560x (CUT-TO-FIT)

Box of 5

- ▶ *Flextend* (extended wear) Skin Barrier
- ▶ Without Tape Border

HCPCS for 15602 & 15603: A4409

HCPCS for 15604: A4410

Stock No	Flange Size	For Stomas Up To
15602	1-3/4" (44 mm)	1-1/4" (32 mm)
15603	2-1/4" (57 mm)	1-3/4" (44 mm)
15604	2-3/4" (70 mm)	2-1/4" (57 mm)

New Image Two-Piece Pouching System

Flat Skin Barriers with Integrated Floating Flange (continued)



1420x (CUT-TO-FIT)

Box of 5

- ▶ *FlexWear* (standard wear) Skin Barrier
- ▶ With Tape Border

HCPCS for 14202, 14203 & 14204: A4414

HCPCS for 14206: A4415

Stock No	Flange Size	For Stomas Up To
14202	1-3/4" (44 mm)	1-1/4" (32 mm)
14203	2-1/4" (57 mm)	1-3/4" (44 mm)
14204	2-3/4" (70 mm)	2-1/4" (57 mm)
14206	4" (102 mm)	3-1/2" (89 mm)



1520x (CUT-TO-FIT)

Box of 5

- ▶ *FlexWear* (standard wear) Skin Barrier
- ▶ Without Tape Border

HCPCS for 15202 & 15203: A4414

HCPCS for 15204: A4415

Stock No	Flange Size	For Stomas Up To
15202	1-3/4" (44 mm)	1-1/4" (32 mm)
15203	2-1/4" (57 mm)	1-3/4" (44 mm)
15204	2-3/4" (70 mm)	2-1/4" (57 mm)



1470x (PRE-SIZED)

Box of 5

- ▶ *Flexend* (extended wear) Skin Barrier
- ▶ With Tape Border

HCPCS: A4409

Stock No	Flange Size	Barrier Opening
14701	1-3/4" (44 mm)	5/8" (16 mm)
14702	1-3/4" (44 mm)	3/4" (19 mm)
14703	1-3/4" (44 mm)	7/8" (22 mm)
14704	1-3/4" (44 mm)	1" (25 mm)
14705	1-3/4" (44 mm)	1-1/8" (29 mm)
14706	1-3/4" (44 mm)	1-1/4" (32 mm)
14707	2-1/4" (57 mm)	1-3/8" (35 mm)
14708	2-1/4" (57 mm)	1-1/2" (38 mm)
14709	2-1/4" (57 mm)	1-3/4" (44 mm)



1610x (PRE-SIZED)

Box of 5

- ▶ *Flexend* (extended wear) Skin Barrier
- ▶ Without Tape Border

HCPCS: A4409

Stock No	Flange Size	Barrier Opening
16101	1-3/4" (44 mm)	5/8" (16 mm)
16102	1-3/4" (44 mm)	3/4" (19 mm)
16103	1-3/4" (44 mm)	7/8" (22 mm)
16104	1-3/4" (44 mm)	1" (25 mm)
16105	1-3/4" (44 mm)	1-1/8" (29 mm)
16106	1-3/4" (44 mm)	1-1/4" (32 mm)
16107	2-1/4" (57 mm)	1-3/8" (35 mm)
16108	2-1/4" (57 mm)	1-1/2" (38 mm)



DID YOU KNOW...

- Each box of *New Image* pouches includes one disposal bag per pouch for discreet disposal
- Cutting guides are provided for all cut-to-fit barriers
- It is not recommended that you cut a pre-sized skin barrier
- It is not recommended to cut past the stop line on a cut-to-fit barrier

New Image Two-Piece Pouching System

Flat Skin Barriers with Integrated Floating Flange (continued)



1430x (PRE-SIZED)

Box of 5

- ▶ *FlexWear* (standard wear) Skin Barrier
- ▶ With Tape Border

HCPCS: A4414

Stock No	Flange Size	Barrier Opening
14301	1-3/4" (44 mm)	5/8" (16 mm)
14302	1-3/4" (44 mm)	3/4" (19 mm)
14303	1-3/4" (44 mm)	7/8" (22 mm)
14304	1-3/4" (44 mm)	1" (25 mm)
14305	1-3/4" (44 mm)	1-1/8" (29 mm)
14306	1-3/4" (44 mm)	1-1/4" (32 mm)
14307	2-1/4" (57 mm)	1-3/8" (35 mm)
14308	2-1/4" (57 mm)	1-1/2" (38 mm)
14309	2-1/4" (57 mm)	1-3/4" (44 mm)



1640x (PRE-SIZED)

Box of 5

- ▶ *FlexWear* (standard wear) Skin Barrier
- ▶ Without Tape Border

HCPCS: A4414

Stock No	Flange Size	Barrier Opening
16402	1-3/4" (44 mm)	3/4" (19 mm)
16403	1-3/4" (44 mm)	7/8" (22 mm)
16404	1-3/4" (44 mm)	1" (25 mm)
16405	1-3/4" (44 mm)	1-1/8" (29 mm)
16406	1-3/4" (44 mm)	1-1/4" (32 mm)
16407	2-1/4" (57 mm)	1-3/8" (35 mm)
16408	2-1/4" (57 mm)	1-1/2" (38 mm)



1410x (SHAPE-TO-FIT)

Box of 5

- ▶ *FormaFlex* (extended wear) Skin Barrier
- ▶ With Tape Border

HCPCS: A4409

Stock No	Flange Size	For Stomas Up To
14102	1-3/4" (44 mm)	1/2"-1-1/4" (13-32 mm)
14103	2-1/4" (57 mm)	5/8"-1-11/16" (16-43 mm)
14104	2-3/4" (70 mm)	5/8"-2-1/4" (16-57 mm)



DID YOU KNOW...

A True Original

The Hollister integrated floating flange—as featured in the New Image two-piece pouching system—offers flexibility and security for you and your patients. The first of its kind, this integrated floating flange helps relieve pressure on the tender postoperative or soft abdomen while offering an easy-to-use, low-profile solution patients can continue to use after they transition home. Included in all flat and convex New Image skin barriers.



New Image Two-Piece Pouching System

Convex Skin Barriers with Integrated Floating Flange



1140x (CUT-TO-FIT)

- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ With Tape Border

Box of 5



HCPCS: A4407

Stock No	Flange Size	For Stomas Up To
11402	1-3/4" (44 mm)	1" (25 mm)
11403	2-1/4" (57 mm)	1-1/2" (38 mm)
11404	2-3/4" (70 mm)	2" (51 mm)



1530x (CUT-TO-FIT)

- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ Without Tape Border

Box of 5



HCPCS for 15302, 15303: A4407
HCPCS for 15304: A4408

Stock No	Flange Size	For Stomas Up To
15302	1-3/4" (44 mm)	1" (25 mm)
15303	2-1/4" (57 mm)	1-1/2" (38 mm)
15304	2-3/4" (70 mm)	2" (51 mm)



1480x (CUT-TO-FIT)

- ▶ *Flexextend* (extended wear) Skin Barrier
- ▶ With Tape Border

Box of 5



HCPCS: A4407

Stock No	Flange Size	For Stomas Up To
14802	1-3/4" (44 mm)	1" (25 mm)
14803	2-1/4" (57 mm)	1-1/2" (38 mm)
14804	2-3/4" (70 mm)	2" (51 mm)



1580x (CUT-TO-FIT)

- ▶ *Flexextend* (extended wear) Skin Barrier
- ▶ Without Tape Border

Box of 5




HCPCS for 15802 & 15803: A4407
HCPCS for 15804: A4408

Stock No	Flange Size	For Stomas Up To
15802	1-3/4" (44 mm)	1" (25 mm)
15803	2-1/4" (57 mm)	1-1/2" (38 mm)
15804	2-3/4" (70 mm)	2" (51 mm)



DID YOU KNOW...

The **CeraPlus** skin barrier is infused with ceramide



- Ceramide is a natural component of the human skin
- Ceramide links the cells of the stratum corneum together to form a waterproof, protective barrier
- Ceramide helps prevent water loss that can lead to skin damage and dryness

Ostomy Care
Ceramide Infused Products

Ostomy Care
New Image Two-Piece Pouching System

Ostomy Care
Premier One-Piece Pouching System

Ostomy Care
Pouchkins Pediatric Pouching System

Ostomy Care
Ostomy Accessories

Continence Care

Bowel Care

Tube Fasteners

Wound Care

Education

New Image Two-Piece Pouching System

Convex Skin Barriers with Integrated Floating Flange (continued)

Upgraded Design Recommended replacement for 1440x series below



1340x (CUT-TO-FIT)

- ▶ *FlexWear* (standard wear) Skin Barrier
- ▶ With Tape Border

Box of 5



Upgraded Design Features:

- ▶ Flexible Release Liner
- ▶ Tapered Barrier Edges

HCPCS: A4373

Stock No	Flange Size	For Stomas Up To
13402	1-3/4" (44 mm)	1" (25 mm)
13403	2-1/4" (57 mm)	1-1/2" (38 mm)
13404	2-3/4" (70 mm)	2" (51 mm)

Original Design

1440x (CUT-TO-FIT)

- ▶ *FlexWear* (standard wear) Skin Barrier
- ▶ With Tape Border

Box of 5



HCPCS: A4373

Stock No	Flange Size	For Stomas Up To
14402	1-3/4" (44 mm)	1" (25 mm)
14403	2-1/4" (57 mm)	1-1/2" (38 mm)
14404	2-3/4" (70 mm)	2" (51 mm)



1540x (CUT-TO-FIT)

- ▶ *FlexWear* (standard wear) Skin Barrier
- ▶ Without Tape Border

Box of 5



HCPCS: A4373

Stock No	Flange Size	For Stomas Up To
15402	1-3/4" (44 mm)	1" (25 mm)
15403	2-1/4" (57 mm)	1-1/2" (38 mm)
15404	2-3/4" (70 mm)	2" (51 mm)



DID YOU KNOW...

Select **New Image** two-piece pouching systems feature two important product upgrades:

- **Lock 'n Roll** microseal closure with ribbed strips to facilitate ease of use, and transparent flap for increased security (as compared to previous design)
- Upgraded **skin barrier design** to help increase comfort and help improve ease of use:
 - Flexible release liner covers both the skin barrier and the tape border to help provide improved ease of use
 - Tapered barrier edges to help provide increased comfort
 - No change in skin barrier application
 - Same HCPCS code

New Image Two-Piece Pouching System

Convex Skin Barriers with Integrated Floating Flange (continued)



1150x (PRE-SIZED)

- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ With Tape Border

Box of 5



HCPCS: A4407

Stock No	Flange Size	Barrier Opening
11503	 1-3/4" (44 mm)	7/8" (22 mm)
11504	 1-3/4" (44 mm)	1" (25 mm)
11505	 2-1/4" (57 mm)	1-1/8" (29 mm)
11506	 2-1/4" (57 mm)	1-1/4" (32 mm)

Upgraded Design Recommended replacement for 149xx series below



139xx (PRE-SIZED)

- ▶ *Flexend* (extended wear) Skin Barrier
- ▶ With Tape Border

Box of 5



HCPCS: A4407

Upgraded Design Features:

- ▶ Flexible Release Liner
- ▶ Tapered Barrier Edge

Stock No	Flange Size	Barrier Opening
13901	 1-3/4" (44 mm)	5/8" (16 mm)
13902	 1-3/4" (44 mm)	3/4" (19 mm)
13903	 1-3/4" (44 mm)	7/8" (22 mm)
13904	 1-3/4" (44 mm)	1" (25 mm)
13905	 2-1/4" (57 mm)	1-1/8" (29 mm)
13906	 2-1/4" (57 mm)	1-1/4" (32 mm)
13907	 2-1/4" (57 mm)	1-3/8" (35 mm)
13908	 2-1/4" (57 mm)	1-1/2" (38 mm)
13909	 2-3/4" (70 mm)	1-5/8" (41 mm)
13910	 2-3/4" (70 mm)	1-3/4" (44 mm)
13911	 2-3/4" (70 mm)	2" (51 mm)

Original Design



149xx (PRE-SIZED)

- ▶ *Flexend* (extended wear) Skin Barrier
- ▶ With Tape Border

Box of 5



HCPCS: A4407

Stock No	Flange Size	Barrier Opening
14901	 1-3/4" (44 mm)	5/8" (16 mm)
14902	 1-3/4" (44 mm)	3/4" (19 mm)
14903	 1-3/4" (44 mm)	7/8" (22 mm)
14904	 1-3/4" (44 mm)	1" (25 mm)
14905	 2-1/4" (57 mm)	1-1/8" (29 mm)
14906	 2-1/4" (57 mm)	1-1/4" (32 mm)
14907	 2-1/4" (57 mm)	1-3/8" (35 mm)
14908	 2-1/4" (57 mm)	1-1/2" (38 mm)
14909	 2-3/4" (70 mm)	1-5/8" (41 mm)
14910	 2-3/4" (70 mm)	1-3/4" (44 mm)
14911	 2-3/4" (70 mm)	2" (51 mm)



DID YOU KNOW...

Convex Skin Barrier Options

The New Image two-piece pouching system delivers several skin barrier options to fit a variety of patient needs. Built-in convexity options provide security for recessed and flush stomas, and peristomal skin with soft muscle tone. Convex skin barriers are also available with and without tape.

New Image Two-Piece Pouching System

Convex Skin Barriers with Integrated Floating Flange (continued)



- 1590x** (PRE-SIZED)
 ▶ *Flexend* (extended wear) Skin Barrier
 ▶ Without Tape Border



HCPCS: A4407

Stock No	Flange Size	Barrier Opening
15903	1-3/4" (44 mm)	7/8" (22 mm)
15904	1-3/4" (44 mm)	1" (25 mm)
15905	2-1/4" (57 mm)	1-1/8" (29 mm)
15906	2-1/4" (57 mm)	1-1/4" (32 mm)

Upgraded Design Recommended replacement for 1450x series below



- 1350x** (PRE-SIZED)
 ▶ *FlexWear* (standard wear) Skin Barrier
 ▶ With Tape Border



Upgraded Design Features:


- ▶ Flexible Release Liner
- ▶ Tapered Barrier Edge

HCPCS: A4373

Stock No	Flange Size	Barrier Opening
13503	1-3/4" (44 mm)	7/8" (22 mm)
13504	1-3/4" (44 mm)	1" (25 mm)
13505	2-1/4" (57 mm)	1-1/8" (29 mm)
13506	2-1/4" (57 mm)	1-1/4" (32 mm)

Original Design

- 1450x** (PRE-SIZED)
 ▶ *FlexWear* (standard wear) Skin Barrier
 ▶ With Tape Border



HCPCS: A4373

Stock No	Flange Size	Barrier Opening
14503	1-3/4" (44 mm)	7/8" (22 mm)
14504	1-3/4" (44 mm)	1" (25 mm)
14505	2-1/4" (57 mm)	1-1/8" (29 mm)
14506	2-1/4" (57 mm)	1-1/4" (32 mm)



DID YOU KNOW...

Just released!



Get FIT with the new digital ostomy skin barrier **Fit Indicator Tool**

Developed with an international panel of experienced stoma care nurses, a new tool from Hollister Education has been designed to provide a simple patient assessment methodology to determine the need for a flat or convex barrier.

Ready to get started? Visit www.hollister.com/education or download the app by searching "Hollister Education" on iTunes or Google Play.

New Image Two-Piece Pouching System

Convex Skin Barriers with Integrated Floating Flange (continued)



1550x (PRE-SIZED)

- ▶ *FlexWear* (standard wear) Skin Barrier
- ▶ Without Tape Border

Box of 5



HCPCS: A4373

Stock No	Flange Size	Barrier Opening
15503	 1-3/4" (44 mm)	7/8" (22 mm)
15504	 1-3/4" (44 mm)	1" (25 mm)
15505	 2-1/4" (57 mm)	1-1/8" (29 mm)
15506	 2-1/4" (57 mm)	1-1/4" (32 mm)

Drainable Pouches



1818x/1819x

- ▶ Drainable Pouch
- ▶ With Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige — Both Sides

Box of 10

HCPCS: A4425

Transparent Stock No	Beige Stock No	Flange Size	Pouch Length
18192	18182	 1-3/4" (44 mm)	12" (30 cm)
18193	18183	 2-1/4" (57 mm)	12" (30 cm)
18194	18184	 2-3/4" (70 mm)	12" (30 cm)



1828x/1829x

- ▶ Drainable Mini-Pouch
- ▶ With Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige — Both Sides

Box of 20

HCPCS: A4425

Transparent Stock No	Beige Stock No	Flange Size	Pouch Length
18292	18282	 1-3/4" (44 mm)	7" (18 cm)
18293	18283	 2-1/4" (57 mm)	7" (18 cm)
18294	18284	 2-3/4" (70 mm)	7" (18 cm)



1811x/1813x/1800x

- ▶ Drainable Pouch
- ▶ Without Integrated Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels
 - Transparent and Ultra-Clear — Body Side Only
 - Beige — Both Sides

Box of 10

HCPCS: A5063

Ultra-Clear Stock No	Transparent Stock No	Beige Stock No	Flange Size	Pouch Length
18002	18132	18112	 1-3/4" (44 mm)	12" (30 cm)
18003	18133	18113	 2-1/4" (57 mm)	12" (30 cm)
18004	18134	18114	 2-3/4" (70 mm)	12" (30 cm)
18006	—	—	 4" (102 mm)	12" (30 cm)

Ostomy Care
Ceramide Infused Products

Ostomy Care
New Image Two-Piece Pouching System

Ostomy Care
Premier One-Piece Pouching System

Ostomy Care
Pouchkins Pediatric Pouching System

Ostomy Care
Ostomy Accessories

Continence Care

Bowel Care

Tube Fasteners

Wound Care

Education

New Image Two-Piece Pouching System

Drainable Pouches (continued)



1814x/1816x

- ▶ Drainable Pouch
- ▶ With Integrated AF300 Filter
- ▶ Clamp Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige — Both Sides

Box of 10
Includes 1 clamp

HCPCS: A4425

Transparent Stock No	Beige Stock No	Flange Size	Pouch Length
18162	18142	1-3/4" (44 mm)	12" (30 cm)
18163	18143	2-1/4" (57 mm)	12" (30 cm)
18164	18144	2-3/4" (70 mm)	12" (30 cm)



1810x/1812x/1817x

- ▶ Drainable Pouch
- ▶ Without Integrated Filter
- ▶ Clamp Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels
 - Transparent and Ultra-Clear — Body Side Only
 - Beige — Both Sides

Box of 10
Includes 1 clamp

HCPCS: A5063

Ultra-Clear Stock No	Transparent Stock No	Beige Stock No	Flange Size	Pouch Length
18172	18102	18122	1-3/4" (44 mm)	12" (30 cm)
18173	18103	18123	2-1/4" (57 mm)	12" (30 cm)
18174	18104	18124	2-3/4" (70 mm)	12" (30 cm)
18176	—	—	4" (102 mm)	12" (30 cm)



1820x

- ▶ Drainable Mini-Pouch
- ▶ Without Integrated Filter
- ▶ Clamp Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels — Both Sides

Box of 10
Includes 1 clamp

HCPCS: A5063

Beige Stock No	Flange Size	Pouch Length
18202	1-3/4" (44 mm)	9" (23 cm)
18203	2-1/4" (57 mm)	9" (23 cm)
18204	2-3/4" (70 mm)	9" (23 cm)



DID YOU KNOW...

- All *New Image* pouches are made up of multiple layers of odor-barrier film
- *New Image* barriers must be used with *New Image* pouches and can't be interchanged with other pouching systems

New Image Two-Piece Pouching System

High Output Drainable Pouches



1801x

Box of 10

- ▶ Drainable Pouch
- ▶ Without Integrated Filter
- ▶ Soft Drain Tap
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panel — Body Side Only
- ▶ No adapter needed to connect to bedside collector

HCPCS: A4412

Ultra-Clear Stock No	Flange Size	Pouch Length
18012	 1-3/4" (44 mm)	12" (30 cm)
18013	 2-1/4" (57 mm)	12" (30 cm)
18014	 2-3/4" (70 mm)	12" (30 cm)
18016	 4" (102 mm)	12" (30 cm)

Closed Pouches



1836x/1837x

Box of 60

- ▶ Closed Pouch
- ▶ With Integrated *AF300* Filter
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige — Both Sides

HCPCS: A4419

Transparent Stock No	Beige Stock No	Flange Size	Pouch Length
18362	18372	 1-3/4" (44 mm)	9" (23 cm)
18363	18373	 2-1/4" (57 mm)	9" (23 cm)
18364	18374	 2-3/4" (70 mm)	9" (23 cm)



1832x

Box of 30

- ▶ Closed Pouch
- ▶ With Integrated *AF300* Filter
- ▶ With Belt Tabs
- ▶ With *QuietWear* Pouch Material — Both Sides

HCPCS: A4419

Beige Stock No	Flange Size	Pouch Length
18322	 1-3/4" (44 mm)	9" (23 cm)
18323	 2-1/4" (57 mm)	9" (23 cm)
18324	 2-3/4" (70 mm)	9" (23 cm)



DID YOU KNOW...

New Image two-piece pouching system features **QuietWear** material for enhanced discretion

Hollister Ostomy offers many choices for patients—including the unique option of QuietWear pouch material on select closed pouches that reduces the noise that can be experienced with other pouches.

New Image Two-Piece Pouching System

Closed Pouches (continued)



1838x/1839x

Box of 60

- ▶ Closed Mini-Pouch
- ▶ With Integrated AF300 Filter
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige — Both Sides

HCPCS: A4419

Transparent Stock No	Beige Stock No	Flange Size	Pouch Length
18382	18392	1-3/4" (44 mm)	7" (18 cm)
18383	18393	2-1/4" (57 mm)	7" (18 cm)
18384	18394	2-3/4" (70 mm)	7" (18 cm)



1833x

Box of 30

1873x

Box of 60

- ▶ Closed Pouch
- ▶ Without Integrated Filter
- ▶ With Belt Tabs
- ▶ With *QuietWear* Pouch Material — Both Sides

HCPCS: A5054

Beige Stock No Box of 30	Beige Stock No Box of 60	Flange Size	Pouch Length
18332	18732	1-3/4" (44 mm)	9" (23 cm)
18333	18733	2-1/4" (57 mm)	9" (23 cm)
18334	18734	2-3/4" (70 mm)	9" (23 cm)



1835x

Box of 30

1875x

Box of 60

- ▶ Closed Mini-Pouch
- ▶ Without Integrated Filter
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels — Both Sides

HCPCS: A5054

Beige Stock No Box of 30	Beige Stock No Box of 60	Flange Size	Pouch Length
18352	18752	1-3/4" (44 mm)	7" (18 cm)
18353	18753	2-1/4" (57 mm)	7" (18 cm)
18354	18754	2-3/4" (70 mm)	7" (18 cm)

See Irrigation Products for Colostomy Irrigations

New Image Two-Piece Pouching System

Urostomy Pouches

Upgraded Design Recommended replacement for 1840x/1841x/1842xx series below



1890x/1891x/1892x

- ▶ Urostomy Pouch
- ▶ Anti-Reflux Valve
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels
 - Transparent and Ultra-Clear — Body Side Only
 - Beige — Both Sides

**Box of 10
Includes 2 adapters**

Upgraded Design Features:

- ▶ Multi-Chamber Design
- ▶ Softer, Streamlined Tap
- ▶ Drain Valve Indicator

HCPCS: A4432

Ultra-Clear Stock No	Transparent Stock No	Beige Stock No	Flange Size	Pouch Length
18922	18902	18912	1-3/4" (44 mm)	9" (23 cm)
18923	18903	18913	2-1/4" (57 mm)	9" (23 cm)
18924	18904	18914	2-3/4" (70 mm)	9" (23 cm)

Original Design



1840x/1841x/1842x

- ▶ Urostomy Pouch
- ▶ Anti-Reflux Valve
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels
 - Transparent and Ultra-Clear — Body Side Only
 - Beige — Both Sides

**Box of 10
Includes 1 adapter**

HCPCS: A4432

Transparent Stock No	Beige Stock No	Ultra-Clear Stock No	Flange Size	Pouch Length
18402	18412	18422	1-3/4" (44 mm)	9" (23 cm)
18403	18413	18423	2-1/4" (57 mm)	9" (23 cm)
18404	18414	18424	2-3/4" (70 mm)	9" (23 cm)



DID YOU KNOW...

- The upgraded 1890x urostomy pouch series comes with two adapters. The original 1840x urostomy pouch series only comes with one
- Please see **Hollister** bedside drainage collection system with anti-reflux valve or visit www.hollister.com for all **Hollister** leg bag systems

Ostomy Care
Ceramide Infused Products

Ostomy Care
New Image Two-Piece Pouching System

Ostomy Care
Premier One-Piece Pouching System

Ostomy Care
Pouchkins Pediatric Pouching System

Ostomy Care
Ostomy Accessories

Continence Care

Bowel Care

Tube Fasteners

Wound Care

Education

New Image Two-Piece Pouching System

Non Sterile Kits — Drainable



1960x Box of 5

Colostomy/Ileostomy Kit includes:

- ▶ (1) 1120x *CeraPlus* (extended wear) Skin Barrier with Tape Border
- ▶ (1) 1800x (*Lock 'n Roll* Microseal Closure) Ultra-Clear
- ▶ Single-Use Package

HCPCS: A5063 + A4414

Ultra-Clear Stock No	Flange Size	For Stomas Up To
19602	1-3/4" (44 mm)	1-1/4" (32 mm)
19603	2-1/4" (57 mm)	1-3/4" (44 mm)
19604	2-3/4" (70 mm)	2-1/4" (57 mm)



1980x Box of 5

Colostomy/Ileostomy Kit includes:

- ▶ (1) 1120x *CeraPlus* (extended wear) Skin Barrier with Tape Border
- ▶ (1) 1817x (Clamp Closure) Ultra-Clear
- ▶ Single-Use Package

HCPCS: A5063 + A4414

Ultra-Clear Stock No	Flange Size	For Stomas Up To
19802	1-3/4" (44 mm)	1-1/4" (32 mm)
19803	2-1/4" (57 mm)	1-3/4" (44 mm)
19804	2-3/4" (70 mm)	2-1/4" (57 mm)



1900x Box of 5

Colostomy/Ileostomy Kit includes:

- ▶ (1) 1420x *FlexWear* (standard wear) Skin Barrier with Tape Border
- ▶ (1) 1800x (*Lock 'n Roll* Microseal Closure) Ultra-Clear
- ▶ (1) 7767 Loop Ostomy Bridge (in 19006 only)
- ▶ Single-Use Package

HCPCS for 19002, 19003 & 19004: A5063 + A4414

HCPCS for 19006: A5063 + A4415

Ultra-Clear Stock No	Flange Size	For Stomas Up To
19002	1-3/4" (44 mm)	1-1/4" (32 mm)
19003	2-1/4" (57 mm)	1-3/4" (44 mm)
19004	2-3/4" (70 mm)	2-1/4" (57 mm)
19006	4" (102 mm)	3-1/2" (89 mm)



1910x Box of 5

Colostomy/Ileostomy Kit includes:

- ▶ (1) 1420x *FlexWear* (standard wear) Skin Barrier with Tape Border
- ▶ (1) 1817x (Clamp Closure) Ultra-Clear
- ▶ (1) 7767 Loop Ostomy Bridge (in 19106 only)
- ▶ Single-Use Package

HCPCS for 19102, 19103 & 19104: A5063 + A4414

HCPCS for 19106: A5063 + A4415

Ultra-Clear Stock No	Flange Size	For Stomas Up To
19102	1-3/4" (44 mm)	1-1/4" (32 mm)
19103	2-1/4" (57 mm)	1-3/4" (44 mm)
19104	2-3/4" (70 mm)	2-1/4" (57 mm)
19106	4" (102 mm)	3-1/2" (89 mm)

New Image Two-Piece Pouching System

Non Sterile Kits — Drainable (continued)



1930x

Box of 5

Colostomy/Ileostomy Kit includes:

- ▶ (1) 1410x *FormaFlex* (extended wear) Skin Barrier with Tape Border
- ▶ (1) 1800x (*Lock 'n Roll* Microseal Closure) Ultra-Clear
- ▶ Single-Use Package

HCPCS: A5063 + A4409

Ultra-Clear Stock No	Flange Size	For Stomas Up To
19303	2-1/4" (57 mm)	1-3/4" (44 mm)
19304	2-3/4" (70 mm)	2-1/4" (57 mm)

Non Sterile Kits — Urostomy



1990x

Box of 5

Urostomy Kit includes:

- ▶ (1) 1120x *CeraPlus* (extended wear) Skin Barrier with Tape Border
- ▶ (1) 1892x with *Comfortwear* panels (Ultra-Clear Pouch Film)
- ▶ (1) 7331 Urostomy Drain Tube Adapter
- ▶ Single-Use Package

HCPCS: A4432 + A4409

Ultra-Clear Stock No	Flange Size	For Stomas Up To
19902	1-3/4" (44 mm)	1-1/4" (32 mm)
19903	2-1/4" (57 mm)	1-3/4" (44 mm)
19904	2-3/4" (70 mm)	2-1/4" (57 mm)



1920x

Box of 5

Urostomy Kit includes:

- ▶ (1) 1460x *Flexend* (extended wear) Skin Barrier with Tape Border
- ▶ (1) 1842x with *Comfortwear* Panels (Ultra-Clear Pouch Film)
- ▶ (1) 7331 Urostomy Drain Tube Adapter
- ▶ Single-Use Package

HCPCS: A4432 + A4409

Ultra-Clear Stock No	Flange Size	For Stomas Up To
19202	1-3/4" (44 mm)	1-1/4" (32 mm)
19203	2-1/4" (57 mm)	1-3/4" (44 mm)
19204	2-3/4" (70 mm)	2-1/4" (57 mm)



1940x

Box of 5

Urostomy Kit includes:

- ▶ (1) 1410x *FormaFlex* (extended wear) Skin Barrier with Tape Border
- ▶ (1) 1842x with *Comfortwear* Panels (Ultra-Clear Pouch Film)
- ▶ (1) 7331 Urostomy Drain Tube Adapter
- ▶ Single-Use Package

HCPCS: A4432 + A4409

Ultra-Clear Stock No	Flange Size	For Stomas Up To
19402	1-3/4" (44 mm)	1-1/4" (32 mm)
19403	2-1/4" (57 mm)	1-3/4" (44 mm)

Ostomy Care
Ceramide Infused Products

Ostomy Care
New Image Two-Piece Pouching System

Ostomy Care
Premier One-Piece Pouching System

Ostomy Care
Pouchkins Pediatric Pouching System

Ostomy Care
Ostomy Accessories

Continence Care

Bowel Care

Tube Fasteners

Wound Care

Education

New Image Two-Piece Pouching System

Sterile Kits — Drainable



1905x

Box of 5

Colostomy/Ileostomy Kit includes:

- ▶ (1) 1420x *FlexWear* (standard wear) Skin Barrier with Tape Border
- ▶ (1) 1800x (*Lock 'n Roll* Microseal Closure) Ultra-Clear with *Comfortwear* Panels
- ▶ (1) 7767 Loop Ostomy Bridge (in 19056 only)
- ▶ Single-Use Package

HCPCS for 19052, 19053 & 19054: A5063 + A4414

HCPCS for 19056: A5063 + A4415

Ultra-Clear Stock No	Flange Size	For Stomas Up To
19052	1-3/4" (44 mm)	1-1/4" (32 mm)
19053	2-1/4" (57 mm)	1-3/4" (44 mm)
19054	2-3/4" (70 mm)	2-1/4" (57 mm)
19056	4" (102 mm)	3-1/2" (89 mm)



1915x

Box of 5
1 clamp per kit

Colostomy/Ileostomy Kit includes:

- ▶ (1) 1420x *FlexWear* (standard wear) Skin Barrier with Tape Border
- ▶ (1) 1817x (Clamp Closure) Ultra-Clear with *Comfortwear* Panels
- ▶ (1) 7767 Loop Ostomy Bridge (in 19156 only)
- ▶ Single-Use Package

HCPCS for 19152, 19153 & 19154: A5063 + A4414

HCPCS for 19156: A5063 + A4415

Ultra-Clear Stock No	Flange Size	For Stomas Up To
19152	1-3/4" (44 mm)	1-1/4" (32 mm)
19153	2-1/4" (57 mm)	1-3/4" (44 mm)
19154	2-3/4" (70 mm)	2-1/4" (57 mm)
19156	4" (102 mm)	3-1/2" (89 mm)



1935x

Box of 5

Colostomy/Ileostomy Kit includes:

- ▶ (1) 1520x *FlexWear* (standard wear) Skin Barrier **without Tape Border**
- ▶ (1) 1800x (*Lock 'n Roll* Microseal Closure) Ultra-Clear with *Comfortwear* Panels
- ▶ Single-Use Package

HCPCS for 19353: A5063 & A4414

HCPCS for 19354: A5063 & A4415

Ultra-Clear Stock No	Flange Size	For Stomas Up To
19353	2-1/4" (57 mm)	1-3/4" (44 mm)
19354	2-3/4" (70 mm)	2-1/4" (57 mm)

Ostomy Care
Ceramide Infused
Products

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Ostomy Accessories

Continence
Care

Bowel Care

Tube
Fasteners

Wound Care

Education

New Image Two-Piece Pouching System

Sterile Kits — Urostomy



1925x

Box of 5

Urostomy Kit includes:

- ▶ (1) 1460x *Flexlend* (extended wear) Skin Barrier with Tape Border
- ▶ (1) 1842x with *Comfortwear* Panels (Ultra-Clear Pouch Film)
- ▶ (1) 7331 Urostomy Drain Tube Adapter
- ▶ Single-Use Package

HCPCS: A4432 + A4409

Ultra-Clear Stock No	Flange Size	For Stomas Up To
19252	 1-3/4" (44 mm)	1-1/4" (32 mm)
19253	 2-1/4" (57 mm)	1-3/4" (44 mm)
19254	 2-3/4" (70 mm)	2-1/4" (57 mm)



1950x

Box of 5

Urostomy Kit includes:

- ▶ (1) 15603 *Flexlend* (extended wear) Skin Barrier **without Tape Border**
- ▶ (1) 18423 Ultra-Clear Pouch Film
- ▶ (1) 7331 Urostomy Drain Tube Adapter
- ▶ Single-Use Package

HCPCS: A4432 & A4409

Ultra-Clear Stock No	Flange Size	For Stomas Up To
19503	 2-1/4" (57 mm)	1-3/4" (44 mm)

Ostomy Care
Ceramide Infused
Products

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Ostomy Accessories

Continence
Care

Bowel Care

Tube
Fasteners

Wound Care

Education

A combination of the right

fit and formulation is the key

to positive outcomes



Premier

One-Piece Pouching System

Ostomy Care
Ceramide-Infused
Products

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Ostomy Accessories

Continence
Care

Bowel Care

Tube
Fasteners

Wound Care

Education

Why choose a one-piece pouching system?

With a one-piece pouching system, the skin barrier and pouch are connected and cannot be separated. Here are a few reasons why you might want to use a one-piece pouching system:

- You may want your pouching system to lay very flat against your body
- You need a pouching system that is flexible
- You may prefer a pouching system that does not require putting pieces together
- You do not wish to change your pouch more often than your skin barrier

Upgraded Product Design

Select Premier One-Piece Pouching Systems feature two important product updates:



- **Lock 'n Roll** microseal closure with **ribbed strips** to facilitate ease of use, and **transparent flap** for increased security
- Upgraded **skin barrier design** to help increase comfort and help improve ease of use:
 - **Flexible release liner** covers both the skin barrier and the tape border to help provide improved ease of use
 - **Tapered barrier edges** to help provide increased comfort
 - No change in skin barrier application
 - Same HCPCS code

Look for the **UPGRADED DESIGN** header above these stock numbers.

For more information, contact your Hollister Ostomy sales representative.



Premier One-Piece Pouching System

Drainable Pouches with Flat Barriers



8901 (CUT-TO-FIT)

Box of 10

- ▶ Drainable Pouch
- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With *ComfortWear* Panels — Both Sides

HCPCS: A5056

Beige w/Viewing Option

Stock No	For Stomas Up To	Pouch Length
8901	2-1/2" (64 mm)	12" (30 cm)



8931 (CUT-TO-FIT)

Box of 10

- ▶ Drainable Pouch
- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ Without Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With *ComfortWear* Panels — Body Side Only

HCPCS: A4388

Ultra-Clear

Stock No	For Stomas Up To	Pouch Length
8931	2-1/2" (64 mm)	12" (30 cm)



8331/8371 (CUT-TO-FIT)

Box of 10

- ▶ Drainable Pouch
- ▶ *Flextend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige with Viewing Option — Both Sides

HCPCS: A5056

Transparent Stock No	Beige w/Viewing Option Stock No	For Stomas Up To	Pouch Length
8331	8371	2-1/2" (64 mm)	12" (30 cm)



8531/8551 (CUT-TO-FIT)

Box of 10

- ▶ Drainable Pouch
- ▶ *Flextend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ Without Integrated Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With *ComfortWear* Panels
 - Beige — Both Sides
 - Ultra-Clear — Body Side Only

HCPCS: A4388

Ultra-Clear Stock No	Beige Stock No	For Stomas Up To	Pouch Length
8531	8551	2-1/2" (64 mm)	12" (30 cm)

Premier One-Piece Pouching System

Drainable Pouches with Flat Barriers (continued)



- 8631/8651 (CUT-TO-FIT)** **Box of 10**
Includes 1 clamp
- ▶ Drainable Pouch
 - ▶ *Flexend* (extended wear) Skin Barrier
 - ▶ With Tape Border
 - ▶ Without Integrated Filter
 - ▶ Clamp Closure
 - ▶ With *ComfortWear* Panels — Body Side Only

HCPCS: A4388

Ultra-Clear Stock No	Transparent Stock No	For Stomas Up To	Pouch Length
8651	8631	2-1/2" (64 mm)	12" (30 cm)



- 8641 (CUT-TO-FIT)** **Box of 10**
Includes 1 clamp
- ▶ Drainable Mini-Pouch
 - ▶ *Flexend* (extended wear) Skin Barrier
 - ▶ With Tape Border
 - ▶ Without Integrated Filter
 - ▶ Clamp Closure
 - ▶ With *ComfortWear* Panels — Body Side Only

HCPCS: A4388

Transparent Stock No	For Stomas Up To	Pouch Length
8641	2-1/2" (64 mm)	9" (23 cm)



- 8081/8181/8381 (CUT-TO-FIT)** **Box of 10**
- ▶ Drainable Pouch
 - ▶ *FlexWear* (standard wear) Skin Barrier
 - ▶ With Tape Border
 - ▶ With Integrated *AF300* Filter
 - ▶ *Lock 'n Roll* Microseal Closure
 - ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige and Beige with Viewing Option — Both Sides


HCPCS: A4424

Transparent Stock No	Beige Stock No	Beige w/Viewing Option Stock No	For Stomas Up To	Pouch Length
8081	8181	8381	2-1/2" (64 mm)	12" (30 cm)



DID YOU KNOW...

The **CeraPlus** skin barrier is infused with ceramide



- Ceramide is a natural component of the human skin
- Ceramide links the cells of the stratum corneum together to form a waterproof, protective barrier
- Ceramide helps prevent water loss that can lead to skin damage and dryness

Premier One-Piece Pouching System

Drainable Pouches with Flat Barriers (continued)



8031 (CUT-TO-FIT)

Box of 10

- ▶ Drainable Pouch
- ▶ *FlexWear* (standard wear) Skin Barrier
- ▶ With Tape Border
- ▶ Without Integrated Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With *ComfortWear* Panels — Body Side Only

HCPCS: A5061

Ultra-Clear Stock No	For Stomas Up To	Pouch Length
8031	2-1/2" (64 mm)	12" (30 cm)



8131 (CUT-TO-FIT)

Box of 10
Includes 1 clamp

- ▶ Drainable Pouch
- ▶ *FlexWear* (standard wear) Skin Barrier
- ▶ With Tape Border
- ▶ Without Integrated Filter
- ▶ Clamp Closure
- ▶ With *ComfortWear* Panels — Body Side Only

HCPCS: A5061

Ultra-Clear Stock No	For Stomas Up To	Pouch Length
8131	2-1/2" (64 mm)	12" (30 cm)



88300/88400 (CUT-TO-FIT)

Box of 10

- ▶ Drainable Pouch
- ▶ *SoftFlex* (standard wear) Skin Barrier
- ▶ **Without Tape Border**
- ▶ With Integrated AF300 Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige — Both Sides

HCPCS: A4424

Transparent Stock No	Beige Stock No	For Stomas Up To	Pouch Length
88400	88300	2-1/8" (55 mm)	12" (30 cm)



88800/88700 (CUT-TO-FIT)

Box of 20

- ▶ Drainable Mini-Pouch
- ▶ *SoftFlex* (standard wear) Skin Barrier
- ▶ **Without Tape Border**
- ▶ With Integrated AF300 Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige — Both Sides

HCPCS: A4424

Transparent Stock No	Beige Stock No	For Stomas Up To	Pouch Length
88800	88700	2-1/8" (55 mm)	7" (18 cm)

Premier One-Piece Pouching System

Drainable Pouches with Flat Barriers (continued)



- 88302/88402** (OVAL CUT-TO-FIT)
- ▶ Drainable Pouch
 - ▶ *SoftFlex* (standard wear) Skin Barrier
 - ▶ **Without Tape Border**
 - ▶ With Integrated *AF300* Filter
 - ▶ *Lock 'n Roll* Microseal Closure
 - ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige — Both Sides

Box of 10

OVAL

HCPCS: A4424

Transparent Stock No	Beige Stock No	For Stomas Up To	Pouch Length
88402	88302	3" x 2.5" (75 mm x 65 mm)	12" (30 cm)



- 88401** (CUT-TO-FIT)
- ▶ Drainable Pouch
 - ▶ *SoftFlex* (standard wear) Skin Barrier
 - ▶ **Without Tape Border**
 - ▶ Without Integrated Filter
 - ▶ *Lock 'n Roll* Microseal Closure
 - ▶ With *ComfortWear* Panels — Body Side Only

Box of 10

HCPCS: A5061

Ultra-Clear Stock No	For Stomas Up To	Pouch Length
88401	2-1/2" (64 mm)	12" (30 cm)



- 88501** (CUT-TO-FIT)
- ▶ Drainable Pouch
 - ▶ *SoftFlex* (standard wear) Skin Barrier
 - ▶ **Without Tape Border**
 - ▶ Without Integrated Filter
 - ▶ Clamp Closure
 - ▶ With *ComfortWear* Panels — Body Side Only

Box of 10
Includes 1 clamp

HCPCS: A5061

Ultra-Clear Stock No	For Stomas Up To	Pouch Length
88501	2-1/2" (64 mm)	12" (30 cm)



- 89xx** (PRE-SIZED)
- ▶ Drainable Pouch
 - ▶ *CeraPlus* (extended wear) Skin Barrier
 - ▶ With Tape Border
 - ▶ With Integrated *AF300* Filter
 - ▶ *Lock 'n Roll* Microseal Closure
 - ▶ With *ComfortWear* Panels — Both Sides

Box of 10

HCPCS: A5056

Beige w/Viewing Option Stock No	Barrier Opening	Pouch Length
8925	1" (25 mm)	12" (30 cm)
8930	1-3/16" (30 mm)	12" (30 cm)
8935	1-3/8" (35 mm)	12" (30 cm)
8940	1-9/16" (40 mm)	12" (30 cm)

Ostomy Care
Ceramide Infused
Products

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Ostomy Accessories

Continence
Care

Bowel Care

Tube
Fasteners

Wound Care

Education

Premier One-Piece Pouching System

Drainable Pouches with Flat Barriers (continued)



853x/855x (PRE-SIZED)

- ▶ Drainable Pouch
- ▶ *Flexextend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ Without Integrated Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige — Both Sides

Box of 10

HCPCS: A4388

Transparent Stock No	Beige Stock No	Barrier Opening	Pouch Length
8532	8552	1" (25 mm)	12" (30 cm)
8538	8558	1-1/4" (32 mm)	12" (30 cm)
8533	8553	1-1/2" (38 mm)	12" (30 cm)
8539	8559	1-3/4" (44 mm)	12" (30 cm)



863x (PRE-SIZED)

- ▶ Drainable Pouch
- ▶ *Flexextend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ Without Integrated Filter
- ▶ Clamp Closure
- ▶ With *ComfortWear* Panels — Body Side Only

**Box of 10
Includes 1 clamp**

HCPCS: A4388

Transparent Stock No	Barrier Opening	Pouch Length
8637	3/4" (19 mm)	12" (30 cm)
8632	1" (25 mm)	12" (30 cm)
8638	1-1/4" (32 mm)	12" (30 cm)
8633	1-1/2" (38 mm)	12" (30 cm)
8639	1-3/4" (44 mm)	12" (30 cm)
8634	2" (51 mm)	12" (30 cm)
8635	2-1/2" (64 mm)	12" (30 cm)



864x (PRE-SIZED)

- ▶ Drainable Mini-Pouch
- ▶ *Flexextend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ Without Integrated Filter
- ▶ Clamp Closure
- ▶ With *ComfortWear* Panels — Both Sides

**Box of 10
Includes 1 clamp**

HCPCS: A4388

Beige Stock No	Barrier Opening	Pouch Length
8647	3/4" (19 mm)	9" (23 cm)
8642	1" (25 mm)	9" (23 cm)
8648	1-1/4" (32 mm)	9" (23 cm)
8643	1-1/2" (38 mm)	9" (23 cm)
8649	1-3/4" (44 mm)	9" (23 cm)

Premier One-Piece Pouching System

Drainable Pouches with Flat Barriers (continued)



80xx/81xx (PRE-SIZED)

Box of 10

- ▶ Drainable Pouch
- ▶ *FlexWear* (standard wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige — Both Sides

HCPCS: A4424

Transparent Stock No	Beige Stock No	Barrier Opening	Pouch Length
8025	8125	1" (25 mm)	12" (30 cm)
8030	8130	1-3/16" (30 mm)	12" (30 cm)
8035	8135	1-3/8" (35 mm)	12" (30 cm)
8040	8140	1-9/16" (40 mm)	12" (30 cm)



883xx/884xx (PRE-SIZED)

Box of 10

- ▶ Drainable Pouch
- ▶ *SoftFlex* (standard wear) Skin Barrier
- ▶ **Without Tape Border**
- ▶ With Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige — Both Sides

HCPCS: A4424

Transparent Stock No	Beige Stock No	Barrier Opening	Pouch Length
88425	88325	1" (25 mm)	12" (30 cm)
88430	88330	1-3/16" (30 mm)	12" (30 cm)
88435	88335	1-3/8" (35 mm)	12" (30 cm)
88440	88340	1-9/16" (40 mm)	12" (30 cm)



887xx (PRE-SIZED)

Box of 20

- ▶ Drainable Mini-Pouch
- ▶ *SoftFlex* (standard wear) Skin Barrier
- ▶ **Without Tape Border**
- ▶ With Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With *ComfortWear* Panels — Both Sides

HCPCS: A4424

Beige Stock No	Barrier Opening	Pouch Length
88725	1" (25 mm)	7" (18 cm)
88730	1-3/16" (30 mm)	7" (18 cm)
88735	1-3/8" (35 mm)	7" (18 cm)
88740	1-9/16" (40 mm)	7" (18 cm)

Ostomy Care
Ceramide Infused
Products

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Ostomy Accessories

Continence
Care

Bowel Care

Tube
Fasteners

Wound Care

Education

Ostomy Care
Ceramide Intused Products

Ostomy Care
New Image Two-Piece Pouching System

Ostomy Care
Premier One-Piece Pouching System

Ostomy Care
Pouchkins Pediatric Pouching System

Ostomy Care
Ostomy Accessories

Continence Care

Bowel Care

Tube Fasteners

Wound Care

Education



A Flexible Fit

Made Secure

Premier Soft Convex Pouching System

Designed to stay secure, no matter what shape the body takes.

When it comes to finding the right solution for patients, fit comes first. Comfort and ease-of-use are also essential to help ensure positive outcomes—in the hospital and after patients leave your care.

Hollister Ostomy presents Premier soft convex—a pouching system that delivers both convex support and flexibility. The skin barrier is designed to conform to the body, without compromising fit. Belt tabs and a tape-bordered barrier help provide an extra sense of security. And the new viewing option offers a convenient way to help position the pouch and observe output.

- Look for the soft convex symbol 
- Available in CeraPlus and Flextend (extended wear) skin barrier formulations



Premier One-Piece Pouching System

Drainable Pouches with Soft Convex Barriers



895xx (CUT-TO-FIT)

- ▶ Drainable Pouch
- ▶ Soft Convex
- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels — Both Sides

Box of 5



HCPCS Code: A5057

Beige w/Viewing Option

Stock No	For Stomas Up To	Pouch Length
8954	1" (25 mm)	12" (30 cm)
8958	1-1/2" (38 mm)	12" (30 cm)
89511	2-1/8" (55 mm)	12" (30 cm)



857x/867x (CUT-TO-FIT)

- ▶ Drainable Pouch
- ▶ Soft Convex
- ▶ *Flextend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels
 - Ultra-Clear — Body Side Only
 - Beige with Viewing Option — Both Sides

Box of 5



HCPCS: A5057

Ultra-Clear Stock No	Beige w/Viewing Option Stock No	For Stomas Up To	Pouch Length
8574	8674	1" (13-25 mm)	12" (30 cm)
8578	8678	1 1/2" (13-38 mm)	12" (30 cm)
85711	86711	2 1/8" (13-55 mm)	12" (30 cm)



896x (PRE-SIZED)

- ▶ Drainable Pouch
- ▶ Soft Convex
- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels — Both Sides

Box of 5



HCPCS: A5057

Beige w/Viewing Option Stock No	Opening	Pouch Length
8960	3/4" (19 mm)	12" (30 cm)
8961	7/8" (22 mm)	12" (30 cm)
8962	1" (25 mm)	12" (30 cm)
8963	1-1/8" (29 mm)	12" (30 cm)

Ostomy Care
Ceramide Infused
Products

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Ostomy Accessories

Continence
Care

Bowel Care

Tube
Fasteners

Wound Care

Education

Premier One-Piece Pouching System

Drainable Pouches with Soft Convex Barriers (continued)



866x (PRE-SIZED)

- ▶ Drainable Pouch
- ▶ Soft Convex
- ▶ *FlexTend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels — Both Sides

Box of 5



HCPCS: A5057

Beige w/Viewing Option Stock No	Barrier Opening	Pouch Length
8661	7/8" (22 mm)	12" (30 cm)
8662	1" (25 mm)	12" (30 cm)
8663	1-1/8" (29 mm)	12" (30 cm)
8664	1-1/4" (32 mm)	12" (30 cm)
8665	1-1/2" (38 mm)	12" (30 cm)
8666	1-3/4" (44 mm)	12" (30 cm)

Drainable Pouches with Convex Barriers



891x (CUT-TO-FIT)

- ▶ Drainable Pouch
- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels — Both Sides

Box of 5



HCPCS Code: A5057

Beige w/Viewing Option Stock No	For Stomas Up To	Pouch Length
8914	1" (25 mm)	12" (30 cm)
8918	1-1/2" (38 mm)	12" (30 cm)
89111	2" (51 mm)	12" (30 cm)



898x (CUT-TO-FIT)

- ▶ Drainable Pouch
- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ Without Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels — Body Side Only

Box of 5



HCPCS Code: A4390

Ultra-Clear Stock No	For Stomas Up To	Pouch Length
8984	1" (25 mm)	12" (30 cm)
8988	1-1/2" (38 mm)	12" (30 cm)
89811	2" (51 mm)	12" (30 cm)

Ostomy Care
Ceramide Infused Products

Ostomy Care
New Image Two-Piece Pouching System

Ostomy Care
Premier One-Piece Pouching System

Ostomy Care
Pouchkins Pediatric Pouching System

Ostomy Care
Ostomy Accessories

Continence Care

Bowel Care

Tube Fasteners

Wound Care

Education

Premier One-Piece Pouching System

Drainable Pouches with Convex Barriers (continued)

Ostomy Care
 Ceramic Infused Products
 New Image Two-Piece Pouching System
 Ostomy Care
 Premier One-Piece Pouching System
 Ostomy Care
 Pouchkins Pediatric Pouching System
 Ostomy Care
 Ostomy Accessories
 Continence Care
 Bowel Care
 Tube Fasteners
 Wound Care
 Education

Upgraded Design Recommended replacement for 858xx series below



827xx/854xx (CUT-TO-FIT)

- ▶ Drainable Pouch
- ▶ *Flexend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Integrated *AF300* Filter
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige with Viewing Option — Both Sides

Box of 5



Upgraded Design Features:

- ▶ *Lock 'n Roll* Microseal Closure with Ribbed Strips and Transparent Flap
- ▶ Flexible Release Liner
- ▶ Tapered Barrier Edges

HCPCS Code: A5057

Transparent Stock No	Beige w/Viewing Option Stock No	For Stomas Up To	Pouch Length
8274	8544	1" (25 mm)	12" (30 cm)
8278	8548	1-1/2" (38 mm)	12" (30 cm)
82711	—	2" (51 mm)	12" (30 cm)

Original Design

858xx (CUT-TO-FIT)

- ▶ Drainable Pouch
- ▶ *Flexend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Integrated *AF300* Filter
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels — Body Side Only
- ▶ *Lock 'n Roll* Closure

Box of 5



HCPCS Code: A5057

Transparent Stock No	For Stomas Up To	Pouch Length
8584	1" (25 mm)	12" (30 cm)
8588	1-1/2" (38 mm)	12" (30 cm)
85811	2" (51 mm)	12" (30 cm)



DID YOU KNOW...

Select **Premier** one-piece pouching systems feature two important product upgrades:

- **Lock 'n Roll** microseal closure with ribbed strips to facilitate ease of use, and transparent flap for increased security (as compared to previous design)
- Upgraded **skin barrier design** to help increase comfort and help improve ease of use:
 - Flexible release liner covers both the skin barrier and the tape border to help provide improved ease of use
 - Tapered barrier edges to help provide increased comfort
 - No change in skin barrier application
 - Same HCPCS code

Premier One-Piece Pouching System

Drainable Pouches with Convex Barriers (continued)

Upgraded Design Recommended replacement for 850xx/852xx series below



821xx/824xx (CUT-TO-FIT)

- ▶ Drainable Pouch
- ▶ *Flexend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ Without Integrated Filter
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige — Both Sides

Box of 5



Upgraded Design Features:

- ▶ *Lock 'n Roll* Microseal Closure with Ribbed Strips and Transparent Flap
- ▶ Flexible Release Liner
- ▶ Tapered Barrier Edges

HCPCS: A4390

Transparent Stock No	For Stomas Up To	Pouch Length
8244	1" (25 mm)	12" (30 cm)
8248	1-1/2" (38 mm)	12" (30 cm)
82411	2" (51 mm)	12" (30 cm)

Original Design

850xx/852xx (CUT-TO-FIT)

- ▶ Drainable Pouch
- ▶ *Flexend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ Without Integrated Filter
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige — Both Sides
- ▶ *Lock 'n Roll* Closure

Box of 5



HCPCS: A4390

Transparent Stock No	Beige Stock No	For Stomas Up To	Pouch Length
8524	8504	1" (25 mm)	12" (30 cm)
8528	8508	1-1/2" (38 mm)	12" (30 cm)
85211	85011	2" (51 mm)	12" (30 cm)

Ostomy Care
Ceramide Infused Products

Ostomy Care
New Image Two-Piece Pouching System

Ostomy Care
Premier One-Piece Pouching System

Ostomy Care
Pouchkins Pediatric Pouching System

Ostomy Care
Ostomy Accessories

Continence Care

Bowel Care

Tube Fasteners

Wound Care

Education

Premier One-Piece Pouching System

Drainable Pouches with Convex Barriers (continued)

Upgraded Design Recommended replacement for 862xx series below



8262xx (CUT-TO-FIT)

- ▶ Drainable Pouch
- ▶ *Flexextend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ Without Integrated Filter
- ▶ Clamp Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels — Body Side Only

**Box of 5
Includes 1 clamp**



Upgraded Design Features:

- ▶ Flexible Release Liner
- ▶ Tapered Barrier Edges

HCPCS: A4390

Transparent Stock No	For Stomas Up To	Pouch Length
82624	1" (25 mm)	12" (30 cm)
82628	1-1/2" (38 mm)	12" (30 cm)
826211	2" (51 mm)	12" (30 cm)

Original Design

862xx (CUT-TO-FIT)

- ▶ Drainable Pouch
- ▶ *Flexextend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ Without Integrated Filter
- ▶ Clamp Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels — Body Side Only

**Box of 5
Includes 1 clamp**



HCPCS: A4390

Transparent Stock No	For Stomas Up To	Pouch Length
8624	1" (25 mm)	12" (30 cm)
8628	1-1/2" (38 mm)	12" (30 cm)
86211	2" (51 mm)	12" (30 cm)

Ostomy Care
 Ceramic Infused Products
Ostomy Care
 New Image Two-Piece Pouching System
Ostomy Care
 Premier One-Piece Pouching System
Ostomy Care
 Pouchkins Pediatric Pouching System
Ostomy Care
 Ostomy Accessories
Continence Care
Bowel Care
Tube Fasteners
Wound Care
Education

Premier One-Piece Pouching System

Drainable Pouches with Convex Barriers (continued)



899x (PRE-SIZED)

- ▶ Drainable Pouch
- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels — Both Sides

Box of 5



HCPCS: A5057

Beige w/Viewing Option

Stock No	Barrier Opening	Pouch Length
8990	3/4" (19 mm)	12" (30 cm)
8991	7/8" (22 mm)	12" (30 cm)
8992	1" (25 mm)	12" (30 cm)
8993	1-1/8" (29 mm)	12" (30 cm)

Ostomy Care
Ceramide Infused
Products

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Ostomy Accessories

Continence
Care

Bowel Care

Tube
Fasteners

Wound Care

Education



DID YOU KNOW...

- **Premier** cut-to-fit products have a 1/2" (13 mm) starter hole and can be used as a pre-sized product for some small stomas (excludes high output pouches, wound drainage collectors and **Pouchkins** products)
- Cutting guides are provided for all cut-to-fit barriers
- Drainable pouches without an integrated closure include one pouch clamp per box
- We do not recommend protective wipes be used with extended wear skin barriers, as they may decrease wear time
- Filtered pouches do not require a cover
- All convex **Premier** pouches have tape border barriers
- All convex **Premier** pouches have belt tabs

Premier One-Piece Pouching System

Drainable Pouches with Convex Barriers (continued)

Ostomy Care
Ceramide Infused
Products

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Ostomy Accessories

Continence
Care

Bowel Care

Tube
Fasteners

Wound Care

Education

Upgraded Design Recommended replacement for 851xx/859xx series below



828x/829x (PRE-SIZED)

Box of 5

- ▶ Drainable Pouch
- ▶ *Flexend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ Without Integrated Filter
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige — Both Sides



Upgraded Design Features:

- ▶ *Lock 'n Roll* Microseal Closure with Ribbed Strips and Transparent Flap
- ▶ Flexible Release Liner
- ▶ Tapered Barrier Edges

HCPCS: A4390

Transparent Stock No	Beige Stock No	Barrier Opening	Pouch Length
8280	8290	3/4" (19 mm)	12" (30 cm)
8281	8291	7/8" (22 mm)	12" (30 cm)
8282	8292	1" (25 mm)	12" (30 cm)
8283	—	1-1/8" (29 mm)	12" (30 cm)
8284	8294	1-1/4" (32 mm)	12" (30 cm)
8285	8295	1-3/8" (35 mm)	12" (30 cm)
8286	—	1-1/2" (38 mm)	12" (30 cm)
8288	—	1-3/4" (44 mm)	12" (30 cm)

Original Design

851xx/859x (PRE-SIZED)

Box of 5

- ▶ Drainable Pouch
- ▶ *Flexend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ Without Integrated Filter
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige — Both Sides
- ▶ *Lock 'n Roll* Closure



HCPCS: A4390

Transparent Stock No	Beige Stock No	Barrier Opening	Pouch Length
8510	8590	3/4" (19 mm)	12" (30 cm)
8511	8591	7/8" (22 mm)	12" (30 cm)
8512	8592	1" (25 mm)	12" (30 cm)
8513	8593	1-1/8" (29 mm)	12" (30 cm)
8514	8594	1-1/4" (32 mm)	12" (30 cm)
8515	8595	1-3/8" (35 mm)	12" (30 cm)
85116	8596	1-1/2" (38 mm)	12" (30 cm)
—	8597	1-5/8" (41 mm)	12" (30 cm)
85118	—	1-3/4" (44 mm)	12" (30 cm)
85119	8599	2" (51 mm)	12" (30 cm)

Premier One-Piece Pouching System

Drainable Pouches with Convex Barriers (continued)

Upgraded Design Recommended replacement for 861xx series below



8261x (PRE-SIZED)

- ▶ Drainable Pouch
- ▶ *Flexextend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ Without Integrated Filter
- ▶ Clamp Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels — Body Side Only

**Box of 5
Includes 1 clamp**



Upgraded Design Features:

- ▶ Flexible Release Liner
- ▶ Tapered Barrier Edges

HCPCS: A4390

Transparent Stock No	Barrier Opening	Pouch Length
82610	3/4" (19 mm)	12" (30 cm)
82611	7/8" (22 mm)	12" (30 cm)
82612	1" (25 mm)	12" (30 cm)
82613	1-1/8" (29 mm)	12" (30 cm)
82614	1-1/4" (32 mm)	12" (30 cm)
82615	1-3/8" (35 mm)	12" (30 cm)
82616	1-1/2" (38 mm)	12" (30 cm)
82617	1-5/8" (41 mm)	12" (30 cm)
82619	2" (51 mm)	12" (30 cm)

Original Design

861x (PRE-SIZED)

- ▶ Drainable Pouch
- ▶ *Flexextend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ Without Integrated Filter
- ▶ Clamp Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels — Body Side Only

**Box of 5
Includes 1 clamp**



HCPCS: A4390

Transparent Stock No	Barrier Opening	Pouch Length
8610	3/4" (19 mm)	12" (30 cm)
8611	7/8" (22 mm)	12" (30 cm)
8612	1" (25 mm)	12" (30 cm)
8613	1-1/8" (29 mm)	12" (30 cm)
8614	1-1/4" (32 mm)	12" (30 cm)
8615	1-3/8" (35 mm)	12" (30 cm)
8616	1-1/2" (38 mm)	12" (30 cm)
8617	1-5/8" (41 mm)	12" (30 cm)
8618	1-3/4" (44 mm)	12" (30 cm)
8619	2" (51 mm)	12" (30 cm)



DID YOU KNOW...

- All high output pouches and wound drainage collectors can be connected directly into a bedside collector without an adapter

Premier One-Piece Pouching System

High Output Drainable Pouches



80070/80110 (CUT-TO-FIT)

Box of 10

- ▶ Non Sterile, Drainable Pouch
- ▶ *Flexextend* (extended wear) Skin Barrier — Flat
- ▶ **Without Tape Border**
- ▶ With Soft Tap
- ▶ With Access Window
- ▶ With *ComfortWear* Panels — Body Side Only
- ▶ Without Starter Hole
- ▶ No adapter needed to connect to bedside collector

HCPCS: A4435

Ultra-Clear Stock No

For Stomas Up To

Pouch Length

80070	2-3/4" (70 mm)	12" (30 cm)
80110	4-1/3" (110 mm)	12" (30 cm)



81070/81110 (CUT-TO-FIT)

Box of 5

- ▶ Sterile, Drainable Pouch
- ▶ *Flexextend* (extended wear) Skin Barrier — Flat
- ▶ **Without Tape Border**
- ▶ With Soft Tap
- ▶ With Access Window
- ▶ With *ComfortWear* Panels — Body Side Only
- ▶ Without Starter Hole
- ▶ No adapter needed to connect to bedside collector

HCPCS: A4435

Ultra-Clear Stock No

For Stomas Up To

Pouch Length

81070	2-3/4" (70 mm)	12" (30 cm)
81110	4-1/3" (110 mm)	12" (30 cm)



970x/977x

- ▶ For Management of Draining Wounds, Fistulae, or Draining Tube Sites
- ▶ Available With and Without *FlexWear* (standard wear) Skin Barriers — Flat
- ▶ Sterile and Non Sterile Drainable Pouches
- ▶ With Access Window
- ▶ Reclosable Cap
- ▶ Without Starter Hole
- ▶ No adapter needed to connect to bedside collector

HCPCS: A6154

With Barrier: Non Sterile Stock No

Description

Quantity/Box

9773	Small (for wounds up to 3")	10
9776	Medium (for wounds up to 3-3/4")	10
9778	Large (for wounds up to 4" x 8")	3

With Barrier: Sterile Stock No

Description

Quantity/Box

9701	Small (for wounds up to 3")	5
9702	Medium (for wounds up to 3-3/4")	5
9703	Large (for wounds up to 4" x 8")	3

Without Barrier: Sterile Stock No

Description

Quantity/Box

9772	Small (for wounds up to 3")	10
9775	Medium (for wounds up to 3-3/4")	10

Premier One-Piece Pouching System

Closed Pouches with Flat Barriers



82300/82400 (CUT-TO-FIT) Box of 30

- ▶ Closed Pouch
- ▶ *SoftFlex* (standard wear) Skin Barrier — Flat
- ▶ **Without Tape Border**
- ▶ With Integrated *AF300* Filter
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige — Both Sides

HCPCS: A4416

Transparent Stock No	Beige Stock No	For Stomas Up To	Pouch Length
82400	82300	2-1/8" (55 mm)	9" (23 cm)



82302/82402 (OVAL CUT-TO-FIT) Box of 30

- ▶ Closed Pouch
- ▶ *SoftFlex* (standard wear) Skin Barrier — Flat
- ▶ **Without Tape Border**
- ▶ With Integrated *AF300* Filter
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige — Both Sides

OVAL

HCPCS: A4416

Transparent Stock No	Beige Stock No	For Stomas Up To	Pouch Length
82402	82302	3" x 2.5" (75 mm x 65 mm)	9" (23 cm)



82100/82500 (CUT-TO-FIT) Box of 30

- ▶ Closed Mini-Pouch
- ▶ *SoftFlex* (standard wear) Skin Barrier — Flat
- ▶ **Without Tape Border**
- ▶ With Integrated *AF300* Filter
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige — Both Sides

HCPCS: A4416

Transparent Stock No	Beige Stock No	For Stomas Up To	Pouch Length
82500	82100	2-1/8" (55 mm)	7" (18 cm)



823xx (PRE-SIZED) Box of 30

- ▶ Closed Pouch
- ▶ *SoftFlex* (standard wear) Skin Barrier — Flat
- ▶ **Without Tape Border**
- ▶ With Integrated *AF300* Filter
- ▶ With *ComfortWear* Panels — Both Sides

HCPCS: A4416

Beige Stock No	Barrier Opening	Pouch Length
82325	1" (25 mm)	9" (23 cm)
82330	1-3/16" (30 mm)	9" (23 cm)
82335	1-3/8" (35 mm)	9" (23 cm)

Premier One-Piece Pouching System

Closed Pouches with Flat Barriers (continued)



821xx (PRE-SIZED)

Box of 30

- ▶ Closed Mini-Pouch
- ▶ *SoftFlex* (standard wear) Skin Barrier — Flat
- ▶ **Without Tape Border**
- ▶ With Integrated *AF300* Filter
- ▶ With *ComfortWear* Panels — Both Sides

HC416: A4416

Beige Stock No	Barrier Opening	Pouch Length
82125	1" (25 mm)	7" (18 cm)
82130	1-3/16" (30 mm)	7" (18 cm)
82135	1-3/8" (35 mm)	7" (18 cm)

Stoma Cap



1796 (PRE-SIZED)

Box of 30

- ▶ Stoma Cap
- ▶ *SoftFlex* (standard wear) Skin Barrier — Flat
- ▶ **Without Tape Border**
- ▶ With Integrated Filter

HCPCS: A5055

Beige Stock No	Barrier Opening	Size
1796	1-15/16" (50 mm)	4" (102 mm)

Premier One-Piece Pouching System

Urostomy Pouches with Flat Barriers

Upgraded Design Recommended replacement for 8440/8450/8460 series below



84490/84590/84690 (CUT-TO-FIT)

Box of 10

Includes 2 adapters

- ▶ Urostomy Pouch
- ▶ Anti-Reflux Valve
- ▶ *Flexextend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With *ComfortWear* Panels
 - Transparent and Ultra-Clear — Body Side Only
 - Beige — Both Sides

Upgraded Design Features:

- ▶ Multi-Chamber Design
- ▶ Softer, Streamlined Tap

HCPCS: A4428

Ultra-Clear Stock No	Transparent Stock No	Beige Stock No	For Stomas Up To	Pouch Length
84590	84690	84490	2-1/2" (64 mm)	9" (23 cm)

Original Design



8440/8450/8460 (CUT-TO-FIT)

Box of 10

Includes 1 adapter

- ▶ Urostomy Pouch
- ▶ Anti-Reflux Valve
- ▶ *Flexextend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With *ComfortWear* Panels
 - Transparent and Ultra Clear — Body Side Only
 - Beige — Both Sides

HCPCS: A4428

Ultra-Clear Stock No	Transparent Stock No	Beige Stock No	For Stomas Up To	Pouch Length
8450	8460	8440	2-1/2" (64 mm)	9" (23 cm)

The upgraded **Hollister** urostomy pouch is available in the **Premier** one-piece pouching system

The upgraded pouch features several innovations designed to deliver more comfort and a greater sense of control (as compared to previous design):

- Skin barrier options to create custom solutions for your patients' unique needs
 - ▶ Flat or convex
 - ▶ Tape border
 - ▶ Flexextend (extended wear) skin barrier
- Enhanced flow-control tap is made with soft material and a streamlined shape
- Multi-chamber design helps to facilitate even distribution of urine
- Improved adapter helps facilitate ease of use
- Belt tab option (convex)

Urostomy Pouches with Flat Barriers (continued)

Upgraded Design Recommended replacement for 844x/846x series below



8449x/8469x (PRE-SIZED)

- ▶ Urostomy Pouch
- ▶ Anti-Reflux Valve
- ▶ *Flexextend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige — Both Sides

Box of 10
Includes 2 adapters

Upgraded Design Features:

- ▶ Multi-Chamber Design
- ▶ Softer, Streamlined Tap

HCPSCS: A4428

Transparent Stock No	Beige Stock No	Barrier Opening	Pouch Length
—	84497	3/4" (19 mm)	9" (23 cm)
84692	—	1" (25 mm)	9" (23 cm)
84698	84498	1-1/4" (32 mm)	9" (23 cm)
84693	—	1-1/2" (38 mm)	9" (23 cm)
84699	84499	1-3/4" (44 mm)	9" (23 cm)
84694	84494	2" (51 mm)	9" (23 cm)

Original Design



844x/846x (PRE-SIZED)

- ▶ Urostomy Pouch
- ▶ Anti-Reflux Valve
- ▶ *Flexextend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige — Both Sides

Box of 10
Includes 1 adapter

HCPSCS: A4428

Transparent Stock No	Beige Stock No	Barrier Opening	Pouch Length
8467	8447	3/4" (19 mm)	9" (23 cm)
8462	8442	1" (25 mm)	9" (23 cm)
8468	8448	1-1/4" (32 mm)	9" (23 cm)
8463	8443	1-1/2" (38 mm)	9" (23 cm)
8469	8449	1-3/4" (44 mm)	9" (23 cm)
8464	—	2" (51 mm)	9" (23 cm)

Premier One-Piece Pouching System

Urostomy Pouches with Convex Barriers

Upgraded Design Recommended replacement for 847x/843x series below



8479xxx/8439xxx (CUT-TO-FIT)

- ▶ Urostomy Pouch
- ▶ Anti-Reflux Valve
- ▶ *Flexextend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels
 - Ultra-Clear — Body Side Only
 - Beige — Both Sides

**Box of 5
Includes 1 adapter**



Upgraded Design Features:

- ▶ Multi-Chamber Design
- ▶ Softer, Streamlined Tap

HCPCS: A4430

Ultra-Clear Stock No	Beige Stock No	For Stomas Up To	Pouch Length
84794	—	1" (25 mm)	9" (23 cm)
84798	—	1-1/2" (38 mm)	9" (23 cm)
8479111	8439111	2" (51 mm)	9" (23 cm)

Original Design



847x/843x (CUT-TO-FIT)

- ▶ Urostomy Pouch
- ▶ Anti-Reflux Valve
- ▶ *Flexextend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige — Both Sides

**Box of 5
Includes 1 adapter**



HCPCS: A4430

Transparent Stock No	Beige Stock No	For Stomas Up To	Pouch Length
8474	8434	1" (25 mm)	9" (23 cm)
8478	8438	1-1/2" (38 mm)	9" (23 cm)
84711	84311	2" (51 mm)	9" (23 cm)



DID YOU KNOW...

- Leaving urine in the urostomy pouch before connecting to a bedside collector/leg bag may minimize the suction effect and allow urine to drain better into the collection system

Premier One-Piece Pouching System

Urostomy Pouches with Convex Barriers (continued)

Upgraded Design Recommended replacement for 848x/849x series below



8489x/8499x (PRE-SIZED)

- ▶ Urostomy Pouch
- ▶ Anti-Reflux Valve
- ▶ *FlexTend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels
 - Ultra-Clear — Body Side Only
 - Beige — Both Sides

**Box of 5
Includes 1 adapter**



Upgraded Design Features:

- ▶ Multi-Chamber Design
- ▶ Softer, Streamlined Tap

HCPCS: A4430

Ultra Clear Stock No	Beige Stock No	Barrier Opening	Pouch Length
84890	—	1/2" (13 mm)	9" (23 cm)
84891	84991	5/8" (16 mm)	9" (23 cm)
84892	84992	3/4" (19 mm)	9" (23 cm)
84893	84993	7/8" (22 mm)	9" (23 cm)
84894	84994	1" (25 mm)	9" (23 cm)
84895	84995	1-1/8" (29 mm)	9" (23 cm)
84896	84996	1-1/4" (32 mm)	9" (23 cm)
84897	84997	1-3/8" (35 mm)	9" (23 cm)
84898	84998	1-1/2" (38 mm)	9" (23 cm)
84899	84999	1-3/4" (44 mm)	9" (23 cm)

Original Design



848x/849x (PRE-SIZED)

- ▶ Urostomy Pouch
- ▶ Anti-Reflux Valve
- ▶ *FlexTend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige — Both Sides

**Box of 5
Includes 1 adapter**



HCPCS: A4430

Transparent Stock No	Beige Stock No	Barrier Opening	Pouch Length
8480	—	1/2" (13 mm)	9" (23 cm)
8481	8491	5/8" (16 mm)	9" (23 cm)
8482	8492	3/4" (19 mm)	9" (23 cm)
8483	8493	7/8" (22 mm)	9" (23 cm)
8484	8494	1" (25 mm)	9" (23 cm)
8485	8495	1-1/8" (29 mm)	9" (23 cm)
8486	8496	1-1/4" (32 mm)	9" (23 cm)
8487	8497	1-3/8" (35 mm)	9" (23 cm)
8488	8498	1-1/2" (38 mm)	9" (23 cm)
8489	8499	1-3/4" (44 mm)	9" (23 cm)

Ostomy Care
Ceramide Infused
Products

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Ostomy Accessories

Continence
Care

Bowel Care

Tube
Fasteners

Wound Care

Education

Premier One-Piece Pouching System

Urostomy Pouches with CeraPlus Soft Convex Barriers

Upgraded Design



8413xx (CUT-TO-FIT)

- ▶ Urostomy Pouch
- ▶ Anti-Reflux Valve
- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels — Body Side Only

Box of 5 includes 1 adapter



NEW

Upgraded Design Features:

- ▶ Multi-Chamber Design
- ▶ Softer, Streamlined Tap
- ▶ Drain Valve Indicator

HCPCS: A4430

Ultra Clear Stock No	For Stomas Up To	Pouch Length
84134	1" (25 mm)	9" (23 cm)
84138	1-1/2" (38 mm)	9" (23 cm)
841311	2-1/8" (55 mm)	9" (23 cm)

Upgraded Design



841x (PRE-SIZED)

- ▶ Urostomy Pouch
- ▶ Anti-Reflux Valve
- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels — Body Side Only

Box of 5 includes 1 adapter



NEW

Upgraded Design Features:

- ▶ Multi-Chamber Design
- ▶ Softer, Streamlined Tap
- ▶ Drain Valve Indicator

HCPCS: A4430

Beige w/Viewing Option Stock No	Barrier Opening	Pouch Length
8412	3/4" (19 mm)	9" (23 cm)
8413	7/8" (22 mm)	9" (23 cm)
8414	1" (25 mm)	9" (23 cm)
8415	1-1/8" (29 mm)	9" (23 cm)

Ostomy Care
Ceramide Infused Products

Ostomy Care
New Image Two-Piece Pouching System

Ostomy Care
Premier One-Piece Pouching System

Ostomy Care
Pouchkins Pediatric Pouching System

Ostomy Care
Ostomy Accessories

Continence Care

Bowel Care

Tube Fasteners

Wound Care

Education

Premier One-Piece Pouching System

Non Sterile Kits — Drainable



89001 (CUT-TO-FIT)

Box of 5 Kits

Colostomy/Ileostomy Kit Includes:

- ▶ (1) 88401 *SoftFlex* (standard wear) Skin Barrier
- ▶ **Without Tape Border**
- ▶ Without Integrated Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ Single-Use Package

HCPCS: A5061

**Ultra-Clear
Stock No**

For Stomas Up To

Pouch Length

89001

2-1/2" (64 mm)

12" (30 cm)



89002 (CUT-TO-FIT)

**Box of 5 Kits
1 clamp per kit**

Colostomy/Ileostomy Kit Includes:

- ▶ (1) 88501 *SoftFlex* (standard wear) Skin Barrier
- ▶ **Without Tape Border**
- ▶ Without Integrated Filter
- ▶ Clamp Closure
- ▶ Single-Use Package

HCPCS: A5061

**Ultra-Clear
Stock No**

For Stomas Up To

Pouch Length

89002

2-1/2" (64 mm)

12" (30 cm)



89003 (CUT-TO-FIT)

**Box of 5 Kits
1 clamp per kit**

Colostomy/Ileostomy Kit Includes:

- ▶ (1) 8651 *Flexend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ Without Integrated Filter
- ▶ Clamp Closure
- ▶ Single-Use Package

HCPCS: A4388

**Ultra-Clear
Stock No**

For Stomas Up To

Pouch Length

89003

2-1/2" (64 mm)

12" (30 cm)



89004 (CUT-TO-FIT)

Box of 5 Kits

Colostomy/Ileostomy Kit Includes:

- ▶ (1) 8531 *Flexend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ Without Integrated Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ Single-Use Package

HCPCS: A4388

**Ultra-Clear
Stock No**

For Stomas Up To

Pouch Length

89004

2-1/2" (64 mm)

12" (30 cm)

Non Sterile Kits — Urostomy



89005 (CUT-TO-FIT)

**Box of 5 Kits
1 adapter per kit**

Urostomy Kit Includes:

- ▶ (1) 8450 *Flexend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ (1) 7331 Urostomy Drain Tube Adapter
- ▶ Single-Use Package

HCPCS: A4428

**Ultra-Clear
Stock No**

For Stomas Up To

Pouch Length

89005

2-1/2" (64 mm)

9" (23 cm)

Education

Wound Care

Tube
Fasteners

Bowel Care

Continence
Care

Ostomy Care
Ostomy Accessories

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Ceramide Infused
Products

A secure fit is necessary for

preventing leakage and

providing patient comfort



Pouchkins

Pediatric Pouching System

Ostomy Care
Ceramide-Infused
Products

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Ostomy Accessories

Continence
Care

Bowel Care

Tube
Fasteners

Wound Care

Education

Education

Wound Care

Tube
Fasteners

Bowel Care

Continence
Care

Ostomy Care
Ostomy Accessories

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Ceramide Infused
Products

Pouchkins Pediatric Pouching System

Premie Pouch



3777 (CUT-TO-FIT)

Box of 15

- ▶ Use as a Premie Closed or Drainable Pouch
- ▶ *SoftFlex* (standard wear) Skin Barrier — Flat
- ▶ **Without Tape Border**
- ▶ Without Starter Hole
- ▶ Without Integrated Filter
- ▶ Soft Wire Closure (if drain opening is desired)

HCPCS: A5052 + A4362

Ultra-Clear Stock No	For Stomas Up To	Pouch Length
3777	5/8" (15 mm)	4-1/2" (11.5 cm)

Newborn Pouch



3778 (CUT-TO-FIT)

Box of 15

- ▶ Newborn Pouch; Use With or Without the Drain Valve
- ▶ *SoftFlex* (standard wear) Skin Barrier — Flat
- ▶ **Without Tape Border**
- ▶ Without Integrated Filter
- ▶ Soft Wire Closure (if larger drain opening is desired)

HCPCS: A5062 + A4362

Ultra-Clear Stock No	For Stomas Up To	Pouch Length
3778	1-3/8" x 7/8" (35 mm x 22 mm)	6-1/2" (16 cm)

One-Piece Pouches



3796 (CUT-TO-FIT)

Box of 10

- ▶ Drainable Pouch
- ▶ *SoftFlex* (standard wear) Skin Barrier — Flat
- ▶ Without Starter Hole
- ▶ **Without Tape Border**
- ▶ Without Integrated Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With *ComfortWear* Panels — Body Side Only

HCPCS: A5061

Transparent Stock No	For Stomas Up To	Pouch Length
3796	1-1/2" (38 mm)	6-1/2" (16 cm)

Pouchkins Pediatric Pouching System

One-Piece Pouches (continued)



3795 (CUT-TO-FIT)

Box of 10

- ▶ Drainable Pouch
- ▶ *SoftFlex* (standard wear) Skin Barrier — Flat
- ▶ Without Starter Hole
- ▶ **Without Tape Border**
- ▶ With Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With *ComfortWear* Panels — Body Side Only

HCPCS: A4424

Transparent Stock No	For Stomas Up To	Pouch Length
3795	2" (51 mm)	7" (18 cm)



3797 (CUT-TO-FIT)

**Box of 10
Includes 1 adapter**

- ▶ Urostomy Pouch
- ▶ Anti-Reflux Valve
- ▶ *SoftFlex* (standard wear) Skin Barrier — Flat
- ▶ Without Starter Hole
- ▶ **Without Tape Border**
- ▶ With *ComfortWear* Panels — Body Side Only

HCPCS: A5071

Transparent Stock No	For Stomas Up To	Pouch Length
3797	1-1/2" (38 mm)	8-3/4" (22 cm)



DID YOU KNOW...

Pouchkins SoftFlex skin barrier

- Designed to be gentle
- No starter hole, which enables the opening to be placed where needed within the cutting guide
- No tape by design
- Flexible and sized for babies and children
- None of the Pouchkins pediatric pouching systems are made with natural rubber latex
- Any 1-3/4" New Image pouch will fit this skin barrier



Pouchkins Pediatric Pouching System

Two-Piece Skin Barrier



3761 (CUT-TO-FIT) **Box of 5**

- ▶ *SoftFlex* (standard wear) Skin Barrier — Flat
- ▶ Without Starter Hole
- ▶ Without Tape Border

HCPCS: A4414

Transparent Stock No	Flange Size	For Stomas Up To
3761	 1-3/4" (44 mm)	1-1/4" (32 mm)

Two-Piece Pouches



3799 **Box of 10**

- ▶ Drainable Pouch
- ▶ Without Integrated Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panel — Body Side Only

HCPCS: A5063

Transparent Stock No	Flange Size	Pouch Length
3799	 1-3/4" (44 mm)	6-1/2" (16 cm)



18292/18282 **Box of 20**

- ▶ Drainable Pouch
- ▶ With Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panel
 - Transparent — Body Side Only
 - Beige — Both Sides

HCPCS: A4425

Transparent Stock No	Beige Stock No	Flange Size	Pouch Length
18292	18282	 1-3/4" (44 mm)	7" (18 cm)



18382/18392 **Box of 60**

- ▶ Closed Mini-Pouch
- ▶ With Integrated *AF300* Filter
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige — Both Sides

HCPCS: A4419

Transparent Stock No	Beige Stock No	Flange Size	Pouch Length
18382	18392	 1-3/4" (44 mm)	7" (18 cm)

Pouchkins Pediatric Pouching System

Two-Piece Pouches (continued)

Upgraded Design Recommended replacement for 1840x/1841x/1842xx series below



1890x/1891x/1892x

- ▶ Urostomy Pouch
- ▶ Anti-Reflux Valve
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels
 - Transparent and Ultra-Clear — Body Side Only
 - Beige — Both Sides

**Box of 10
Includes 2 adapters**

Upgraded Design Features:

- ▶ Multi-Chamber Design
- ▶ Softer, Streamlined Tap
- ▶ Drain Valve Indicator

HCPCS: A4432

Ultra-Clear Stock No	Transparent Stock No	Beige Stock No	Flange Size	Pouch Length
18922	18902	18912	 1-3/4" (44 mm)	9" (23 cm)

Original Design



18402/18412/18422

- ▶ Urostomy Pouch
- ▶ Anti-Reflux Valve
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels
 - Transparent and Ultra-Clear — Body Side Only
 - Beige — Both Sides

**Box of 10
Includes 1 adapter**

HCPCS: A4432

Ultra-Clear Stock No	Transparent Stock No	Beige Stock No	Flange Size	Pouch Length
18422	18402	18412	 1-3/4" (44 mm)	9" (23 cm)

Pediatric Ostomy Belt



3774

- ▶ Ostomy Belt
- White Color

Single Unit

HCPCS: A4367

Stock No	Length
3774	10"- 17" (25 cm-43 cm)



DID YOU KNOW...

Shadow Buddies dolls

- Exclusive to **Hollister**, Shadow Buddies dolls are special friends for children who have had or will have ostomy surgery
- Great teaching tool to help with education and care giving
- Huggable dolls with heart-shaped eyes and a stoma provide comfort and support



To order, refer to Shadow Buddies information in the Ostomy Education Resources section.

Shadow Buddies is a trademark of The Shadow Buddies Foundation

Ostomy Care
Ceramide-Infused
Products

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Ostomy Accessories

Continence
Care

Bowel Care

Tube
Fasteners

Wound Care

Education

Accessories can enhance

performance of the skin barrier

and promote healthy skin



Ostomy Accessories

Why use an ostomy accessory?

Ostomy accessories can be used together with your pouching system. You may use one or more of these products. Here are a few reasons why you might want to use an ostomy accessory:

- It may help the performance of your pouching system. For example, an accessory could help increase wear time, protect your skin, or make the removal of your adhesive easier
- It may help solve a problem such as leakage or skin irritation
- It may help improve the quality of your life by eliminating odor

Count on the Right Fit With

Adapt Barrier Rings



Hollister understands that no two people or stomas are alike — and that barrier rings can play a critical role in helping ensure the best outcome for your patients

Whether you choose flat or convex, slim or regular, round or oval— Adapt barrier rings:

- Are designed to be stretched or shaped
- Help prevent ostomy output from getting under the pouch seal, which may reduce leakage and skin irritation
- Can be cut or stacked together to improve the fit of the pouching system
- Are alcohol free, so there is no sting from alcohol when applied to broken or irritated skin

Talk to your Hollister ostomy representative about the full line of Adapt barrier rings, now available with and without ceramide.



Ostomy Accessories

Adapt Skin Barrier Rings



88xx

- ▶ *CeraRing* Barrier Rings (Flat)
- ▶ *Flexlend M* (extended wear) with Ceramide

Box of 10

HCPCS: A4385

Stock No	Outer Diameter	Thickness
8805	2" (48 mm)	4.5 mm
8815	2" (48 mm)	2.3 mm SLIM



78xx

- ▶ Barrier Rings (Flat)
- ▶ *Flexlend M* (extended wear)

Box of 10

HCPCS: A4385

Stock No	Outer Diameter	Thickness
7805	2" (48 mm)	4.5 mm
7815	2" (48 mm)	2.3 mm SLIM
7806	4" (98 mm)	2.3 mm SLIM



895xx

- ▶ *CeraRing* Convex Barrier Rings
- ▶ *Flexlend* (extended wear) with Ceramide

Box of 10

CONVEX

NEW

HCPCS: A4411

Stock No	Inner Diameter
89520	13/16" (20 mm) — can be stretched to 1" (25 mm)
89530	1-3/16" (30 mm) — can be stretched to 1-3/8" (35 mm)
89540	1-9/16" (40 mm) — can be stretched to 1-3/4" (45 mm)



795xx

- ▶ Convex Barrier Rings
- ▶ *Flexlend* (extended wear)

Box of 10

CONVEX

HCPCS: A4411

Stock No	Inner Diameter
79520	13/16" (20 mm) — can be stretched to 1" (25 mm)
79530	1-3/16" (30 mm) — can be stretched to 1-3/8" (35 mm)
79540	1-9/16" (40 mm) — can be stretched to 1-3/4" (45 mm)



896xx

- ▶ *CeraRing* Oval Convex Barrier Rings
- ▶ *Flexlend* (extended wear) with Ceramide

Box of 10

CONVEX

NEW

HCPCS: A4411

Stock No	Inner Diameter
89601	7/8" x 1-1/2" (22 mm x 38 mm) — can be stretched to 1-1/8" x 1-3/4" (27 mm x 43 mm)
89602	1-3/16" x 1-7/8" (30 mm x 48 mm) — can be stretched to 1-3/8" x 2-1/8" (35 mm x 53 mm)
89603	1-1/2" x 2-3/16" (38 mm x 56 mm) — can be stretched to 1-3/4" x 2-3/8" (43 mm x 61 mm)



796xx

- ▶ Oval Convex Barrier Rings
- ▶ *Flexlend* (extended wear)

Box of 10

CONVEX

HCPCS: A4411

Stock No	Inner Diameter
79601	7/8" x 1-1/2" (22 mm x 38 mm) — can be stretched to 1-1/8" x 1-3/4" (27 mm x 43 mm)
79602	1-3/16" x 1-7/8" (30 mm x 48 mm) — can be stretched to 1-3/8" x 2-1/8" (35 mm x 53 mm)
79603	1-1/2" x 2-3/16" (38 mm x 56 mm) — can be stretched to 1-3/4" x 2-3/8" (43 mm x 61 mm)



A Sense of Security

Means Everything



Three strips create a custom frame around your skin barrier.

Adapt barrier extenders

At Hollister Ostomy Care, we know a sense of security and leakage prevention are top priorities for you.

Our Adapt barrier extenders offer an additional layer of protection for those concerned about lifting around the edges of their skin barrier. With a thin and flexible hydrocolloid formulation, Adapt barrier extenders conform easily to body contours. Three convenient strips make it easy to create a custom frame around the skin barrier. The end result is enhanced security, when it matters most.



Ostomy Care
Ceramide-Infused
Products

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Ostomy Accessories

Continence
Care

Bowel Care

Tube
Fasteners

Wound Care

Education

Ostomy Accessories

Adapt Skin Barrier Strips



79400

- ▶ Skin Barrier Strips
- ▶ *Flexextend* (extended wear)

HCPCS: A4385

Box of 10

Stock No	Size	Length
79400	60 g	5-1/8" (13 cm)

Adapt Barrier Extenders



79402

- ▶ Barrier Extenders
- ▶ Hydrocolloid Formulation

HCPCS: A4362*

Box of 30

Stock No
79402

* 3 extenders = 1 unit, Box of 30 = 10 billable units

Adapt Skin Barrier Paste



79300/79301

- ▶ Non-Adhesive

HCPCS: A4406

Stock No	Description	Quantity
79300	2.1 oz (60 g) tube	Box of 1
79301	0.5 oz (14 g) tube	Box of 20

Adapt Stoma Powder



7906

Box of 1

HCPCS: A4371

Stock No	Description
7906	1 oz (28.3 g) puff bottle

Adapt Odor Eliminator and Lubricant



78500/78501

HCPCS Code: A4394

Stock No	Description	Quantity
78500	8 oz (236 mL) bottle	Box of 1
78501	0.27 oz (8 mL) packet	50 Packets/Box

Adapt Ostomy Belts



7300/7299

Box of 10

HCPCS: A4367

Beige Stock No	Size
7300	23"-43" (58 cm-109 cm) Medium
7299	34"-65" (86 cm-165 cm) Large

Ostomy Care
Ceramide Infused Products

Ostomy Care
New Image Two-Piece Pouching System

Ostomy Care
Premier One-Piece Pouching System

Ostomy Care
Pouchkins Pediatric Pouching System

Ostomy Care
Ostomy Accessories

Continence Care

Bowel Care

Tube Fasteners

Wound Care

Education

Adapt Skin Prepping Product



- 7917**
- ▶ Skin Protective Wipes
 - ▶ Non-Adhesive
 - ▶ Alcohol Free

HCPCS: A5120

**1 Box
50 Wipes Each**

NO STING

Stock No

7917

*Not recommended for use with extended wear skin barriers, as it may decrease wear time

Adapt Adhesives and Removers



- 7760**
- ▶ Universal Remover Wipes
 - ▶ For Adhesives and Barriers
 - ▶ Silicone Based
 - ▶ Alcohol Free

HCPCS: A4456

**1 Box
50 Wipes Each**

NO STING

Stock No

7760



- 7737**
- ▶ Medical Adhesive Remover
 - ▶ Silicone Based
 - ▶ Alcohol Free

HCPCS: A4455

Box of 1

NO STING

Stock No

7737

Description

1.7 oz (50 mL) 360 degree spray can



- 7730**
- ▶ Medical Adhesive
 - ▶ Improves Adhesive Contact Between Skin and Barrier

HCPCS: A4364

Box of 1

Stock No

7730

Description

3.2 oz (90 g) spray can

Flextend (extended wear) Skin Barrier



- 8800**
8801
- HCPCS for 8800: A4385**
HCPCS for 8801: A5122

Stock No

8800

4 x 4" (10 cm x 10 cm)

Box of 5

8801

8 x 8" (20 cm x 20 cm)

Box of 3

Ostomy Accessories

Odor Eliminators



773x

- ▶ *m9* Odor Eliminator Spray

HCPCS: No Code

Stock No (unscented)	Description	Quantity
7732	2 oz (60 mL) pump spray	Box of 12
7733	8 oz (240 mL) pump spray	Box of 6

Stock No (Green Apple Scented)	Description	Quantity
7734	2 oz (60 mL) pump spray	Box of 12
7735	8 oz (240 mL) pump spray	Box of 6



771x

- ▶ *m9* Odor Eliminator Drops

HCPCS: A4394

Stock No	Description	Quantity
7715	1 oz (30 mL) bottle	Box of 12
7717	8 oz (240 mL) bottle	Box of 6

Cleaner/Decrystallizer



7736

- ▶ *m9* Cleaner/Decrystallizer
- ▶ Cleans Urinary Drainage Systems
- ▶ pH Balanced and Non-Acidic
- ▶ Free Wash Bottle Included

Box of 12

HCPCS: A5131

Stock No	Description
7736	16 oz (480 mL) bottle

Drainable Pouch Clamp



8770

HCPCS: A4363

Box of 20

Stock No	Color
8770	Beige

Urostomy Drain Tube Adapter



7331

HCPCS: A4421

Box of 10

Stock No
7331

Hollister Universal Catheter Access Port



9779

- ▶ Graduated nipple and punch allows introduction of catheter or drain through pouch
- ▶ Sterile

Box of 8

Stock No
9779

Ostomy Accessories

Other Accessories

Loop Ostomy Bridge



7767
▶ Sterile

Box of 20

HCPCS: A4421

Stock No

7767

Bedside Drainage Collection System with Anti-Reflux Valve



9839
▶ Sample Port
▶ Sterile Fluid Pathway

Box of 20

HCPCS: A4357

Stock No

9839

Description

68 oz (2000 mL) Capacity

Skin Barrier Pastes and Powders



7910
▶ Karaya 5 Paste

Box of 12

HCPCS: A4405

Stock No

7910

Description

4.5 oz (135 mL) tube



7905
▶ Karaya 5 Powder

Box of 12

HCPCS: A4371

Stock No

7905

Description

2.5 oz (75 mL) puff bottle

Skin Barriers



7800
▶ Premium (standard wear) Skin Barrier

Box of 5

HCPCS: A4362

Stock No

7800

Size

4 x 4" (10 cm x 10 cm)



7700
▶ *HolliHesive* (standard wear) Skin Barrier

Box of 5

HCPCS: A4362

Stock No

7700

Size

4 x 4" (10 cm x 10 cm)

Replacement Filters



7766
▶ For use with 366x, 383x and 389x
▶ 1 pack contains 20 Filters and 2 Plugs

Box of 5 Packs

HCPCS: A4368


Stock No

7766

Irrigation Products

General Information

- ▶ Similar to an enema with the water being placed into the stoma.
- ▶ Routine irrigation is recommended to be done on doctor's orders (to regulate emptying the bowel at a scheduled time).
- ▶ Appropriate for some people with sigmoid or descending colostomies.
- ▶ Not appropriate for people with urostomy or ileostomy.
- ▶ Allow approximately one hour for irrigation procedure.



7721

Box of 1

- ▶ Cone Irrigator Kit Includes:
 - ▶ (1) Irrigator
 - ▶ (1) Stoma Cone with Connector
 - ▶ (1) Cleaning Brush

HCPCS: A4398 + A4399

Stock No

7721



7723

Box of 10

- ▶ Stoma Cone with Connector
- ▶ Replacement Unit

HCPCS: A4399

Stock No

7723



1460x (CUT-TO-FIT)

Box of 5

- ▶ *Flexextend* (extended wear) Skin Barrier
- ▶ With Tape Border

HCPCS: A4409

Stock No

Flange Size

For Stomas Up To

14602	 1-3/4" (44 mm)	1-1/4" (32 mm)
14603	 2-1/4" (57 mm)	1-3/4" (44 mm)
14604	 2-3/4" (70 mm)	2-1/4" (57 mm)



1815x

Box of 5

- ▶ Drainable Irrigator Sleeve to be used with *New Image* Two-Piece System
- ▶ With Belt Tabs
- ▶ Includes Clamp

HCPCS: A4397

Transparent Stock No

Flange Size

Sleeve Length

18152	 1-3/4" (44 mm)	30" (76 cm)
18153	 2-1/4" (57 mm)	30" (76 cm)
18154	 2-3/4" (70 mm)	30" (76 cm)



7728/7724

Box of 20

- ▶ Irrigator Sleeve with Belt Tabs
- ▶ Includes Clamp

HCPCS: A4397

Stock No	Opening Size	Sleeve Length
7728	2" (51 mm)	22" (56 cm)
7724	3" (76 mm)	22" (56 cm)



318x (PRE-SIZED)

Box of 30

- ▶ Stoma Cap, Tape Adhesive, Integrated Filter, Odor-Barrier Film, Without *ComfortWear* Panels

HCPCS: A5055

Beige Stock No	Opening	Size
3184	2" (51 mm)	4-1/4" (11 cm)
3186	3" (76 mm)	4-1/4" (11 cm)



7300/7299

Box of 10

- ▶ Ostomy Belts

HCPCS: A4367

Beige Stock No	Size
7300	23"-43" (58 cm-109 cm) Medium
7299	34"-65" (86 cm-165 cm) Large



7740

Box of 10

- ▶ *Adapt* Stoma Lubricant
- ▶ For Use with Stoma Irrigator

HCPCS: A4402

Stock No	Description
7740	4 oz (120 mL) bottle

Upgrade to the New Image Two-Piece Pouching System!

CenterPointLock Two-Piece Pouching System

At Hollister, we are dedicated to providing a choice of quality, differentiated products, with features that meet the individual needs of clinicians and patients.

Please be advised that the CenterPointLock two-piece pouching system will be discontinued effective December 31, 2018.

The CenterPointLock two-piece ostomy system was developed and launched over three decades ago. Since then, our products have undergone advancements for the benefit of our customers, including improvements in both skin barrier and pouching technology.

We are encouraging people who are using the CenterPointLock two-piece ostomy system to upgrade to the New Image two-piece system so they can experience the security and quality of the new generation Hollister two-piece system!

[Click below to see conversion charts](#)

[Skin Barrier Conversion](#)

[Pouch Conversion](#)





Integrated Floating Flange

- Audible clicks signal a secure attachment
- Discreet low profile attachment of pouch to barrier



CeraPlus Skin Barrier

The CeraPlus skin barrier with Remois technology* is infused with ceramide, the skin's naturally occurring protection against dryness. It is designed to maintain adhesive properties, and features a formulation infused with ceramide to help:

- Protect the skin's natural moisture barrier
- Maintain good peristomal skin health
- Decrease transepidermal water loss (TEWL) from damaged or eroded skin

**Remois technology is a technology of Alcare Co., Ltd.*



New Image Pouches

- Waterproof, integrated AF300 filter is designed to allow gas and air to escape, but keeps odor and output in the pouch. Liquid cannot penetrate the filter, therefore does not require a protective sticker
- Choice of secure pouch closure options – clamp or integrated closure
- Comfortwear panels provide enhanced discretion and comfort
- Numerous pouch size options to meet individual preferences

CenterPointLock Two-Piece Pouching System

Letter Matching System

To help you select compatible CenterPointLock two-piece system components, every CenterPointLock system product box label features a letter designation: Ⓜ, Ⓛ, Ⓤ, Ⓚ, or Ⓦ. Once you've selected a CenterPointLock two-piece system skin barrier, select a pouch, irrigator drain, or stoma cap with a matching letter on the box label.

Flat Skin Barriers with Floating Flange



872x (CUT-TO-FIT)

Box of 5

▶ *Flexend* (extended wear) Skin Barrier, Tape Border

HCPCS: A4409

Stock No	Letter Match	Flange Size	For Stomas Up To
8727	Ⓜ	1-1/2" (38 mm)	1" (25 mm)
8722	Ⓛ	1-3/4" (44 mm)	1-1/4" (32 mm)
8723	Ⓤ	2-1/4" (57 mm)	1-3/4" (44 mm)
8724	Ⓚ	2-3/4" (70 mm)	2-1/4" (57 mm)



875x (CUT-TO-FIT)

Box of 5

▶ *Flexend* (extended wear) Skin Barrier, Without Tape Border

HCPCS FOR 8752, 8753 & 8757: A4409

HCPCS FOR 8754: A4410

Stock No	Letter Match	Flange Size	For Stomas Up To
8752	Ⓛ	1-3/4" (44 mm)	1-1/4" (32 mm)
8753	Ⓤ	2-1/4" (57 mm)	1-3/4" (44 mm)
8754	Ⓚ	2-3/4" (70 mm)	2-1/4" (57 mm)



372x (CUT-TO-FIT)

Box of 5

▶ *FlexWear* (standard wear) Skin Barrier, Tape Border

HCPCS FOR 3722-3724 & 3727: A4414

HCPCS FOR 3726: A4415

Stock No	Letter Match	Flange Size	For Stomas Up To
3727	Ⓜ	1-1/2" (38 mm)	1" (25 mm)
3722	Ⓛ	1-3/4" (44 mm)	1-1/4" (32 mm)
3723	Ⓤ	2-1/4" (57 mm)	1-3/4" (44 mm)
3724	Ⓚ	2-3/4" (70 mm)	2-1/4" (57 mm)
3726	Ⓦ	4" (102 mm)	3-1/2" (89 mm)



376x (CUT-TO-FIT)

Box of 5

▶ *SoftFlex* (standard wear) Skin Barrier, Without Tape Border

HCPCS FOR 3762, 3763 & 3767: A4414

HCPCS FOR 3764 & 3766: A4415

Stock No	Letter Match	Flange Size	For Stomas Up To
3767	Ⓜ	1-1/2" (38 mm)	1" (25 mm)
3762	Ⓛ	1-3/4" (44 mm)	1-1/4" (32 mm)
3763	Ⓤ	2-1/4" (57 mm)	1-3/4" (44 mm)
3764	Ⓚ	2-3/4" (70 mm)	2-1/4" (57 mm)
3766	Ⓦ	4" (102 mm)	3-1/2" (89 mm)

Note: **New Image** and **CenterPointLock** pouching systems are not interchangeable.

Ostomy Care
Ceramide Infused Products
Ostomy Care
New Image Two-Piece Pouching System
Ostomy Care
Premier One-Piece Pouching System
Ostomy Care
Pouchkins Pediatric Pouching System
Ostomy Care
Ostomy Accessories
Continence Care
Bowel Care
Tube Fasteners
Wound Care
Education

CenterPointLock Two-Piece Pouching System

Flat Skin Barriers with Floating Flange (continued)



370x (CUT-TO-FIT)

Box of 5

▶ *HolliHesive* (standard wear) Skin Barrier, Tape Border

HCPCS FOR 3702-3704 & 3707: A4414

HCPCS FOR 3706: A4415

Stock No	Letter Match	Flange Size	For Stomas Up To
3707	Ⓟ	1-1/2" (38 mm)	up to 3/4" (19 mm)
3702	Ⓛ	1-3/4" (44 mm)	up to 1" (25 mm)
3703	Ⓧ	2-1/4" (57 mm)	up to 1-1/2" (38 mm)
3704	Ⓢ	2-3/4" (70 mm)	up to 2" (51 mm)
3706	Ⓛ	4" (102 mm)	up to 3-1/4" (82 mm)



874x (PRE-SIZED)

Box of 5

▶ *Flexend* (extended wear) Skin Barrier, Tape Border

HCPCS: A4409

Stock No	Letter Match	Flange Size	Barrier Opening
8743	Ⓟ	1-1/2" (38 mm)	7/8" (22 mm)
8744	Ⓟ	1-1/2" (38 mm)	1" (25 mm)
8745	Ⓛ	1-3/4" (44 mm)	1-1/8" (29 mm)
8746	Ⓛ	1-3/4" (44 mm)	1-1/4" (32 mm)
8747	Ⓧ	2-1/4" (57 mm)	1-3/8" (35 mm)
8748	Ⓧ	2-1/4" (57 mm)	1-1/2" (38 mm)



374x (PRE-SIZED)

Box of 5

▶ *FlexWear* (standard wear) Skin Barrier, Tape Border

HCPCS: A4414

Stock No	Letter Match	Flange Size	Barrier Opening
3747	Ⓟ	1-1/2" (38 mm)	3/4" (19 mm)
3742	Ⓟ	1-1/2" (38 mm)	1" (25 mm)
3748	Ⓛ	1-3/4" (44 mm)	1-1/4" (32 mm)
3743	Ⓧ	2-1/4" (57 mm)	1-1/2" (38 mm)
3749	Ⓧ	2-1/4" (57 mm)	1-3/4" (44 mm)

Convex Skin Barriers with Floating Flange



877x (CUT-TO-FIT)

Box of 5

▶ *Flexend* (extended wear) Skin Barrier, Convex, Tape Border



HCPCS: A4407

Stock No	Letter Match	Flange Size	Barrier Opening
8774	Ⓛ	1-3/4" (44 mm)	1" (25 mm)
8778	Ⓧ	2-1/4" (57 mm)	1-1/2" (38 mm)
87711	Ⓢ	2-3/4" (70 mm)	2" (51 mm)



873x (PRE-SIZED)

Box of 5

▶ *Flexend* (extended wear) Skin Barrier, Convex, Tape Border



HCPCS: A4407

Stock No	Letter Match	Flange Size	Barrier Opening
8731	Ⓟ	1-1/2" (38 mm)	5/8" (16 mm)
8732	Ⓟ	1-1/2" (38 mm)	3/4" (19 mm)
8733	Ⓛ	1-3/4" (44 mm)	7/8" (22 mm)
8734	Ⓛ	1-3/4" (44 mm)	1" (25 mm)
8735	Ⓧ	2-1/4" (57 mm)	1-1/8" (29 mm)
8736	Ⓧ	2-1/4" (57 mm)	1-1/4" (32 mm)
8737	Ⓧ	2-1/4" (57 mm)	1-3/8" (35 mm)
8738	Ⓧ	2-1/4" (57 mm)	1-1/2" (38 mm)

CenterPointLock Two-Piece Pouching System

Convex Skin Barriers with Floating Flange (continued)



373x (PRE-SIZED)

- ▶ *FlexWear* (standard wear) Skin Barrier, Convex, Tape Border

Box of 5



HCPCS: A4373

Stock No	Letter Match	Flange Size	Barrier Opening
3730	Ⓜ	1-1/2" (38 mm)	1/2" (13 mm)
3731	Ⓜ	1-1/2" (38 mm)	5/8" (16 mm)
3732	Ⓜ	1-1/2" (38 mm)	3/4" (19 mm)
3733	Ⓛ	1-3/4" (44 mm)	7/8" (22 mm)
3734	Ⓛ	1-3/4" (44 mm)	1" (25 mm)
3735	Ⓛ	2-1/4" (57 mm)	1-1/8" (29 mm)
3736	Ⓛ	2-1/4" (57 mm)	1-1/4" (32 mm)
3737	Ⓛ	2-1/4" (57 mm)	1-3/8" (35 mm)
3738	Ⓛ	2-1/4" (57 mm)	1-1/2" (38 mm)
3739	Ⓚ	2-3/4" (70 mm)	1-5/8" (41 mm)
37310	Ⓚ	2-3/4" (70 mm)	1-3/4" (44 mm)
37311	Ⓚ	2-3/4" (70 mm)	2" (51 mm)

Drainable Pouches



389x/383x

Box of 10

- ▶ Drainable Pouch, Replaceable Filter, Clamp Closure, Belt Tabs, Without *ComfortWear* Panels. Includes 1 Bag of 20 *Hollister* Replacement Filters and 2 Plugs, #7766

HCPCS: A4425 + A4368

Transparent Stock No	Beige Stock No	Letter Match	Flange Size	Pouch Length
3897	3837	Ⓜ	1-1/2" (38 mm)	12" (30 cm)
3892	3832	Ⓛ	1-3/4" (44 mm)	12" (30 cm)
3893	3833	Ⓛ	2-1/4" (57 mm)	12" (30 cm)
3894	3834	Ⓚ	2-3/4" (70 mm)	12" (30 cm)



881x

Box of 10

- ▶ Drainable Pouch, Without Integrated Filter, Clamp Closure, Belt Tabs, *ComfortWear* Panels — Both Sides

HCPCS: A5063

Beige Stock No	Letter Match	Flange Size	Pouch Length
8817	Ⓜ	1-1/2" (38 mm)	12" (30 cm)
8812	Ⓛ	1-3/4" (44 mm)	12" (30 cm)
8813	Ⓛ	2-1/4" (57 mm)	12" (30 cm)
8814	Ⓚ	2-3/4" (70 mm)	12" (30 cm)



380x/381x

Box of 10

- ▶ Drainable Pouch, Without Integrated Filter, Clamp Closure, Belt Tabs, Without *ComfortWear* Panels

HCPCS: A5063

Transparent Stock No	Beige Stock No	Letter Match	Flange Size	Pouch Length
3807	3817	Ⓜ	1-1/2" (38 mm)	11" (28 cm)
3802	3812	Ⓛ	1-3/4" (44 mm)	11" (28 cm)
3803	3813	Ⓛ	2-1/4" (57 mm)	11" (28 cm)
3804	3814	Ⓚ	2-3/4" (70 mm)	11" (28 cm)
3806	—	Ⓛ	4" (102 mm)	11" (28 cm)

CenterPointLock Two-Piece Pouching System

Drainable Pouches (continued)



387x Box of 10

- ▶ High Output Drainable Pouch, Without Integrated Filter, Clamp Closure, Belt Tabs, Without *ComfortWear* Panels

HCPCS: A4412

Transparent Stock No	Letter Match	Flange Size	Pouch Length
3872	Ⓡ	1-3/4" (44 mm)	16" (40 cm)
3873	Ⓢ	2-1/4" (57 mm)	16" (40 cm)
3874	Ⓣ	2-3/4" (70 mm)	16" (40 cm)



384x Box of 10

- ▶ Drainable Mini-Pouch, Without Integrated Filter, Clamp Closure, Belt Tabs, Without *ComfortWear* Panels

HCPCS: A5063

Beige Stock No	Letter Match	Flange Size	Pouch Length
3847	Ⓡ	1-1/2" (38 mm)	9" (23 cm)
3842	Ⓢ	1-3/4" (44 mm)	9" (23 cm)
3843	Ⓣ	2-1/4" (57 mm)	9" (23 cm)
3844	Ⓤ	2-3/4" (70 mm)	9" (23 cm)

Closed Pouches



834x Box of 15

- ▶ Closed Pouch, Integrated Filter, Odor-Barrier Rustle-Free Film, Belt Tabs, *ComfortWear* Panels — Both Sides

HCPCS: A4419

Beige Stock No	Letter Match	Flange Size	Pouch Length
8347	Ⓡ	1-1/2" (38 mm)	7-1/2" (19 cm)
8342	Ⓢ	1-3/4" (44 mm)	7-1/2" (19 cm)
8343	Ⓣ	2-1/4" (57 mm)	7-1/2" (19 cm)
8344	Ⓤ	2-3/4" (70 mm)	7-1/2" (19 cm)



334x Box of 15

- ▶ Closed Pouch, Integrated Filter, Belt Tabs, Without *ComfortWear* Panels

HCPCS: A4419

Beige Stock No	Letter Match	Flange Size	Pouch Length
3347	Ⓡ	1-1/2" (38 mm)	9" (23 cm)
3342	Ⓢ	1-3/4" (44 mm)	9" (23 cm)
3343	Ⓣ	2-1/4" (57 mm)	9" (23 cm)
3344	Ⓤ	2-3/4" (70 mm)	9" (23 cm)

Urostomy Pouches



890x Box of 10

- ▶ Urostomy Pouch, Belt Tabs, *ComfortWear* Panels — Both Sides

HCPCS: A4432

Beige Stock No	Letter Match	Flange Size	Pouch Length
8902	Ⓡ	1-3/4" (44 mm)	9" (23 cm)
8903	Ⓢ	2-1/4" (57 mm)	9" (23 cm)

CenterPointLock Two-Piece Pouching System

Urostomy Pouches (continued)



390x

- ▶ Urostomy Pouch, Belt Tabs, Without *ComfortWear* Panels

Box of 10

HCPCS: A4432

Stock No	Transparent	Letter Match	Flange Size	Pouch Length
3907		H	1-1/2" (38 mm)	9" (23 cm)
3902		I	1-3/4" (44 mm)	9" (23 cm)
3903		J	2-1/4" (57 mm)	9" (23 cm)
3904		K	2-3/4" (70 mm)	9" (23 cm)

Stoma Cap



340x

- ▶ Stoma Cap, Integrated Filter, Without *ComfortWear* Panels

Box of 25

HCPCS: A5055

Beige Stock No	Letter Match	Flange Size	Pouch Length
3407	H	1-1/2" (38 mm)	4-1/2" (11 cm)
3402	I	1-3/4" (44 mm)	4-1/2" (11 cm)
3403	J	2-1/4" (57 mm)	4-1/2" (11 cm)
3404	K	2-3/4" (70 mm)	4-1/2" (11 cm)

Irrigator Sleeve



382x

- ▶ Irrigator Sleeve, Belt Tabs
- ▶ Includes clamp

Box of 5

HCPCS: A4397

Transparent Stock No	Letter Match	Flange Size	Pouch Length
3827	H	1-1/2" (38 mm)	30" (76 cm)
3822	I	1-3/4" (44 mm)	30" (76 cm)
3823	J	2-1/4" (57 mm)	30" (76 cm)
3824	K	2-3/4" (70 mm)	30" (76 cm)
3826	L	4" (102 mm)	30" (76 cm)

Ostomy Care
Ceramide Infused
Products

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Ostomy Accessories

Continence
Care

Bowel Care

Tube
Fasteners

Wound Care

Education

Karaya 5 One-Piece Pouching System

Drainable Pouches



366x (PRE-SIZED)

Box of 10

- ▶ Drainable Pouch, *Karaya 5* (standard wear)
Skin Barrier — Convex, Tape Border, Replaceable Filter, Pouch Clamp Closure, Belt Tabs, Without *ComfortWear* Panels.
Includes 1 Bag of 20 *Hollister* Replacement Filters and 2 Plugs, #7766

HPCPS: A4389 + A4368



Transparent Stock No	Pouch Size	Barrier Opening	Pouch Length
3662	1" (25 mm)	5/8" (16 mm)	12" (30 cm)
3668	1-1/4" (32 mm)	7/8" (22 mm)	12" (30 cm)
3663	1-1/2" (38 mm)	1-1/8" (29 mm)	12" (30 cm)
3669	1-3/4" (44 mm)	1-3/8" (35 mm)	12" (30 cm)
3664	2" (51 mm)	1-1/2" (38 mm)	12" (30 cm)
3665	2-1/2" (64 mm)	2" (51 mm)	12" (30 cm)



360x (PRE-SIZED)

Box of 10

- ▶ Drainable Pouch, *Karaya 5* (standard wear)
Skin Barrier — Convex, Tape Border, Without Integrated Filter, Pouch Clamp Closure, Belt Tabs, Without *ComfortWear* Panels

HPCPS: A4389



Transparent Stock No	Pouch Size	Barrier Opening	Pouch Length
3608	1-1/4" (32 mm)	7/8" (22 mm)	12" (30 cm)
3603	1-1/2" (38 mm)	1-1/8" (29 mm)	12" (30 cm)
3609	1-3/4" (44 mm)	1-3/8" (35 mm)	12" (30 cm)
3604	2" (51 mm)	1-1/2" (38 mm)	12" (30 cm)
3605	2-1/2" (64 mm)	2" (51 mm)	12" (30 cm)
3606	3" (76 mm)	2-1/2" (64 mm)	12" (30 cm)



323x (PRE-SIZED)

Box of 30

- ▶ Drainable Mini-Pouch, *Karaya 5* (standard wear)
Skin Barrier — Convex, Tape Border, Without Integrated Filter, Pouch Clamp Closure, Belt Tabs, Without *ComfortWear* Panels

HPCPS: A4389



Transparent Stock No	Pouch Size	Barrier Opening	Pouch Length
3232	1" (25 mm)	5/8" (16 mm)	9" (23 cm)
3238	1-1/4" (32 mm)	7/8" (22 mm)	9" (23 cm)
3233	1-1/2" (38 mm)	1-1/8" (29 mm)	9" (23 cm)
3239	1-3/4" (44 mm)	1-3/8" (35 mm)	9" (23 cm)
3234	2" (51 mm)	1-1/2" (38 mm)	9" (23 cm)



322x (PRE-SIZED)

Box of 30

- ▶ Drainable Pouch, *Karaya 5* (standard wear)
Skin Barrier — Convex, Tape Border, Without Integrated Filter, Pouch Clamp Closure, Belt Tabs, Without *ComfortWear* Panels

HPCPS: A4389



Transparent Stock No	Pouch Size	Barrier Opening	Pouch Length
3222	1" (25 mm)	5/8" (16 mm)	12" (30 cm)
3228	1-1/4" (32 mm)	7/8" (22 mm)	12" (30 cm)
3223	1-1/2" (38 mm)	1-1/8" (29 mm)	12" (30 cm)
3229	1-3/4" (44 mm)	1-3/8" (35 mm)	12" (30 cm)
3224	2" (51 mm)	1-1/2" (38 mm)	12" (30 cm)
3225	2-1/2" (64 mm)	2" (51 mm)	12" (30 cm)
3226	3" (76 mm)	2-1/2" (64 mm)	12" (30 cm)

Karaya 5 One-Piece Pouching System

Drainable Pouches (continued)



327x (PRE-SIZED)

Box of 30

- ▶ Drainable Pouch, *Karaya 5* (standard wear)
Skin Barrier — Convex, Tape Border, Without Integrated Filter,
Pouch Clamp Closure, Belt Tabs, Without *ComfortWear* Panels

HCPSCS: A4389



Transparent Stock No	Pouch Size	Barrier Opening	Pouch Length
3278	1-1/4" (32 mm)	7/8" (22 mm)	16" (40 cm)
3273	1-1/2" (38 mm)	1-1/8" (29 mm)	16" (40 cm)
3279	1-3/4" (44 mm)	1-3/8" (35 mm)	16" (40 cm)
3274	2" (51 mm)	1-1/2" (38 mm)	16" (40 cm)
3275	2-1/2" (64 mm)	2" (51 mm)	16" (40 cm)
3276	3" (76 mm)	2-1/2" (64 mm)	16" (40 cm)

Closed Pouches



332x (PRE-SIZED)

Box of 30

- ▶ Closed Pouch, *Karaya 5* (standard wear)
Skin Barrier — Convex, Tape Border, Integrated Filter,
Belt Tabs, Without *ComfortWear* Panels

HCPSCS: A4417



Transparent Stock No	Pouch Size	Barrier Opening	Pouch Length
3328	1-1/4" (32 mm)	7/8" (22 mm)	9" (23 cm)
3323	1-1/2" (38 mm)	1-1/8" (29 mm)	9" (23 cm)
3329	1-3/4" (44 mm)	1-3/8" (35 mm)	9" (23 cm)
3324	2" (51 mm)	1-1/2" (38 mm)	9" (23 cm)
3325	2-1/2" (64 mm)	2" (51 mm)	9" (23 cm)
3326	3" (76 mm)	2-1/2" (64 mm)	9" (23 cm)



716x (PRE-SIZED)

Box of 30

- ▶ Closed Pouch, *Karaya 5* (standard wear)
Skin Barrier — Convex, Without Tape Border, Without
Integrated Filter, Belt Tabs, Without *ComfortWear* Panels
- ▶ Requires a Belt

HCPSCS: A4387



Transparent Stock No	Pouch Size	Barrier Opening	Pouch Length
7168	1-1/4" (32 mm)	7/8" (22 mm)	9" (23 cm)
7163	1-1/2" (38 mm)	1-1/8" (29 mm)	9" (23 cm)
7169	1-3/4" (44 mm)	1-3/8" (35 mm)	9" (23 cm)
7164	2" (51 mm)	1-1/2" (38 mm)	9" (23 cm)
7165	2-1/2" (64 mm)	2" (51 mm)	9" (23 cm)
7166	3" (76 mm)	2-1/2" (64 mm)	9" (23 cm)

Urostomy Pouches



143x (PRE-SIZED)

Box of 10

- ▶ Urostomy Pouch, *Karaya 5* (standard wear)
Skin Barrier — Convex, Tape Border, Belt Tabs,
Without *ComfortWear* Panels

HCPSCS: A4429



Transparent Stock No	Pouch Size	Barrier Opening	Pouch Length
1432	1" (25 mm)	5/8" (16 mm)	9" (23 cm)
1438	1-1/4" (32 mm)	7/8" (22 mm)	9" (23 cm)
1433	1-1/2" (38 mm)	1-1/8" (29 mm)	9" (23 cm)
1439	1-3/4" (44 mm)	1-3/8" (35 mm)	9" (23 cm)
1434	2" (51 mm)	1-1/2" (38 mm)	9" (23 cm)

Ostomy Care
Ceramide-Infused
Products

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Ostomy Accessories

Continence
Care

Bowel Care

Tube
Fasteners

Wound Care

Education

Continence Care

Ostomy Care
Ceramide Infused
Products

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Ostomy Accessories

Continence
Care

Bowel Care


Tube
Fasteners

Wound Care

Education

Hollister Secure Start Services

Personalized Patient Support



Support Throughout the Continuum of Care

We are proud to offer dedicated support for each and every patient we serve along the continuum of care. Hollister Secure Start services provide a lifetime of personalized support for our members.

Our team is focused on helping patients with:

- Finding the right product fit
- Navigating insurance coverage
- Identifying supplier options
- Providing condition specific education
- And more!

To learn more, call us at **1.888.808.7456**.



The Hollister logo and Secure Start are trademarks of Hollister Incorporated.
©2017 Hollister Incorporated.

Medicare Utilization Chart

Intermittent Catheters

VaPro Plus Touch Free Hydrophilic Intermittent Catheter

VaPro Plus Pocket Touch Free Hydrophilic Intermittent Catheter

VaPro Touch Free Hydrophilic Intermittent Catheter, Coudé Tip

Onli Ready to Use Hydrophilic Intermittent Catheter

Advance Plus Touch Free Intermittent Catheter System

Advance Plus Pocket Touch Free Intermittent Catheter System

Advance Plus Touch Free Intermittent Catheter Kit

Apogee HC Hydrophilic Intermittent Catheter

Apogee Intermittent Catheters

Apogee Plus System

Apogee Plus System Kit

Male External Catheters

InView Male External Catheters

Extended Wear Male External Catheter

Everyday Male External Catheter

Leg Bag Systems

Vented Urinary Leg Bag Combination Pack

Urinary Leg Bag

Urinary Leg Bag Combination Pack

Leg Bag Accessories

Fabric Leg Bag Straps

Extension Tubing

Collection Systems

Bedside Drainage Collection System with Anti-Reflux Valve

Male Urinary Pouch External Collection Device

Female Urinary Pouch External Collection Device

Ostomy Care
Ceramide Infused
Products

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Ostomy Accessories

Continence
Care

Bowel Care

Tube
Fasteners

Wound Care

Education

Ostomy Care
Ceramide Infused Products

Ostomy Care
New Image Two-Piece Pouching System

Ostomy Care
Premier One-Piece Pouching System

Ostomy Care
Pouchkins Pediatric Pouching System

Ostomy Care
Ostomy Accessories

Continence Care

Bowel Care

Tube Fasteners

Wound Care


Education



VaPro Plus Pocket

Touch Free Hydrophilic Intermittent Catheter

We're for the one who makes their dreams a reality.



We want to help you get there.

The VaPro Plus Pocket hydrophilic intermittent catheter provides the right balance of ease of use and protection without compromising convenience or discretion. Its ready to use, touch free catheter with integrated collection bag is designed for easy handling and protection against bacterial contamination.

Request a sample: call 888.740.8999 or visit www.hollister.com/vapro

RX Only. Federal (USA) Law restricts this device to sale by or on the order of a physician or other healthcare practitioner licensed under state law to order this product.



Continence Care
People First

Medicare's Usual Maximum Quantity for Urological Supplies

The quantity of urological supplies needed by a patient is determined by many factors. There will be variation according to individual need and their needs may vary over time. The table below lists the maximum number of items/units of service that are medically necessary. The actual quantity needed for a particular patient may be more or less than the amount listed. The medical necessity for use of a greater quantity of supplies than the amounts listed must be clearly documented in the patient's medical record, and must be available upon request. The HCPCS and reimbursement information provided herein is intended to provide general information concerning coding of Hollister products only. Please refer to the **HCPCS Code Disclaimer** for more information.

Per Medicare LCD for Urological Supplies (L33803), effective 1-1-2017.

HCPCS	Long Description	Quantity
A4217	Sterile water/saline, 500 mL	*
A4310	Insert tray without drainage bag and without catheter (accessories only)	1 tray/insertion
A4311	Insertion tray without drainage bag with indwelling catheter, foley type, two-way latex with coating (Teflon, silicone, silicone elastomer, or hydrophilic, etc.)	1 tray/insertion
A4312	Insertion tray without drainage bag with indwelling catheter, foley type, two-way all silicone	1 tray/insertion
A4313	Insertion tray without drainage bag with indwelling catheter, foley type, three-way, for continuous irrigation	1 tray/insertion
A4314	Insertion tray with drainage bag with indwelling catheter, foley type, two-way latex with coating (Teflon, silicone, silicone elastomer, or hydrophilic, etc.)	1 ea/mo
A4315	Insertion tray with drainage bag with indwelling catheter, foley type, two-way, all silicone	1 ea/mo
A4316	Insertion tray with drainage bag with indwelling catheter, foley type, three-way, for continuous irrigation	1 ea/mo
A4320	Irrigation tray with bulb or piston syringe, any purpose	*
A4321	Therapeutic agent for urinary catheter irrigation	*
A4322	Irrigation syringe, bulb, or piston, each	*
A4326	Male external catheter with integral collection chamber, any type, each	*
A4327	Female external urinary collection device; meatal cup, each	1 ea/wk
A4328	Female external urinary collection device; pouch, each	1 ea/day
A4331	Extension drainage tubing, any type, any length, with connector/adaptor, for use with urinary leg bag or urostomy pouch, each	*
A4332	Lubricant, individual sterile packet, each	200 ea/mo
A4333	Urinary catheter anchor device, adhesive skin attachment, each	3 ea/wk
A4334	Urinary catheter anchoring device, leg strap, each	1 ea/mo
A4335	Incontinence supply; miscellaneous	*
A4336	Incontinence supply, urethral insert, any type, each	*
A4338	Indwelling catheter; foley type, two-way latex with coating (Teflon, silicone, silicone elastomer, or hydrophilic, etc.), each	1 ea/mo
A4340	Indwelling catheter; specialty type, (e.g., coudé, mushroom, wing, etc.), each	1 ea/mo
A4344	Indwelling catheter, foley type, two-way, all silicone, each	1 ea/mo
A4346	Indwelling catheter; foley type, three-way for continuous irrigation, each	1 ea/mo
A4349	Male external catheter, with or without adhesive, disposable, each	35 ea/mo

* No established usual maximum, refer to Medicare Utilization Guidelines.

Medicare Utilization Chart

HCPCS	Long Description	Quantity
A4351	Intermittent urinary catheter; straight tip, with or without coating (Teflon, silicone, silicone elastomer, or hydrophilic, etc.), each	200 ea/mo
A4352	Intermittent urinary catheter; coudé (curved) tip, with or without coating (Teflon, silicone, silicone elastomeric, or hydrophilic, etc.), each	200 ea/mo
A4353	Intermittent urinary catheter, with insertion supplies	200 ea/mo
A4354	Insertion tray with drainage bag but without catheter	1 ea/mo
A4355	Irrigation tubing set for continuous bladder irrigation through a three-way indwelling foley catheter, each	1 ea/day
A4356	External urethral clamp or compression device (not to be used for catheter clamp), each	1 ea/3 mo
A4357	Bedside drainage bag, day or night, with or without anti-reflux device, with or without tube, each	2 ea/mo
A4358	Urinary drainage bag, leg or abdomen, vinyl, with or without tube, with straps, each	2 ea/mo
A4360	Disposable external urethral clamp or compression device, with pad and/or pouch, each	Not Covered
A4402	Lubricant, per ounce	*
A4450	Tape, non-waterproof, per 18 square inches	10 units (18 sq in)/mo
A4452	Tape, waterproof, per 18 square inches	10 units (18 sq in)/mo
A4455	Adhesive remover or solvent (for tape, cement, or other adhesive), per ounce	*
A4520	Incontinence garment, any type (e.g., brief, diaper), each	*
A4456	Adhesive remover, wipes, any type, each	Not Covered
A4554	Disposable underpads, all sizes	*
A5102	Bedside drainage bottle with or without tubing, rigid or expandable, each	1 ea/3 mo
A5105	Urinary suspensory with leg bag, with or without tube, each	*
A5112	Urinary drainage bag, leg or abdomen, latex, with or without tube, with straps, each	1 ea/mo
A5113	Leg strap; latex, replacement only, per set	*
A5114	Leg strap, foam or fabric, replacement only, per set	*
A5131	Appliance cleaner, incontinence and ostomy appliances, per 16 ounce	1 (16 oz)/mo
A5200	Percutaneous catheter/tube anchoring device, adhesive skin attachment	*
A9270	Non-covered item or service	

* No established usual maximum, refer to Medicare Utilization Guidelines.

Note: Intermittent catheterization using a sterile intermittent catheter kit (A4353)

This product is covered when the patient requires catheterization and the beneficiary meets one of the following criteria (1-5):

- 1) The beneficiary resides in a nursing facility
- 2) The beneficiary is immunosuppressed, for example (not all-inclusive): on a regimen of immunosuppressive drugs post-transplant, on cancer chemotherapy, has AIDS, has a drug-induced state such as chronic oral corticosteroid use
- 3) The beneficiary has radiologically documented vesico-urethral reflux while on a program of intermittent catheterization
- 4) The beneficiary is a spinal-cord-injured-female with neurogenic bladder who is pregnant (for duration of pregnancy only)
- 5) The beneficiary has had distinct, recurrent urinary tract infections, while on a program of sterile intermittent catheterization with A4351/A4352 and sterile lubricant A4332, twice within the 12 months prior to the initiation of sterile intermittent catheter kits.

Per Medicare LCD for Urological Supplies (L33803), effective 1-1-2017.

Ostomy Care
Ceramide Infused Products

Ostomy Care
New Image Two-Piece Pouching System

Ostomy Care
Premier One-Piece Pouching System

Ostomy Care
Pouchkins Pediatric Pouching System

Ostomy Care
Ostomy Accessories

Continence Care

Bowel Care

Tube Fasteners

Wound Care

Education

Onli

Ready to Use Hydrophilic Intermittent Catheter

Whether their day is ordinary or extraordinary, the Onli catheter is there for them



The Onli catheter fits into your patient's day — whatever the day may bring. This easy to open hydrophilic intermittent catheter is designed for smooth insertion with reduced spills and no staining.** This ready to use catheter is as easy to learn as it is to use, so you can be confident that he will be able to use the Onli catheter the right way at the right time.

For more information or a sample, call 888.740.8999 or visit hollister.com/onli.

Ready. Simple.

**On denim, linen, or polyester. Data on file.

RX Only. Federal (USA) Law restricts this device to sale by or on the order of a physician or other healthcare practitioner licensed under state law to order this product.

 **Hollister.**

Continence Care
People First

Ostomy Care
Ceramide Infused
Products

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Ostomy Accessories

Continence
Care

Bowel Care

Tube
Fasteners

Wound Care

Education

Intermittent Catheters

VaPro Plus Touch Free Hydrophilic Intermittent Catheter



74xxx

Box of 30

- ▶ Protective Tip and Sleeve allows for no touch catheterization
- ▶ Pre-Hydrated and Ready to Use
- ▶ Integrated 1000 mL Collection Bag
- ▶ Smooth hydrophilic coating
- ▶ Smooth Eyelets
- ▶ Not Made with Natural Rubber Latex

HCPCS: A4353

Stock No	Length (inches)	Size	Description
74122	8	12 Fr	Straight
74142	8	14 Fr	Straight
74124	16	12 Fr	Straight
74144	16	14 Fr	Straight

VaPro Plus Pocket Touch Free Hydrophilic Intermittent Catheter



71xxx-30

Box of 30

- ▶ All the attributes of *VaPro Plus* Catheter in a Pocket Size Package

HCPCS: A4353

Stock No	Length (inches)	Size	Description
71082-30	8	8 Fr	Straight
71102-30	8	10 Fr	Straight
71122-30	8	12 Fr	Straight
71142-30	8	14 Fr	Straight
71084-30	16	8 Fr	Straight
71104-30	16	10 Fr	Straight
71124-30	16	12 Fr	Straight
71144-30	16	14 Fr	Straight
71164-30	16	16 Fr	Straight

VaPro Touch Free Hydrophilic Intermittent Catheter, Coudé Tip



73xxx

Box of 30

- ▶ Protective Tip and Sleeve allows for no touch catheterization
- ▶ Pre-Hydrated and Ready to Use
- ▶ Coudé Tip
- ▶ Smooth Hydrophilic Coating
- ▶ Smooth Eyelets
- ▶ Not Made with Natural Rubber Latex

HCPCS: A4352

Stock No	Length (inches)	Size	Description
73124	16	12 Fr	Coudé
73144	16	14 Fr	Coudé
73164	16	16 Fr	Coudé

Onli Ready to Use Hydrophilic Intermittent Catheter



82xxx-30

Box of 30

- ▶ Pre-Hydrated and Ready to Use
- ▶ Smooth Hydrophilic Coating
- ▶ Smooth Eyelets
- ▶ Ergonomic Gripper on 16" Catheter
- ▶ Not Made with Natural Rubber Latex

HCPCS: A4351

Stock No	Length (inches)	Size	Description
82081-30	7	8 Fr	Straight
82101-30	7	10 Fr	Straight
82121-30	7	12 Fr	Straight
82141-30	7	14 Fr	Straight
82084-30	16	8 Fr	Straight
82104-30	16	10 Fr	Straight
82124-30	16	12 Fr	Straight
82144-30	16	14 Fr	Straight
82164-30	16	16 Fr	Straight

RX Only. Federal (USA) Law restricts this device to sale by or on the order of a physician or other healthcare practitioner licensed under state law to order this product.

Intermittent Catheters

Advance Plus Touch Free Intermittent Catheter System



94xxx/95xxx

Box of 100

- ▶ 16" length
- ▶ Gel Reservoir
- ▶ Protective Tip
- ▶ Integrated 1500 mL Collection Bag with Large Handle
- ▶ Smooth Eyelets
- ▶ Not Made with Natural Rubber Latex

HCPCS: A4353

Stock No	Size	Description
94064	6 Fr	Straight
94084	8 Fr	Straight
94104	10 Fr	Straight
94124	12 Fr	Straight
94144	14 Fr	Straight
94164	16 Fr	Straight
94184	18 Fr	Straight
95124	12 Fr	Coudé
95144	14 Fr	Coudé
95164	16 Fr	Coudé

Advance Plus Pocket Touch Free Intermittent Catheter System



93xxx

Box of 100

- ▶ All the attributes of the full length **Advance Plus** Catheter in a Pocket Size Package

HCPCS: A4353

Stock No	Size	Description
93104	10 Fr	Pocket
93124	12 Fr	Pocket
93144	14 Fr	Pocket
93164	16 Fr	Pocket

Advance Plus Touch Free Intermittent Catheter Kit



96xxx/97xxx

Box of 100

- ▶ **Advance Plus** Catheter
- ▶ Kit Contains Gloves, Underpad Drape, and Benzalkonium Chloride Swabsticks
- ▶ Not Made with Natural Rubber Latex

HCPCS: A4353

Stock No	Size	Description
96064	6 Fr	Straight
96084	8 Fr	Straight
96104	10 Fr	Straight
96124	12 Fr	Straight
96144	14 Fr	Straight
96164	16 Fr	Straight
96184	18 Fr	Straight
97124	12 Fr	Coudé
97144	14 Fr	Coudé
97164	16 Fr	Coudé

RX Only. Federal (USA) Law restricts this device to sale by or on the order of a physician or other healthcare practitioner licensed under state law to order this product.

Intermittent Catheters

Apogee HC Hydrophilic Intermittent Catheter



90xxx/91xxx

Box of 30

- ▶ Smooth Insertion
- ▶ Simple To Use
- ▶ Not Made with Natural Rubber Latex
- ▶ Sliding sleeve to aid insertion for 16 inch product
- ▶ Burstable sterile water pack that activates coating

HCPCS: A4351 for all except

HCPCS: A4352 for 91226, 91426, 91626

Stock No	Length	Size	Description
91006	6	10 Fr	Straight
91206	6	12 Fr	Straight
91406	6	14 Fr	Straight
90810	10	8 Fr	Straight
91010	10	10 Fr	Straight
91216	16	12 Fr	Straight
91416	16	14 Fr	Straight
91616	16	16 Fr	Straight
91816	16	18 Fr	Straight
91226	16	12 Fr	Coudé
91426	16	14 Fr	Coudé
91626	16	16 Fr	Coudé

RX Only. Federal (USA) Law restricts this device to sale by or on the order of a physician or other healthcare practitioner licensed under state law to order this product.

Intermittent Catheters

Apogee Family of Intermittent Catheters



10xx/1xxxx

- ▶ Smooth Eyelets
- ▶ Not Made with Natural Rubber Latex (unless noted)
- ▶ Non-lubricated

Box of 30

Box of 50

Apogee IC Intermittent Catheters

HCPCS: A4351

Stock No	Length (inches)	Size	Description
10806	6	8 Fr	Firm
11006	6	10 Fr	Firm
11206	6	12 Fr	Firm
11406	6	14 Fr	Firm
11446*	6	14 Fr	Firm
10610	10	6 Fr	Firm
10810	10	8 Fr	Firm
11010	10	10 Fr	Firm
10616	16	6 Fr	Firm
10816	16	8 Fr	Firm
11016	16	10 Fr	Firm
11216	16	12 Fr	Firm
11416	16	14 Fr	Firm
11616	16	16 Fr	Firm
11816	16	18 Fr	Firm

*No funnel Qty: 30/box

Apogee IC Intermittent Catheters Coudé Tip

HCPCS: A4352

Stock No	Length (inches)	Size	Description
10826	16	8 Fr	Firm
11026	16	10 Fr	Firm
11226	16	12 Fr	Firm
11426	16	14 Fr	Firm
11626	16	16 Fr	Firm
11826	16	18 Fr	Firm

Qty: 30/box

Apogee Intermittent Catheters, Other Styles

HCPCS: A4351 for all except

HCPCS: A4352 for 1087

Stock No	Length (inches)	Size	Description
1036	16	10 Fr	Soft
1046	16	12 Fr	Soft
1061	16	14 Fr	Soft
1083	16	16 Fr	Soft
1049	16	12 Fr	Firm, Curved Catheter
1065	16	14 Fr	Firm, Curved Catheter
1074A	16	14 Fr	Soft, Curved Catheter
1087	16	16 Fr	Firm, Coudé, Curved Catheter
1051	16	12 Fr	Red Rubber Latex*
1067	16	14 Fr	Red Rubber Latex*
1084	16	16 Fr	Red Rubber Latex*

Qty: 30/box except curved catheters: 50/box

***Caution:** This product contains natural rubber latex which may cause allergic reactions.

RX Only. Federal (USA) Law restricts this device to sale by or on the order of a physician or other healthcare practitioner licensed under state law to order this product.

Intermittent Catheters

Apogee Plus System



Bxxx

Case of 100

- ▶ 16" length
- ▶ Protective Tip
- ▶ Pre-lubricated
- ▶ Smooth Eyelets
- ▶ Integrated 1500 mL Collection Bag
- ▶ Not Made with Natural Rubber Latex (unless noted)

HCPCS: A4353

Stock No	Size	Description
B8F	8 Fr	Firm
B10F	10 Fr	Firm
B12F	12 Fr	Firm
B14F	14 Fr	Firm
B16F	16 Fr	Firm
B18F	18 Fr	Firm
B12C	12 Fr	Coudé
B14C	14 Fr	Coudé
B16C	16 Fr	Coudé
B12S	12 Fr	Soft
B14S	14 Fr	Soft
B14R	14 Fr	Red Rubber Latex*
B16R	16 Fr	Red Rubber Latex*
B14F-Female	14 Fr, 6"	Firm

***Caution:** This product contains natural rubber latex which may cause allergic reactions.

Apogee Plus System Kit



BxxxB

Case of 100

- ▶ **Apogee Plus Catheter**
- ▶ Kit Contains Gloves, Underpad Drape, Gauze, and PVP Swabsticks
- ▶ Not Made with Natural Rubber Latex (unless noted)

HCPCS: A4353

Stock No	Size	Description
B6FB	6 Fr	Firm**
B8FB	8 Fr	Firm**
B10FB	10 Fr	Firm
B12FB	12 Fr	Firm
B14FB	14 Fr	Firm
B16FB	16 Fr	Firm
B18FB	18 Fr	Firm
B12CB	12 Fr	Coudé
B14CB	14 Fr	Coudé
B16CB	16 Fr	Coudé
B12SB	12 Fr	Soft
B14SB	14 Fr	Soft
B12RB	12 Fr	Red Rubber Latex*
B14RB	14 Fr	Red Rubber Latex*
B16RB	16 Fr	Red Rubber Latex*
B14FB-Female	14 Fr, 6"	Firm

⚠ WARNING: The gloves in this kit can expose you to chemicals including Diisononyl phthalate (DINP) which is known to the State of California to cause cancer. For more information go to www.P65Warnings.ca.gov.

***Caution:** This product contains natural rubber latex which may cause allergic reactions.

******The protective tip on this product is closer to 10mm.

RX Only. Federal (USA) Law restricts this device to sale by or on the order of a physician or other healthcare practitioner licensed under state law to order this product.

Ostomy Care
Ceramide Infused Products

Ostomy Care
New Image Two-Piece Pouching System

Ostomy Care
Premier One-Piece Pouching System

Ostomy Care
Pouchkins Pediatric Pouching System

Ostomy Care
Ostomy Accessories

Continence Care

Bowel Care

Tube Fasteners

Wound Care

Education

Male External Catheters

InView Male External Catheters



974xx/975xx/976xx
974xx-100/975xx-100/976xx-100

Box of 30
Box of 100

- ▶ Self-Adhesive
- ▶ Clear Material Enables Monitoring of Skin
- ▶ Not Made with Natural Rubber Latex

HCPCS: A4349

SPECIAL: shorter length sheath



Stock No Box of 30	Stock No Box of 100	Color	Size
97425	97425-100	Green	25 mm
97429	97429-100	Purple	29 mm
97432	97432-100	Brown	32 mm
97436	97436-100	Blue	36 mm
97441	97441-100	Orange	41 mm

STANDARD: the daily catheter for general use



Stock No Box of 30	Stock No Box of 100	Color	Size
97525	97525-100	Green	25 mm
97529	97529-100	Purple	29 mm
97532	97532-100	Brown	32 mm
97536	97536-100	Blue	36 mm
97541	97541-100	Orange	41 mm

EXTRA: 70% more adhesive than Standard for added security



Stock No Box of 30	Stock No Box of 100	Color	Size
97625	97625-100	Green	25 mm
97629	97629-100	Purple	29 mm
97632	97632-100	Brown	32 mm
97636	97636-100	Blue	36 mm
97641	97641-100	Orange	41 mm

Extended Wear Male External Catheter



- 920x**
- ▶ Self-Adhesive
 - ▶ Unique Inner Flap
 - ▶ Easy Applicator Collar
 - ▶ Contains Natural Rubber Latex

Box of 30

HCPCS: A4349

Stock No	Color	Size
9206	Green	22 mm-25 mm
9207	Purple	26 mm-30 mm
9208	Brown	31 mm-35 mm
9209	Blue	36 mm-39 mm

Caution: This product contains natural rubber latex which may cause allergic reactions.

Everyday Male External Catheter



- 910x**
- ▶ Self-Adhesive
 - ▶ Easy Applicator Collar
 - ▶ Contains Natural Rubber Latex

Box of 30

HCPCS: A4349

Stock No	Color	Size
9106	Green	22 mm-25 mm
9107	Purple	26 mm-30 mm
9108	Brown	31 mm-35 mm
9109	Blue	36 mm-39 mm

Caution: This product contains natural rubber latex which may cause allergic reactions.

Leg Bag Systems

Vented Urinary Leg Bag Combination Pack



9645/9655

Box of 10

- ▶ Includes:
 - (1) Vented Urinary Leg Bag — Medium or Large
 - (1) Extension Tubing with Vacuum Relief Valve, (18" Long, Oval, 11/32" Effective Inside Diameter), with Connector
 - (1) Pair, Fabric Leg Bag Straps — Medium (23" Long, 1-1/8" Wide)
- ▶ Not Made with Natural Rubber Latex, Contains Sterile Components
- ▶ Anti-reflux valve to help prevent back flow of urine into tube

HCPCS: A4358

Stock No	Description	Length	Width
9645	Medium — 18 oz (540 mL)	10"	4-1/4"
9655	Large — 30 oz (900 mL)	11-1/8"	5"

Urinary Leg Bag



9814/9805

Box of 10

- ▶ Includes:
 - (1) Urinary Leg Bag — Medium or Large
- ▶ Not Made with Natural Rubber Latex
- ▶ Anti-reflux valve to help prevent back flow of urine into tube

HCPCS: A9270

HCPCS: A4358 (if straps are added)

Stock No	Description	Length	Width
9814	Medium — 18 oz (540 mL)	10"	4-1/4"
9805	Large — 30 oz (900 mL)	11-1/8"	5"

Urinary Leg Bag Combination Pack



9348/9349

Box of 12

- ▶ Includes:
 - (1) Urinary Leg Bag — Medium or Large
 - (1) Extension Tubing (18" Long, Oval, 11/32" Effective Inside Diameter) with Connector
 - (1) Pair, Fabric Leg Bag Straps — Medium (23" Long, 1-1/8" Wide)
- ▶ Not Made with Natural Rubber Latex; Contains Sterile Components
- ▶ Anti-reflux valve to help prevent back flow of urine into tube

HCPCS: A4358

Stock No	Description	Length	Width
9348	Medium — 18 oz (540 mL)	10"	4-1/4"
9349	Large — 30 oz (900 mL)	11-1/8"	5"

Ostomy Care
Ceramide Infused
Products

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Ostomy Accessories

Continence
Care

Bowel Care

Tube
Fasteners

Wound Care

Education

Leg Bag Accessories

Fabric Leg Bag Straps



9342/9343

10 Pair

- ▶ Not Made with Natural Rubber Latex

HCPCS: A5114

Stock No	Description	Length	Width
9342	Small	15"	1-1/8"
9343	Medium	23"	1-1/8"

Extension Tubing



9345/9346

Box of 10

- ▶ Kink-resistant with Connector
- ▶ Oval, 11/32" Effective Inside Diameter
- ▶ Thermoplastic Elastomer, Not Made with Natural Rubber Latex

HCPCS: A4331

Stock No	Description	Length
9345	Non Sterile	18"
9346	Sterile	18"

Ostomy Care
Ceramide Infused
Products

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Ostomy Accessories

Continence
Care

Bowel Care

Tube
Fasteners

Wound Care

Education

Collection Systems

Bedside Drainage Collection System with Anti-Reflux Valve



9839

Box of 20

- ▶ Sample Port
- ▶ Sterile Fluid Pathway
- ▶ Not Made with Natural Rubber Latex

HCPCS: A4357

Stock No	Description	Length	Width
9839	68 oz (2000 mL)	12"	9-3/4"

Male Urinary Pouch External Collection Device*



9873

Box of 10

- ▶ Cut-to-Fit, *Flexend* Skin Barrier
- ▶ Thin and Flexible
- ▶ No Additional Tape Needed
- ▶ Odor-Barrier Film
- ▶ Convenience Drain
- ▶ Not Made with Natural Rubber Latex

HCPCS: A4326

Stock No	Description	Length	Width
9873	Male Urinary Pouch	7-1/2"	4-1/2"



9811

Box of 10

- ▶ Cut-to-Fit, *SoftFlex* Skin Barrier
- ▶ Odor-Barrier Film
- ▶ Convenience Drain
- ▶ Tape Strips and Skin Gel Wipe Included
- ▶ Not Made with Natural Rubber Latex

HCPCS: A4326

Stock No	Description	Length	Width
9811	Male Urinary Pouch	7-1/2"	4-1/2"

Female Urinary Pouch External Collection Device*



9840

Box of 10

- ▶ Cut-to-Fit, *SoftFlex* Skin Barrier
- ▶ Odor-Barrier Film
- ▶ Convenience Drain
- ▶ Alternative to Indwelling Catheter for Bedridden Females
- ▶ *Adapt* Paste, Tape Strips, and Skin Gel Wipe Included
- ▶ Not Made with Natural Rubber Latex

HCPCS: A4328

Stock No	Description	Length	Width
9840	Female Urinary Pouch	7-1/2"	4-1/2"

* Meant to be connected to leg bag or gravity drainage

Ostomy Care
Ceramide-Infused
Products

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Ostomy Accessories

Continence
Care

Bowel Care

Tube
Fasteners

Wound Care

Education

Bowel Care

Education

Wound Care

Tube
Fasteners

Bowel Care

Continence
Care

Ostomy Care
Ostomy Accessories

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Ceramide Infused
Products

Table of Contents • Bowel Care

Fecal Collectors

Fecal Collectors

Ostomy Care
Ceramide Infused
Products

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Ostomy Accessories

Continence
Care

Bowel Care

Tube
Fasteners

Wound Care

Education

Fecal Collectors

Fecal Collectors



9880

Box of 10

- ▶ Cut-to-Fit *Flectend* (extended wear) Tapered Skin Barrier
- ▶ Odor-Barrier Pouch Film
- ▶ Convenience Drain

HCPCS: A4330

Stock No	Description	Pouch Length
9880	Large	12" (30 cm)



9823

Box of 10

- ▶ Cut-to-Fit *Flectend* (extended wear) Tapered Skin Barrier
- ▶ Ultra-Clear Pouch Film
- ▶ Convenience Drain

HCPCS: A4330

Stock No	Description	Pouch Length
9823	Medium	10" (25 cm)



9822/9821

Box of 10

- ▶ Cut-to-Fit *SoftFlex* (standard wear) Skin Barrier
- ▶ Odor-Barrier Pouch Film
- ▶ Convenience Drain

HCPCS: A4330

Stock No	Description	Pouch Length
9822	Medium	10" (25 cm)
9821	Large	12" (30 cm)



DID YOU KNOW...

To view video instructions for product application, visit

<http://www.hollister.com/us/bm/fc/>

Ostomy Care
Ceramide-Infused
Products

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Ostomy Accessories

Continence
Care

Bowel Care

Tube
Fasteners

Wound Care

Education

Tube Fasteners

Education

Wound Care

Tube
Fasteners

Bowel Care

Continence
Care

Ostomy Care
Ostomy Accessories

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Ceramide Infused
Products

Table of Contents • Tube Fasteners

Tube Fasteners

- AnchorFast Guard** Oral Endotracheal Tube Fastener
- AnchorFast** Oral Endotracheal Tube Fastener
- Horizontal Tube Attachment Device (HTAD): Sterile
- Vertical Tube Attachment Device (VTAD): Sterile
- Drain Tube Attachment Device (DTAD): Sterile
- Feeding Tube Attachment Device (FTAD)
- Suction Tube Attachment Device (STAD)

Ostomy Care
Caramide Infused
Products

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Ostomy Accessories

Continence
Care

Bowel Care

Tube
Fasteners

Wound Care

Education

Ostomy Care
Ceramide Infused
Products

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Ostomy Accessories


Continance
Care

Bowel Care

Tube
Fasteners

Wound Care

Education



*When your choice
is too important
to take a chance.*

AnchorFast

Oral Endotracheal Tube Fastener

**With a long history of quality and innovation,
you can trust the AnchorFast product.**

- Used at over 3,700 hospitals in the United States
- 20 years of demonstrated quality
- Ongoing innovation driven by customer feedback

**Learn more at www.hollister.com/anchorfast
or call 888.740.8999**



CAUTION: Federal (USA) Law restricts this device to sale by or on the order of a physician. Prior to use of the AnchorFast oral endotracheal tube fastener, be sure to read the entire product Instructions for Use package insert that accompanies the product.

The Hollister logo and AnchorFast are trademarks of Hollister Incorporated.

© 2016 Hollister Incorporated

AnchorFast Guard Oral Endotracheal Tube Fastener



9800

- ▶ Integrated Tube Protection Sleeve
- ▶ Secures Oral Endotracheal Tubes (5 mm-10 mm inner diameter)
- ▶ Helps Optimize Oral Care
- ▶ Helps Safeguard Skin Health
- ▶ Individually Packaged
- ▶ Not Made with Natural Rubber Latex

Stock No	HCPCS	Description	Quantity/Box
9800	A9900* / A9999**	AnchorFast Guard Oral Endotracheal Tube Fastener	12

**to be used with a rented ventilator **to be used with a patient-owned ventilator*

AnchorFast Oral Endotracheal Tube Fastener



9799

- ▶ Secures Oral Endotracheal Tubes (5 mm-10 mm inner diameter)
- ▶ Helps Optimize Oral Care
- ▶ Helps Safeguard Skin Health
- ▶ Individually Packaged
- ▶ Not Made with Natural Rubber Latex

Stock No	HCPCS	Description	Quantity/Box
9799	A9900* / A9999**	AnchorFast Oral Endotracheal Tube Fastener	12

**to be used with a rented ventilator **to be used with a patient-owned ventilator*

AnchorFast Replacement Neck Band Straps

Stock No	HCPCS	Description	Quantity/Box
9788	—	Replacement Neck Band Straps	4

Caution: Federal (USA) Law restricts this device to sale by or on the order of a physician. Prior to use of the AnchorFast Guard oral endotracheal tube fastener, be sure to read the entire product Instructions for Use package insert that accompanies the product.



DID YOU KNOW...

For more information about **AnchorFast** oral endotracheal tube fasteners, visit www.hollister.com/anchorfast

Tube Fasteners

Horizontal Tube Attachment Device (HTAD): Sterile



9781

- ▶ Secures Drains, Tubes, or Catheters (5 Fr-40 Fr) Parallel to Patient's Skin
- ▶ Tubes Can be Secured at Exit Site, or Distal to Exit Site
- ▶ Helps Safeguard Skin Health
- ▶ Sterile, Not Made with Natural Rubber Latex

Stock No	HCPCS	Description	Quantity/Box
9781	A5200	Horizontal Tube Attachment Device (HTAD): Sterile	5

Vertical Tube Attachment Device (VTAD): Sterile



9782

- ▶ Secures Drains, Tubes, or Catheters (5 Fr-40 Fr) Vertical to Patient's Skin
- ▶ Secures Tube at Exit Site
- ▶ Helps Safeguard Skin Health
- ▶ Sterile, Not Made with Natural Rubber Latex

Stock No	HCPCS	Description	Quantity/Box
9782	A5200	Vertical Tube Attachment Device (VTAD): Sterile	5

Drain Tube Attachment Device (DTAD): Sterile



9780

- ▶ Secures a Variety of Drains or Tubes (10 Fr-24 Fr) Vertical to Patient's Skin
- ▶ Sterile, Not Made with Natural Rubber Latex

Stock No	HCPCS	Description	Quantity/Box
9780	A5200	Drain Tube Attachment Device (DTAD): Sterile	5

Feeding Tube Attachment Device (FTAD)



9786

- ▶ Secures Nasogastric/Nasointestinal Tubes (5 Fr-18 Fr)
- ▶ Butterfly-Shaped *FlexWear* Skin Barrier Conforms to Patient's Nose
- ▶ Unique Swivel Clamp Helps Eliminate Pressure on the Nares
- ▶ Not Made with Natural Rubber Latex

Stock No	HCPCS	Description	Quantity/Box
9786	B9998	Feeding Tube Attachment Device (FTAD)	20

Suction Tube Attachment Device (STAD)



9785

- ▶ Secures Nasogastric Tubes (12 Fr-18 Fr)
- ▶ Butterfly-Shaped *FlexWear* Skin Barrier Conforms to Patient's Nose
- ▶ Unique Swivel Clamp Helps Eliminate Pressure on the Nares
- ▶ Not Made with Natural Rubber Latex

Stock No	HCPCS	Description	Quantity/Box
9785	A9999	Suction Tube Attachment Device (STAD)	20

Ostomy Care
Ceramide-Infused
Products

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Ostomy Accessories

Continence
Care

Bowel Care

Tube
Fasteners

Wound Care

Education

Wound Care

Education

Wound Care

Tube
Fasteners

Bowel Care

Continence
Care

Ostomy Care
Ostomy Accessories

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Ceramide Infused
Products

Table of Contents • Wound Care

Medicare Utilization Chart

FastCast OLS Total Contact Cast System

Foam Dressings

TRIACT Foam with Silicone Border

Restore Foam with TRIACT Technology Non-Border

Restore LITE Foam with TRIACT Technology Non-Border

Restore Silver Foam with TRIACT Technology Non-Border

Restore Foam with Silicone Border

Restore Foam with Silicone Non-Border

Contact Layers with **TRIACT** Technology

Restore Contact Layer FLEX

Restore Contact Layer, Silver

Restore to **CalciCare** Conversion Chart

CalciCare and **Restore** Calcium Alginates

CalciCare Calcium Alginate

CalciCare Calcium Alginate, Silver

Restore Calcium Alginate

Restore Calcium Alginate, Silver

Restore Hydrogel

Restore Wound Cleanser and Skin Care

Restore Wound Cleanser

Restore Skin Cleanser

Restore Cleanser & Moisturizer

Restore Skin Conditioning Creme

Restore Hydrocolloids

Restore Hydrocolloid: Sterile

Restore Extra Thin Hydrocolloid: Sterile

Restore Hydrocolloid with Foam Backing: Sterile

Ostomy Care
Ceramide Infused
Products

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Ostomy Accessories

Continence
Care

Bowel Care

Tube
Fasteners

Wound Care

Education

Medicare Usual Dressing Change Frequency for Wound Care Supplies

The quantity and type of dressings dispensed at any one time must take into account the current status of the wound(s), the likelihood of change, and the recent use of dressings. Dressing needs may change frequently (e.g., weekly) in the early phases of wound treatment and/or with heavily draining wounds. Surgical dressings must be tailored to the specific needs of an individual patient. The table below lists some coverage guidelines and usual dressing change frequencies. The medical necessity for more frequent dressing changes must be documented in the patient's medical record and submitted with the claim. The HCPCS and reimbursement information provided herein is intended to provide general information concerning coding of Hollister products only. Please refer to the **HCPCS Code Disclaimer** for more information.

Per Medicare LCD for Surgical Dressings (L33831), effective 7-24-2017.

HCPCS	Long Description	Frequency of Change	Wound Condition
A4450	Tape, non-waterproof, per 18 square inches	based on drsg size	
A4452	Tape, waterproof, per 18 square inches	based on drsg size	
A4461	Surgical dressing holder, non-reusable, each	*	
A4463	Surgical dressing holder, reusable, each	*	
A4465	Non-elastic binder for extremity	*	
A4490	Surgical stockings above knee length, each	*	
A4495	Surgical stockings thigh length, each	*	
A4500	Surgical stockings below knee length, each	*	
A4510	Surgical stockings full length, each	*	
A4649	Surgical supply; miscellaneous	*	
A6025	Gel Sheet for dermal or epidermal application, (e.g., silicone, hydrogel, other), each	*	
A6154	Wound pouch, each	3/wk	
A6196	Alginate or other fiber gelling dressing, wound cover, sterile, pad size 16 sq in or less, each dressing	1/day	mod-heavy exudate
A6197	Alginate or other fiber gelling dressing, wound cover, sterile, pad size more than 16 sq in but less than or equal to 48 sq in, each dressing	1/day	mod-heavy exudate
A6198	Alginate or other fiber gelling dressing, wound cover, sterile, pad size more than 48 sq in, each dressing	1/day	mod-heavy exudate
A6199	Alginate or other fiber gelling dressing, wound filler, sterile, per 6 inches	2 units (12 in)/day	mod-heavy exudate

* No established usual maximum, refer to Medicare Utilization Guidelines.

Medicare Utilization Chart

HCPDS	Long Description	Frequency of Change	Wound Condition
A6203	Composite dressing, sterile, pad size 16 sq in or less, with any size adhesive border, each dressing	3/wk	
A6204	Composite dressing, sterile, pad size more than 16 sq in but less than or equal to 48 sq in, with any size adhesive border, each dressing	3/wk	
A6205	Composite dressing, sterile, pad size more than 48 sq in, with any size adhesive border, each dressing	3/wk	
A6206	Contact layer, sterile, 16 sq in or less, each dressing	1/wk	line the wound
A6207	Contact layer, sterile, more than 16 sq in but less than or equal to 48 sq in, each dressing	1/wk	line the wound
A6208	Contact layer, sterile, more than 48 sq in, each dressing	1/wk	line the wound
A6209	Foam dressing, wound cover, sterile, pad size 16 sq in or less, without adhesive border, each dressing	3/wk	mod-heavy exudate
A6210	Foam dressing, wound cover, sterile, pad size more than 16 sq in but less than or equal to 48 sq in, without adhesive border, each dressing	3/wk	mod-heavy exudate
A6211	Foam dressing, wound cover, sterile, pad size more than 48 sq in, without adhesive border, each dressing	3/wk	mod-heavy exudate
A6212	Foam dressing, wound cover, sterile, pad size 16 sq in or less, with any size adhesive border, each dressing	3/wk	mod-heavy exudate
A6213	Foam dressing, wound cover, sterile, pad size more than 16 sq in but less than or equal to 48 sq in, with any size adhesive border, each dressing	3/wk	mod-heavy exudate
A6214	Foam dressing, wound cover, sterile, pad size more than 48 sq in, with any size adhesive border, each dressing	3/wk	mod-heavy exudate
A6216	Gauze, non-impregnated, non sterile, pad size 16 sq in or less, without adhesive border, each dressing	3/day	
A6217	Gauze, non-impregnated, non sterile, pad size more than 16 sq in but less than or equal to 48 sq in, without adhesive border, each dressing	3/day	
A6218	Gauze, non-impregnated, non sterile, pad size more than 48 sq in, without adhesive border, each dressing	3/day	
A6219	Gauze, non-impregnated, sterile, pad size 16 sq in or less, with any size adhesive border, each dressing	1/day	

* No established usual maximum, refer to Medicare Utilization Guidelines.

Ostomy Care
Ceramide Infused
Products

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Ostomy Accessories

Continence
Care

Bowel Care

Tube
Fasteners

Wound Care

Education

Medicare Utilization Chart

HCPCS	Long Description	Frequency of Change	Wound Condition
A6220	Gauze, non-impregnated, sterile, pad size more than 16 sq in but less than or equal to 48 sq in, with any size adhesive border, each dressing	1/day	
A6221	Gauze, non-impregnated, sterile, pad size more than 48 sq in, with any size adhesive border, each dressing	1/day	
A6222	Gauze, impregnated with other than water, normal saline, or hydrogel, sterile pad size 16 sq in or less, without adhesive border, each dressing	1/day	
A6223	Gauze, impregnated with other than water, normal saline, or hydrogel, sterile pad size more than 16 sq in but less than or equal to 48 sq in, without adhesive border, each dressing	1/day	
A6224	Gauze, impregnated with other than water, normal saline, or hydrogel, sterile pad size more than 48 sq in, without adhesive border, each dressing	1/day	
A6228	Gauze, impregnated with water or normal saline, sterile, pad size 16 sq in or less, without adhesive border, each dressing		* Non-Covered
A6229	Gauze, impregnated with water or normal saline, sterile, pad size more than 16 sq in but less than or equal to 48 sq in, without adhesive border, each dressing		* Non-Covered
A6230	Gauze, impregnated with water or normal saline, sterile, pad size more than 48 sq in, without adhesive border, each dressing		* Non-Covered
A6234	Hydrocolloid dressing, wound cover, sterile, pad size 16 sq in or less, without adhesive border, each dressing	3/wk	light-mod exudate
A6235	Hydrocolloid dressing, wound cover, sterile, pad size more than 16 sq in but less than or equal to 48 sq in, without adhesive border, each dressing	3/wk	light-mod exudate
A6236	Hydrocolloid dressing, wound cover, sterile, pad size more than 48 sq in, without adhesive border, each dressing	3/wk	light-mod exudate
A6237	Hydrocolloid dressing, wound cover, sterile, pad size 16 sq in or less, with any size adhesive border, each dressing	3/wk	light-mod exudate
A6238	Hydrocolloid dressing, wound cover, sterile, pad size more than 16 sq in but less than or equal to 48 sq in, with any size adhesive border, each dressing	3/wk	light-mod exudate
A6239	Hydrocolloid dressing, wound cover, sterile, pad size more than 48 sq in, with any size adhesive border, each dressing	3/wk	light-mod exudate

* No established usual maximum, refer to Medicare Utilization Guidelines.

Medicare Utilization Chart

HPCS	Long Description	Frequency of Change	Wound Condition
A6240	Hydrocolloid dressing, wound filler, paste, sterile, per ounce	3/wk	light-mod exudate
A6241	Hydrocolloid dressing, wound filler, dry form, sterile, per gram	3/wk	light-mod exudate
A6242	Hydrogel dressing, wound cover, sterile, pad size 16 sq in or less, without adhesive border, each dressing	1/day	min-no exudate
A6243	Hydrogel dressing, wound cover, sterile, pad size more than 16 sq in but less than or equal to 48 sq in, without adhesive border, each dressing	1/day	min-no exudate
A6244	Hydrogel dressing, wound cover, sterile, pad size more than 48 sq in, without adhesive border, each dressing	1/day	min-no exudate
A6245	Hydrogel dressing, wound cover, sterile, pad size 16 sq in or less, with any size adhesive border, each dressing	3/wk	min-no exudate
A6246	Hydrogel dressing, wound cover, sterile, pad size more than 16 sq in but less than or equal to 48 sq in, with any size adhesive border, each dressing	3/wk	min-no exudate
A6247	Hydrogel dressing, wound cover, sterile, pad size more than 48 sq in, with any size adhesive border, each dressing	3/wk	min-no exudate
A6248	Hydrogel dressing, wound filler, gel, per fluid ounce	1/day up to 3 oz/mo	min-no exudate
A6250	Skin sealants, protectants, moisturizers, ointments, any type, any size	*	
A6251	Specialty absorptive dressing, wound cover, sterile, pad size 16 sq in or less without adhesive border, each dressing	1/day	mod-heavy exudate
A6252	Specialty absorptive dressing, wound cover, sterile, pad size more than 16 sq in but less than or equal to 48 sq in, without adhesive border, each dressing	1/day	mod-heavy exudate
A6253	Specialty absorptive dressing, wound cover, sterile, more than 48 sq in, without adhesive border, each dressing	1/day	mod-heavy exudate
A6254	Specialty absorptive dressing, wound cover, sterile, pad size 16 sq in or less with any size adhesive border, each dressing	every other day	mod-heavy exudate
A6255	Specialty absorptive dressing, wound cover, sterile, pad size more than 16 sq in but less than or equal to 48 sq in, with any size adhesive border, each dressing	every other day	mod-heavy exudate
A6256	Specialty absorptive dressing, wound cover, sterile, more than 48 sq in, with any size adhesive border, each dressing	every other day	mod-heavy exudate

* No established usual maximum, refer to Medicare Utilization Guidelines.

Medicare Utilization Chart

HCPCS	Long Description	Frequency of Change	Wound Condition
A6257	Transparent film, sterile, 16 sq in or less, each dressing	3/wk	min exudate
A6258	Transparent film, sterile, more than 16 sq in but less than or equal to 48 sq in, each dressing	3/wk	min exudate
A6259	Transparent film, sterile, more than 48 sq in, each dressing	3/wk	min exudate
A6260	Wound cleanser, any type, any size	*	
A6261	Wound filler, gel/paste, per fluid ounce, not otherwise classified	1/day	
A6262	Wound filler, dry form, per gram, not otherwise classified	1/day	
A6266	Gauze, impregnated, other than water, normal saline, or zinc paste, sterile, any width, per linear yard	1/day	
A6402	Gauze, non-impregnated, sterile, pad size 16 sq in or less, without adhesive border, each dressing	3/day	
A6403	Gauze, non-impregnated, sterile, pad size more than 16 sq in but less than or equal to 48 sq in, without adhesive border, each dressing	3/day	
A6404	Gauze, non-impregnated, sterile, pad size more than 48 sq in, without adhesive border, each dressing	3/day	
A6407	Packing strips, non-impregnated, sterile up to 2 inches in width, per linear yard	3/day	
A6410	Eye pad, sterile, each	*	
A6411	Eye pad, non sterile, each	*	
A6412	Eye patch, occlusive, each	*	
A6413	Adhesive bandage, first-aid type, any size each	*	
A6441	Padding bandage, non-elastic, non-woven/non-knitted, width greater than or equal to three inches and less than five inches, per yard	same as drsg chg freq	
A6442	Conforming bandage, non-elastic, knitted/woven, non sterile, width less than three inches, per yard	same as drsg chg freq	
A6443	Conforming bandage, non-elastic, knitted/woven, non sterile, width greater than or equal to three inches and less than five inches, per yard	same as drsg chg freq	
A6444	Conforming bandage, non-elastic, knitted/woven, non sterile, width greater than or equal to five inches, per yard	same as drsg chg freq	

* No established usual maximum, refer to Medicare Utilization Guidelines.

Ostomy Care
Ceramide Infused
Products

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Ostomy Accessories

Continence
Care

Bowel Care

Tube
Fasteners

Wound Care

Education

Medicare Utilization Chart

HCPCS	Long Description	Frequency of Change	Wound Condition
A6445	Conforming bandage, non-elastic, knitted/woven, sterile, width less than three inches, per yard	same as drsg chg freq	
A6446	Conforming bandage, non-elastic, knitted/woven, sterile, width greater than or equal to three inches and less than five inches, per yard	same as drsg chg freq	
A6447	Conforming bandage, non-elastic, knitted/woven, sterile, width greater than or equal to five inches, per yard	same as drsg chg freq	
A6448	Light compression bandage, elastic, knitted/woven, width less than three inches, per yard	1/wk	
A6449	Light compression bandage, elastic, knitted/woven, width greater than or equal to three inches and less than five inches, per yard	1/wk	
A6450	Light compression bandage, elastic, knitted/woven, width greater than or equal to five inches, per yard	1/wk	
A6451	Moderate-compression bandage, elastic, knitted/woven, load resistance of 1.25 to 1.34 foot pounds at 50% maximum stretch, width greater than or equal to three inches and less than five inches, per yard	1/wk	
A6452	High-compression bandage, elastic, knitted/woven, load resistance greater than or equal to 1.35 foot pounds at 50% maximum stretch, width greater than or equal to three inches and less than five inches, per yard	1/wk	
A6453	Self-adherent bandage, elastic, non-knitted/non-woven, width less than three inches, per yard	1/wk	
A6454	Self-adherent bandage, elastic, non-knitted/non-woven, width greater than or equal to three inches and less than five inches, per yard	1/wk	
A6455	Self-adherent bandage, elastic, non-knitted/non-woven, width greater than or equal to five inches, per yard	1/wk	
A6456	Zinc paste-impregnated bandage, non-elastic, knitted/woven, width greater than or equal to three inches and less than five inches, per yard	1/wk	
A6457	Tubular dressing with or without elastic, any width, per linear yard	*	
A6501	Compression burn garment, bodysuit (head to foot), customer fabricated	*	

* No established usual maximum, refer to Medicare Utilization Guidelines.

Medicare Utilization Chart

HCPCS	Long Description	Frequency of Change	Wound Condition
A6502	Compression burn garment, chin strap, customer fabricated	*	
A6503	Compression burn garment, facial hood, custom fabricated	*	
A6504	Compression burn garment, glove to wrist, custom fabricated	*	
A6505	Compression burn garment, glove to elbow, custom fabricated	*	
A6506	Compression burn garment, glove to axilla, custom fabricated	*	
A6507	Compression burn garment, foot to knee length, custom fabricated	*	
A6508	Compression burn garment, foot to thigh length, custom fabricated	*	
A6509	Compression burn garment, upper trunk to waist including arm openings (vest), custom fabricated	*	
A6510	Compression burn garment, trunk, including arms down to leg openings (leotard), custom fabricated	*	
A6511	Compression burn garment, lower trunk including leg openings (panty), custom fabricated	*	
A6512	Compression burn garment, not otherwise classified	*	
A6513	Compression burn mask, face and/or neck, plastic or equal, custom fabricated	*	
A6530	Gradient compression stocking, below knee, 18-30 mm HG, each	*	open venous ulcer
A6531	Gradient compression stocking, below knee, 30-40 mm HG, each	*	open venous ulcer
A6532	Gradient compression stocking, below knee, 40-50 mm HG, each	*	open venous ulcer
A6533	Gradient compression stocking, thigh length, 18-30 mm HG, each	*	open venous ulcer
A6534	Gradient compression stocking, thigh length, 30-40 mm HG, each	*	open venous ulcer
A6535	Gradient compression stocking, thigh length, 40-50 mm HG, each	*	open venous ulcer
A6536	Gradient compression stocking, full length/chap style, 18-30 mm HG, each	*	open venous ulcer
A6537	Gradient compression stocking, full length/chap style, 30-40 mm HG, each	*	open venous ulcer
A6538	Gradient compression stocking, full		open venous ulcer

* No established usual maximum, refer to Medicare Utilization Guidelines.

Ostomy Care
Ceramide Infused Products

Ostomy Care
New Image Two-Piece Pouching System

Ostomy Care
Premier One-Piece Pouching System

Ostomy Care
Pouchkins Pediatric Pouching System

Ostomy Care
Ostomy Accessories

Continence Care

Bowel Care

Tube Fasteners

Wound Care

Education

Medicare Utilization Chart

HCP/CS	Long Description	Frequency of Change	Wound Condition
	length/chap style, 40-50 mm HG, each		
A6539	Gradient compression stocking, waist length, 18-30 mm HG, each	*	open venous ulcer
A6540	Gradient compression stocking, waist length, 30-40 mm HG, each	*	open venous ulcer
A6541	Gradient compression stocking, waist length, 40-50 mm HG, each	*	open venous ulcer
A6544	Gradient compression stocking, garter belt	*	open venous ulcer
A6545	Gradient compression wrap, non-elastic, below knee, 30-50 mm HG, each	1/6 mo	open venous ulcer
A6549	Gradient compression stocking/sleeve, not otherwise specified	*	open venous ulcer
A9720	Non-covered item or service		

* No established usual maximum, refer to Medicare Utilization Guidelines.

Ostomy Care
Ceramide Infused
Products

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Ostomy Accessories

Continence
Care

Bowel Care

Tube
Fasteners

Wound Care

Education

FastCast OLS Total Contact Cast System

FastCast OLS Total Contact Cast*



- ▶ **Vascular Management** with low-level compression sock
- ▶ **Infection Management and Prevention** with the included *Hydrofera Blue Ready* foam and *Endoform* dermal template dressings
- ▶ **Pressure Relief** with *FastCast OLS* total contact cast system

Stock No	Product Description
FC0001	FastCast OLS 50" - 1 Cast Kit
FC0001-8	FastCast OLS 50" - 8 Cast Kits
FC0002	FastCast OLS 50" - 5 Cast Kits & 1 Med Boot
FC0003	FastCast OLS 50" - 5 Cast Kits & 1 Lg Boot
FC0006	FastCast OLS 56" - 1 Cast Kit
FC0006-8	FastCast OLS 56" - 8 Cast Kits
FC0007	FastCast OLS 56" - 5 Cast Kits & 1 Lg Boot
FC0008	FastCast OLS 56" - 5 Cast Kits & 1 XLg Boot
FC0011	FastCast OLS Small CAM Boot
FC0012	FastCast OLS Medium CAM Boot
FC0013	FastCast OLS Large CAM Boot
FC0014	FastCast OLS Extra Large CAM Boot
FC0024	FastCast OLS 4" Cast Roll Replacement

*Caution: Federal (USA) law restricts this device to sale by or on the order of a physician or licensed healthcare professional. Refer to Instructions for Use for contraindications, warnings, precautions, and possible complications.

Ostomy Care
Ceramide Infused
Products

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Ostomy Accessories

Continence
Care

Bowel Care

Tube
Fasteners

Wound Care

Education

TRIACT Technology Wound Interface

TRIACT technology is an innovative lipido-colloid formulation in select *Restore* wound dressings. The unique composition of *TRIACT* technology provides a moist wound interface and allows virtually pain-free removal. It also reduces trauma to delicate tissue upon removal. Dressings with *TRIACT* technology are designated with the *TRIACT* logo on the packaging.

- ▶ Helps to provide a moist, wound-healing environment
- ▶ Helps to minimize pain at dressing change¹
- ▶ Helps to prevent damage to newly formed tissue during dressing removal

TRIACT Foam with Silicone Border



- ▶ *TRIACT* Technology helps create a moist, wound-healing interface
- ▶ *TRIACT* Technology helps prevent damage to newly formed tissue during dressing removal
- ▶ Silicone border offers virtually pain-free removal

Stock No	HCPCS	Description	Quantity/Box
550762	A6212	4" x 4" with pad size 2.5" x 2.5"	10
550763	A6212	6" x 6" with pad size 4" x 4"	10
550764	A6213	6" x 8" with pad size 4" x 6"	10
550765	A6213	8" x 8" Sacral Pad 5" x 5.5"	10

Restore Foam with TRIACT Technology Non-Border



- ▶ Absorbent Foam Dressing with *TRIACT* Wound Interface Layer
- ▶ Requires Secondary Fixation
- ▶ Can be Cut to Size
- ▶ Indicated for Moderate to Heavy Exuding Acute and Chronic Wounds

Stock No	HCPCS	Description	Quantity/Box
509380	A6209	2.5" x 2.5"	10
509381	A6209	4" x 4"	10
509382	A6210	6" x 6"	10



509384	A6210	Heel 4.7" x 7.5"	10
--------	-------	------------------	----

Restore LITE Foam with TRIACT Technology Non-Border



- ▶ Thin, Conformable Non-Adhesive Absorbent Dressing with *TRIACT* Technology
- ▶ Manages Light to Moderate Exuding Acute and Chronic Wounds

Stock No	HCPCS	Description	Quantity/Box
509393	A6209	2.5" x 2.5"	10
509394	A6210	4" x 5"	10
509395	A6210	6" x 6"	10
509396	A6210	6" x 8"	10

¹ Meaume S. et al, The importance of pain reduction through dressing selection in routine wound management: the MAPP study. Journal of Wound Care November 2004, 13(10)

Foam Dressings

Restore Silver Foam with TRIACT Technology Non-Border*



- ▶ Non-Adhesive, Non-Occlusive Antimicrobial Absorbent Dressing with TRIACT
- ▶ May Help to Reduce Infection in Moderately to High Exuding Partial- and Full-Thickness Wounds
- ▶ Super-Absorbent Foam Pad Ensures Drainage of Exudate and Helps Protect the Skin Around the Lesion from Maceration
- ▶ Soft, Semi-Permeable Backing
- ▶ Contains Silver Sulfate

Stock No	HCPCS	Description	Quantity/Box
509345	A6209	4" x 4"	10
509346	A6210	6" x 8"	10

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician or licensed healthcare professional. Refer to Instructions for Use for contraindications, warnings, precautions, and possible complications.

Silicone Wound Interface

Restore Foam with Silicone Border



- ▶ Absorption of moderate to heavy exudate
- ▶ Protection of difficult wound-sites, peri-wound skin and bony prominences
- ▶ Simple application, easy removal & patient comfort

Stock No	HCPCS	Description	Quantity/Box
520003	A6212	3.5" x 3.5" with pad size 2.5" x 2.5"	10
520005	A6212	5" x 6.2" with pad size 3.4" x 4.6"	10
520000	A6213	7" x 7" with pad size 5.5" x 5.5"	10
520001	A6213	7.3" x 7.4" Sacral Pad 6.1" x 5.7"	10

Restore Foam with Silicone Non-Border



- ▶ Absorption of moderate to heavy exudate
- ▶ Protection of difficult wound-sites, peri-wound skin and bony prominences
- ▶ Simple application, easy removal & patient comfort

Stock No	HCPCS	Description	Quantity/Box
520023	A6209	2.5" x 2.5"	10
520025	A6210	4" x 5"	10
520027	A6211	6.5" x 8"	10



520020	A6210	Heel 8" x 5"	10
--------	-------	--------------	----

Restore Contact Layer FLEX



- ▶ TRIACT Formulation Promotes Non-Adherence, and Helps Prevent Disruption of the Healing Process During Dressing Removal
- ▶ Can Be Shaped to Fit Contours of Deep or Irregular Wounds, and in Difficult Anatomical Locations
- ▶ Virtually Pain-Free Dressing Removal for Enhanced Patient Comfort
- ▶ Forms a Gel in the Presence of Exudate to Help Maintain a Moist Wound Interface
- ▶ Dressing's Gentle Tack Helps Hold Dressing in Place and Minimizes Trauma to the Wound Bed
- ▶ Can Be Used on a Wide Variety of Wound Types, and in Conjunction with NPWT

Stock No	HCPCS	Description	Quantity/Box
506487	A6206	2" x 2"	10
506488	A6207	4" x 5"	10
506489	A6207	6" x 8"	10

Restore Contact Layer, Silver*



- ▶ Non-Adhesive, Non-Occlusive Antimicrobial Wound Contact Dressing with TRIACT Technology
- ▶ May Help to Reduce Infection in Low to Moderately Exuding Partial- and Full-Thickness Wounds
- ▶ Shown to be Effective Against Bacteria Most Frequently Associated with Wound Infections
- ▶ Apply Directly to the Wound and Cover with a Secondary Dressing
- ▶ Contains Silver Sulfate

Stock No	HCPCS	Description	Quantity/Box
509341	A6207	4" x 5"	10
509342	A6207	6" x 8"	10

*Caution: Federal (USA) law restricts this device to sale by or on the order of a physician or licensed healthcare professional. Refer to Instructions for Use for contraindications, warnings, precautions, and possible complications.

Restore to CalciCare Conversion Chart

Restore		CalciCare	
Description	SKUs	Description	SKUs
Restore Calcium Alginate 2"x 2"	529938	CalciCare Calcium Alginate 2"x 2"	529938R
Restore Calcium Alginate 4"x 4"	529937	CalciCare Calcium Alginate 4"x 4.75"	529937R
Restore Calcium Alginate 4"x 8"	529939	CalciCare Calcium Alginate 4"x 8"	529939R
Restore Calcium Alginate, 12" rope	529940	CalciCare Calcium Alginate, 18" rope	529940R
Restore Calcium Alginate, Silver 2"x 2"	529967	CalciCare Calcium Alginate, Silver 2"x 2"	529967R
Restore Calcium Alginate, Silver 4"x 4.75"	529968	CalciCare Calcium Alginate, Silver 4"x 4.75"	529968R
Restore Calcium Alginate, Silver, 12" rope	529969	CalciCare Calcium Alginate, Silver 18" rope	529969R

CalciCare and Restore Calcium Alginates

CalciCare Calcium Alginate



- ▶ Composed of a High G (Guluronic Acid) Calcium/Sodium Alginate
- ▶ Soft and Conformable, Highly Absorbent, and Fast Gelling
- ▶ Forms a Moist Environment for Optimal Wound Healing
- ▶ Non-Woven Pad is Reinforced for Intact Removal

Stock No	HCPCS	Description	Quantity/Box
529938R	A6196	2" x 2"	10
529937R	A6197	4" x 4.75"	10
529939R	A6197	4" x 8"	5
529940R	A6199	1" x 18" Rope	5

CalciCare Calcium Alginate, Silver*



- ▶ Composed of a High G (Guluronic Acid) Calcium/Sodium Alginate
- ▶ The Dressings are Soft and Conformable, Highly Absorbent, and Fast Gelling, as Wound Fluid is Absorbed
- ▶ Forms a Moist Environment for Optimal Wound Healing
- ▶ Non-Woven Pad is Reinforced for Intact Removal
- ▶ Releases Silver Ions in the Presence of Wound Exudate
- ▶ Effective Barrier to Bacterial Penetration

Stock No	HCPCS	Description	Quantity/Box
529967R	A6196	2" x 2"	10
529968R	A6197	4" x 4.75"	10
529969R	A6199	1" x 18" Rope	5

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician or licensed healthcare professional. Refer to Instructions for Use for contraindications, warnings, precautions, and possible complications.

Restore Calcium Alginate



- ▶ Highly Absorptive
- ▶ Made From Calcium Sodium Alginate, Creating Conformable Protective Gel When in Contact with Exudate
- ▶ Forms Moist Wound Interface
- ▶ Remains Intact When Saturated so it is Removed in One Piece

Stock No	HCPCS	Description	Quantity/Box
529938	A6196	2" x 2"	10
529937	A6196	4" x 4"	10
529939	A6197	4" x 8"	5
529940	A6199	12" Rope	5

Restore Calcium Alginate, Silver*



- ▶ Releases Silver Ions in the Presence of Wound Exudate
- ▶ Effective Antimicrobial Agent Against a Broad Spectrum of Microorganisms
- ▶ Indicated for Moderate to Heavily Exuding Partial- to Full-Thickness Wounds
- ▶ Assists in Maintaining a Moist Environment and Allows Intact Removal
- ▶ Contains Silver Sodium Hydrogen Zirconium Phosphate

Stock No	HCPCS	Description	Quantity/Box
529967	A6196	2" x 2"	10
529968	A6197	4" x 4.75"	10
529969	A6199	1" x 12" Rope	5

**Caution: Federal (USA) law restricts this device to sale by or on the order of a physician or licensed healthcare professional. Refer to Instructions for Use for contraindications, warnings, precautions, and possible complications.*

Restore Hydrogel



- ▶ Maintains Moist Healing Environment in Partial- and Full-Thickness Wounds
- ▶ Clear, Viscous Gel
- ▶ Up to 72-hour Wear Time

Stock No	HCPCS	Description	Quantity/Box
529974	A6248	Amorphous Gel in Tube 3 oz — Non Sterile	12

Ostomy Care
Ceramide Infused
Products

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Ostomy Accessories

Continence
Care

Bowel Care

Tube
Fasteners

Wound Care

Education

Restore Wound Cleanser and Skin Care

Restore Wound Cleanser



- ▶ Helps Remove Particulate Matter and Other Debris From Wound Bed
- ▶ No-Rinse Alternative to Saline and Harsh Antiseptic Cleansers

Stock No	HCPCS	Description	Quantity/Box
529975	A6260	8 oz pump spray	12
529976	A6260	12 oz pump spray	12

Restore Skin Cleanser



- ▶ For Use on Intact Skin
- ▶ Cleanses Peristomal and Perineal Skin of Urinary and Fecal Discharge
- ▶ Gentle, Non-Irritating
- ▶ Must Be Rinsed

Stock No	HCPCS	Description	Quantity/Box
517210	A9270	8 oz (237 mL) Pump Spray	12

Restore Cleanser & Moisturizer



- ▶ Effective Skin Protection for People with Urinary and/or Fecal Incontinence
- ▶ Cleanser, Moisturizer, Protective Breathable Barrier for Skin
- ▶ Dimethicone-Based
- ▶ No Rinse

Stock No	HCPCS	Description	Quantity/Box
529977	A6250	7.5 oz (222 mL) Pump Spray	12
529978	A6250	11.5 oz (340 mL) Pump Spray	12

Restore Skin Conditioning Creme



- ▶ Moisturizes and Soothes Dry, Flaky Skin
- ▶ May be Used Under Adhesive Products

Stock No	HCPCS	Description	Quantity/Box
517220	A6250	4 oz (120 mL) Bottle	12

Ostomy Care
Ceramide Infused
Products

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Ostomy Accessories

Continence
Care

Bowel Care

Tube
Fasteners

Wound Care

Education

Restore Hydrocolloids

Restore Hydrocolloid: Sterile



- ▶ Conformable Semi-occlusive Dressing with Flexible Outer Layer
- ▶ Maintains Moist Wound Environment
- ▶ Manages Minimal to Moderate Exudate in Partial- and Full-Thickness Wounds
- ▶ Available With or Without Tapered Edges, Specialty Shapes
- ▶ Includes Disposable Wound Measuring Guide

Stock No	HCPCS	Description	Quantity/Box
519953	A6234	Without Tapered Edges 4" x 4"	5
519954	A6235	Without Tapered Edges 6" x 8"	3
519955	A6236	Without Tapered Edges 8" x 8"	3
519956	A6234	With Tapered Edges 4" x 4"	5
519963	A6235	With Tapered Edges 6" x 6"	5
519957	A6235	With Tapered Edges 6" x 8"	3
519958	A6236	With Tapered Edges 8" x 8"	3
519959	A6235	With Tapered Edges/Triangle Shaped 17 sq in	5
519965	A6235	With Tapered Edges/Triangle Shaped 26.5 sq in	5

Restore Extra Thin Hydrocolloid: Sterile



- ▶ Thinner Semi-occlusive Dressing with Flexible Backing
- ▶ Maintains Moist Wound Environment
- ▶ Manages Minimal Exudate in Shallow Wounds
- ▶ Helps Protect Skin from Friction Injury
- ▶ Includes Disposable Wound Measuring Guide

Stock No	HCPCS	Description	Quantity/Box
519921	A6234	4" x 4"	5
519923	A6235	6" x 8"	3
519925	A6236	8" x 8"	3

Restore Hydrocolloid with Foam Backing: Sterile



- ▶ Occlusive Dressing with White Foam Backing
- ▶ For Low Friction Areas with Moderate Exudate
- ▶ Includes Disposable Wound Measuring Guide

Stock No	HCPCS	Description	Quantity/Box
519930	A6234	4" x 4"	5
519932	A6235	6" x 8"	3
519935	A6236	8" x 8"	3



DID YOU KNOW...

Disposable Wound Guides



Disposable Wound Measuring Guides are included with all **Restore** hydrocolloid dressings. To order additional Wound Measuring Guides, please refer to the Educational Resources section of this catalog.

Education

Wound Care

Tube
Fasteners

Bowel Care

Continence
Care

Ostomy Care
Ostomy Accessories

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Ceramide Infused
Products

Ostomy Care
Ceramide-Infused
Products

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Ostomy Accessories

**Continence
Care**

Bowel Care

**Tube
Fasteners**

Wound Care

Education

Education

Education

Wound Care

**Tube
Fasteners**

Bowel Care

**Continence
Care**

Ostomy Care
Ostomy Accessories

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Ceramide Infused
Products

AF300 Filter: “AF” stands for air flow and “300” represents the 300 mm² air flow surface area on every AF300 filter.

Ceramide: A natural component of the human skin that links the cells of the stratum corneum together to form a waterproof protective barrier. This helps prevent water loss that can lead to skin damage and dryness.

CeraPlus Skin Barrier: A skin barrier that has all the basic characteristics of an ostomy skin barrier (adhesion, absorption, and erosion resistance) and is infused with ceramide. *CeraPlus* is an example of an extended wear skin barrier.

Closed Pouch: A pouch without an opening at the bottom. Quick and easy to discard. It must be removed to be emptied.

Colostomy: An ostomy (surgical opening) created in the colon, part of the large intestine.

ComfortWear (pouch) Panels: Provide a soft cloth between the pouch film and the skin.

Convex Skin Barrier: A skin barrier shaped to press more deeply around the stoma. May be used when the stoma does not protrude from the abdomen.

Cut-to-Fit Skin Barrier: A skin barrier that can be customized to a needed size and shape. Ideal for stomas that are not round or stomas that are still changing in size.

Drainable Pouch: A pouch with an opening at the bottom. An integrated closure system or plastic clamp is used to keep the pouch closed until it is time to empty it.

Filter: Releases gas but not odor out of the pouch. The key performance factors of any filter are air flow, deodorization, and liquid protection.

Flange: The plastic ring in a two-piece ostomy system that is used to connect the skin barrier and pouch together.

Flat Skin Barrier: May be used when the stoma protrudes from the abdomen.

FlexTend Skin Barrier: A skin barrier designed to be durable and resistant to liquid discharge. *FlexTend* is an example of an extended wear skin barrier.

Ostomy Education • Glossary

FlexWear Skin Barrier: A skin barrier designed to be worn more than one day. A *FlexWear* skin barrier is designed to be more resistant to breakdown than a *SoftFlex* skin barrier. *FlexWear* is an example of a standard wear skin barrier.

Floating Flange: Type of skin barrier flange found on Hollister two-piece pouching systems that allows you to put your fingers under the skin barrier flange while attaching the pouch.

FormaFlex Skin Barrier: A skin barrier that can be stretched to fit the size and shape of the stoma without the use of scissors. It is a barrier option in the New Image two-piece pouching system. *FormaFlex* is an example of an extended wear skin barrier.

Ileostomy: A surgically created opening in the last portion of the small intestine.

Irrigation: A procedure sometimes used by individuals with a colostomy to control bowel elimination. It is not for everyone with a colostomy and requires special equipment and instruction by a knowledgeable healthcare provider.

Lock 'n Roll Microseal Closure: Integrated system for closing the bottom of drainable pouches. Utilizes plastic interlocking fasteners to create a secure closure. Found on *New Image*, *Premier*, and *Pouchkins* drainable pouches.

Medical Adhesive-Related Skin Injury (MARS): Occurrences in which erythema and/or other manifestation such as blister, tear or other injury persists more than 30 minutes after removal of medical adhesive.

Odor-Barrier Pouch: Is highly resistant to letting odors pass through the pouch film. The odor barrier quality increases confidence when wearing a pouch.

Peristomal Moisture-Associated Skin Damage (PMASD): Inflammation and erosion of the skin adjacent to the stoma, associated with exposure to effluent such as urine or stool.

Pre-Sized Skin Barrier: A skin barrier that is already cut to the proper size.

QuietWear Pouch Material: Material comprised of a pouch film and *ComfortWear* panel that are laminated together. This special design results in the reduction of the pouch rustling noise.

SoftFlex Skin Barrier: A skin barrier that provides gentle adhesive attachment to the skin. *SoftFlex* is designed to limit skin damage with frequent removal and can be changed daily (or more often). *SoftFlex* is an example of a standard wear skin barrier.

Stoma: A surgically created opening in the gastrointestinal or urinary tract. Also known as an ostomy. The stoma is pink and moist and doesn't have feeling.

Urostomy: A surgically created opening into the urinary system; drains urine and sometimes mucus.

Urostomy Pouches: Pouch with a drain tap at the bottom so urine can be emptied quickly and easily.

Wear Time: The length of time a product can be worn before failure. Wear time varies widely based upon a variety of factors.

Key Goals for Nursing Management

- **Protect the skin** – The output from a fistula can be very irritating to the skin. Broken skin contributes to pain and discomfort. A skin barrier can protect the skin from contact with the drainage.
- **Contain the output** – If the output is minimal, gauze dressings may be sufficient for containment. However, if the volume of the output from the fistula is significant or frequently overwhelms the dressing, a pouch may be more appropriate. Using a pouch can be more efficient and can protect healthcare providers from exposure to bodily fluids.
- **Measure the output** – With a high output fistula, loss of fluids and electrolytes can be a major concern. Measuring the drainage is easier using a pouching system. Sometimes connecting a pouch to a bedside collector can relieve the stress on the pouch, so the wear time improves.
- **Manage costs** – Care of fistulas can be costly in terms of supplies and human resources. Using a pouching system and appropriate accessories can be more cost-effective and manageable than frequent dressing changes.
- **Control odor** – Containment of the drainage in a pouching system, even in small quantities, can manage odor effectively. Other accessory products can also assist with odor control when the pouch is changed or emptied.
- **Provide patient comfort** – Effectively managing a fistula can contribute to patient dignity. Protecting the skin, collecting the discharge, containing costs, and controlling odor all contribute to patient comfort.

All of the Hollister educational materials listed below are also available on our website: www.hollister.com

Educational DVDs	Literature No
DVD 1: A Guide to Your Ostomy Surgery (English)	921856
DVD 1: A Guide to Your Ostomy Surgery (Spanish)	922224
DVD 2: A Guide to Living with Your Ostomy (English)	921857
DVD 2: A Guide to Living with Your Ostomy (Spanish)	922225

Online Education Program *available at www.hollister.com*

Case Studies

For a complete description and current listing, visit the Clinical Evaluations and Case Studies section found in the Professional Resources portion of our website, www.hollister.com. For more information, contact your Hollister Ostomy Sales Representative.

Stoma Measuring Guide	903515
------------------------------	--------

Educational Booklets	Literature No
Understanding Your Colostomy	923054
Understanding Your Colostomy (Spanish) <i>(online only)</i>	923234
Understanding Your Ileostomy	923055
Understanding Your Ileostomy (Spanish) <i>(online only)</i>	923235
Understanding Your Urostomy	923056
Understanding Your Urostomy (Spanish) <i>(online only)</i>	923236
Can They Still Wiggle and Giggle? Coloring Book	907617
The Day I Got My Stoma <i>(online only)</i>	910497
Caring for Your Child With an Ostomy – Pediatric (English)	923057
Caring for Your Child With an Ostomy – Pediatric (Spanish) <i>(online only)</i>	923257
Caring for your Loved One with an Ostomy	923058
Caring for your Loved One with an Ostomy (Spanish) <i>(online only)</i>	923258

Ostomy Care Tips

<i>Adapt</i> Barrier Rings	923059
<i>Adapt</i> Convex Barrier Rings	923060
<i>Adapt</i> Lubricating Deodorant	923061
<i>Adapt</i> Medical Adhesive Spray/ <i>Adapt</i> Remover Spray & Wipes	923073
<i>Adapt</i> No-Sting Skin Protective Wipe	923079
<i>Adapt</i> Ostomy Belt	923062
<i>Adapt</i> Paste	923063
<i>Adapt</i> Stoma Powder	923064
Colostomy Irrigation (<i>online only</i>)	923065
Using Convexity	923066
Fistula Management (<i>online only</i>)	923067
High Output	923121
Ileostomy Dietary Guidelines	923117
<i>Lock 'n Roll</i> Microseal Closure	923068
Medicare Reimbursement Guide	923118
One-Piece Pouching System (English)	923069
One-Piece Pouching System (Spanish)	923070
Ostomy Educational Theatre (<i>online only</i>)	923115
Patient Anatomy Guide	923072
Peristomal Skin Care	923085
<i>Pouchkins</i> Newborn Pouching System (<i>online only</i>)	923119
<i>Pouchkins</i> One-Piece Pouching System (<i>online only</i>)	923074
<i>Pouchkins</i> Premie Pouching System (<i>online only</i>)	923076
<i>Pouchkins</i> Two-Piece Pouching System (<i>online only</i>)	923075
Routine Care of Your Ostomy (English)	923077
Routine Care of Your Ostomy (Spanish)	923078
Stoma Site Selection (<i>online only</i>)	923080
Two-Piece Pouching System (English)	923081
Two-Piece Pouching System (Spanish)	923082
Two-Piece Pouching System with <i>FormaFlex</i> Skin Barrier (<i>online only</i>)	923083
Types of Ostomies	923084
Universal Catheter Access Port (UCAP) (<i>online only</i>)	923114
Where a Soft Convex Skin Barrier May Improve Fit	923116

Living with an Ostomy Booklets

	Literature No
Travel	923123
Healthy Eating	923124
Sex & Parenthood	923125
Home & Work Life	923126
Sports & Fitness	923127

Professional Organizations:

Wound, Ostomy and Continence Nurses Society (WOCN)
www.wocn.org
1.888.224.9626

Shadow Buddies	Item No
Light Skin Tone Girl	905439
Light Skin Tone Boy	905440
Medium Skin Tone Girl	905441
Medium Skin Tone Boy	905442
Dark Skin Tone Girl	905443
Dark Skin Tone Boy	905444

Additional Resources	Item No
Stoma Model (White Plastic)	858674

Consumer Resources

Hollister Secure Start Service
1.888.808.7456
www.securestartservices.com

United Ostomy Associations of America, Inc. (UOAA)
info@uoaa.org
www.ostomy.org
1.800.826.0826

Crohn's and Colitis Foundation of America (CCFA)
www.ccfa.org
1.800.932.2423

Shadow Buddies
www.shadowbuddies.org
1.913.642.4646

Contenance Care Education • Glossary

Anti-Reflux: Refers to a tube or collapsible material within a urine collection device to help prevent urine from re-entering the tubing

Applicator Collar: Found on some Hollister male external catheters; a plastic guide with notches for the thumb and forefinger to assist proper placement against the tip of the penis

Bladder: A balloon-like muscular organ which lies in the pelvis and functions to store and expel urine

Bladder Catheterization: A procedure in which a catheter is passed through the urethra or stoma into the bladder for the purpose of draining urine and performing diagnostic tests of bladder or urethral function

Benzalkonium Chloride (BZK): An antimicrobial solution used for disinfecting the urethral opening prior to inserting an intermittent catheter. Does not stain skin or clothing. Less likely to cause sensitivity

Catheter (Urinary): A hollow tube inserted through the urethra or a stoma to the bladder to withdraw urine or administer medication

Closed System: Refers to a "touchless" gel catheter located within a urine collection bag. The enclosed catheter is manipulated through the bag, guided through the protective tip as it is inserted into the urethra or stoma without being directly touched. For VaPro Plus, the catheter is placed outside the collection bag and the funnel of the catheter attaches to the collection bag port.

Closed System Kit: Closed system catheter includes the catheter located within the urine collection bag, and usually includes gloves, PVP or BZK swab sticks, and underpad drape

Collection Bag: A leg bag or bedside bag designed to collect urine by gravity that is connected by tubing to an indwelling or external catheter, urinary pouch, or also to wound drainage systems

Color Coded: On Hollister intermittent catheters, the funnel ends are color and size coordinated so catheters may be consistently ordered by size or funnel color

Connector: A device that secures the extension tubing to the catheter or urinary pouch system

Coudé Tip: A slight bend manufactured in the tip of the catheter which makes insertion past the prostate easier for some men. A notch at the funnel end is a guide for alignment during insertion

External Condom Catheter: Device that is secured externally with various adhesives to the penis shaft for the purpose of urine collection. Must be connected with tubing to a urine collection bag

Extension Tubing: Tubing that provides the connection between an external or indwelling catheter and a urine collection bag or leg bag. May also serve as a connection between an ostomy pouch or wound drainage collector. Made of latex or vinyl

Eyelet: An oval-shaped hole in the insertion end of the catheter to facilitate drainage of urine from the bladder. Eyelets are usually two in number, and may be across from each other or offset. Smooth eyelets help protect the urethra from insertion trauma

Female Urinary Pouch: A cut-to-fit, one-piece, vinyl pouch with a flexible synthetic barrier which is attached externally to the female urethral/vaginal area. The pouch is usually connected to a bedside collection bag. Recommended for the more immobile female

Firm Catheter: Refers to lesser pliability of a urethral catheter

Flexible: Refers to greater pliability of a urethral catheter or of a skin barrier found on the external female and male urinary pouches

Flextend Barrier: A skin barrier with special additives that achieve a stronger adhesive seal and are more resistant to breakdown from fluids. Skin-prepping agents are not recommended under *Flextend* Skin Barriers

French Size: Abbreviated Fr, the measuring gauge for the outer diameter of a straight or indwelling catheter

Funnel: Colored non-insertion end of an intermittent catheter allows for ease of fluid control. A catheter tipped syringe may be attached at the funnel for fluid instillation into the bladder. Color corresponds to specific French size

Gel Reservoir: A small, lubricant-filled flexible device unique to Hollister built into the *Advance Plus* touch free intermittent catheter system. The catheter is evenly coated with gel as it passes through the gel reservoir

Grip Section: A finely ribbed area on the catheter in a closed catheter system which assists with catheter control for continuous catheter advancement

Hydrophilic Catheter: A coated catheter that is lubricated with water

Incontinence: Loss of natural control of bowel and/or bladder functions

Indwelling Catheter: A flexible tube that remains in the bladder continuously to drain urine. May be referred to as a Foley catheter

Inner Flap: A thin, latex extra membrane inside Hollister *Extended Wear* latex male external catheters. The inner flap rests against the tip of the penis to prevent urine backwash and undermining of the catheter adhesive

Intermittent Catheter: A flexible tube that is used to empty the bladder on a regular schedule; used for self-catheterization. For single use only.

Latex: A material made from natural rubber. Red rubber catheters and some male external catheters are made from latex

Leg Bag Straps: Hold the leg bag in place

Luer End: No funnel on non insertion end. A syringe may be connected to the Luer end for bladder irrigation or insertion of medication

Lubricant: A water soluble jelly applied to a catheter prior to insertion which allows for easier insertion and helps minimize urethral irritation

Male External Urinary Pouch: A cut-to-fit, one-piece vinyl pouch with a synthetic adhesive backing that may be used for the very short penis, or the flush or retracted penis. Usually attached to a bedside collection bag, and is recommended more for those with limited mobility

Meatus: The opening of the urethra in both men and women

Neurogenic Bladder (also known as, Neurogenic Lower Urinary Tract Dysfunction): A loss of normal neural control of the bladder and sphincter function associated with a neurological condition, such as diabetes, stroke, multiple sclerosis, spina bifida or spinal cord injury

Non Latex: A material not made with Natural Rubber Latex; usually vinyl or silicone for catheters or plastic for collection bags. Minimizes risk of allergic reaction

Pre-Lubricated Catheter: Refers to a catheter which is lubricated by passing through a gel reservoir, or a lubricated catheter in a closed system, or a catheter in packaging which contains a lubricant packet which is ruptured to lubricate the catheter prior to opening the packaging, or a pre-lubricated catheter with hydrophilic technology used in the packaging allows the catheter to be lubricious and ready to use out of the pack

Protective Tip: A specially designed system on some intermittent catheters which functions to prevent urethra bacterial contamination of the catheter. It is inserted into the urethra. The 15 mm length on the Hollister protective tip helps decrease the risk of introducing bacteria into the urethra. The tapered tip with multiple slits is designed for comfortable insertion

Ring Cap: Cap on the protective tip protects the catheter tip. Ring assists with cap removal especially for those with limited dexterity. Oval cap design prevents roll away

Sizing Guide: Used for external catheters, these semicircular cutouts help determine the diameter of the penis in order to guide selection for correct fit

Soft Catheter: A gentle, flexible catheter with limited firmness

SoftFlex Barrier: A flexible synthetic material that will conform to rounded and some irregular surfaces. Gentle enough to be removed in 24 hours or less

Straight Tip: Refers to the tapered, rounded insertion end of a catheter

Urethra: The natural tube through which urine exits the bladder

Vacuum Relief Valve: A feature on the leg bag tubing that helps prevent collapse of the tubing, as well as collapse of an external urinary pouch or catheter

Vented Leg Bag: A special vacuum relief mechanism featured on the oval kink-resistant tubing and the companion pouch which also features an air vent. These combined features help to minimize a vacuum in the leg bag as well as in the tubing. A vacuum is created as the urine cools outside the body. The vacuum will stall urine drainage through the tubing and/or prevent the leg bag from emptying completely

Visit www.hollister.com/continencecare



SEARCH

PRODUCTS OSTOMY CARE CONTINENCE CARE WOUND & SKIN CARE CRITICAL CARE SECURE START SERVICES ABOUT

Continence Care Products

A listing of bladder management products

Learning Center

Guidance beyond product information

Professional Resources

Educational resources for healthcare professionals



Hollister Continence Care

People are at the center of everything we do

Ostomy Care
Ceramide Infused
Products

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Ostomy Accessories

Continence
Care

Bowel Care

Tube
Fasteners

Wound Care

Education

Educational Materials

See hollister.com for a current listing of educational materials.

Video Resources

See hollister.com for a current listing of educational materials.

Specialty Organizations

American Association for Respiratory Care (AARC)

www.aarc.org

1.972.243.2272

Wound, Ostomy and Continence Nurses Society (WOCN)

www.wocn.org

1.888.224.9626

Society of Critical Care Medicine (SCCM)

www.sccm.org

1.847.827.6869

National Pressure Ulcer Advisory Panel

www.npuap.org

1.202.521.6789

American Association of Critical Care Nurses (AACN)

www.aacn.org

1.800.394.5995

<http://www.hollister.com/en/criticalcare/criticalcareprofessionalresources/clinicalevaluationscasestudies>

Abrasion: Wearing away of the skin through some mechanical process (friction or trauma)

Abscess: Accumulation of pus formed in tissue as a result of infection

Alginate: A highly absorptive dressing derived from brown seaweed

Cellulitis: Inflammation of the tissues indicating a local infection; characterized by redness, edema, and tenderness

Collagen: Protein found in skin and connective tissue

Debridement: Removal of foreign material and devitalized or contaminated tissue from a wound

Dehiscence: Separation of wound edges

Denude: Removal or loss of superficial skin layers

Dressing: Any of various material used for covering and protecting wounds

Edema: Swelling

Epidermis: Outermost layer of the skin

Erythema: Diffuse redness of the skin

Eschar: Dead tissue

Excoriation: A scratch or abrasion on the skin

Exudate: Wound fluid or drainage

Friction: Rubbing that causes mechanical trauma to the skin or tissue, layers, mucosa

Full-thickness: Tissue destruction extending through the dermis to involve subcutaneous level and possibly muscle, fascia, or bone

Granulation: Formation of connective tissue and many new capillaries in a full-thickness wound; typically appears as red and cobblestoned

Hydrocolloid dressing: A category of wound dressings composed of materials, such as gelatin, pectin, and carboxymethylcellulose, that provide a moist healing environment and adhere to the skin around the wound

Hydrogel: Water- or glycerin-based gels, impregnated gauzes, or sheet dressings. Hydrogels maintain a moist healing environment and absorb a minimal amount of wound exudate

Hydrophilic: Attracting moisture

Infection: Invasion and multiplication of microorganisms in the body tissues

Maceration: Softening of tissue by soaking in fluids

Necrotic: Devitalized tissue; may appear yellow and moist, gray, or dark and leathery

Partial-thickness: Wounds that extend through the epidermis and may involve the dermis; these wounds heal by re-epithelialization

Peri-wound: The area immediately around the wound

Pressure Injury Stages:

Pressure Injury: A pressure injury is localized damage to the skin and underlying soft tissue usually over a bony prominence or related to a medical or other device. The injury can present as intact skin or an open ulcer and may be painful. The injury occurs as a result of intense and/or prolonged pressure or pressure in combination with shear. The tolerance of soft tissue for pressure and shear may also be affected by microclimate, nutrition, perfusion, co-morbidities and condition of the soft tissue.

Stage 1 Pressure Injury: Non-blanchable erythema of intact skin:

Intact skin with a localized area of non-blanchable erythema, which may appear differently in darkly pigmented skin. Presence of blanchable erythema or changes in sensation, temperature, or firmness may precede visual changes. Color changes do not include purple or maroon discoloration; these may indicate deep tissue pressure injury.

Stage 2 Pressure Injury: Partial-thickness skin loss with exposed dermis

Partial-thickness loss of skin with exposed dermis. The wound bed is viable, pink or red, moist, and may also present as an intact or ruptured serum-filled blister. Adipose (fat) is not visible and deeper tissues are not visible. Granulation tissue, slough and eschar are not present. These injuries commonly result from adverse microclimate and shear in the skin over the pelvis and shear in the heel. This stage should not be used to describe moisture associated skin damage (MASD) including incontinence associated dermatitis (IAD), intertriginous dermatitis (ITD), medical adhesive related skin injury (MARS), or traumatic wounds (skin tears, burns, abrasions).

Stage 3 Pressure Injury: Full-thickness skin loss

Full-thickness loss of skin, in which adipose (fat) is visible in the ulcer and granulation tissue and epibole (rolled wound edges) are often present. Slough and/or eschar may be visible. The depth of tissue damage varies by anatomical location; areas of significant adiposity can develop deep wounds. Undermining and tunneling may occur. Fascia, muscle, tendon, ligament, cartilage and/or bone are not exposed. If slough or eschar obscures the extent of tissue loss this is an Unstageable Pressure Injury.

Stage 4 Pressure Injury: Full-thickness skin and tissue loss

Full-thickness skin and tissue loss with exposed or directly palpable fascia, muscle, tendon, ligament, cartilage or bone in the ulcer. Slough and/or eschar may be visible. Epibole (rolled edges), undermining and/or tunneling often occur. Depth varies by anatomical location. If slough or eschar obscures the extent of tissue loss this is an Unstageable Pressure Injury.

Unstageable Pressure Injury: Obscured full-thickness skin and tissue loss

Full-thickness skin and tissue loss in which the extent of tissue damage within the ulcer cannot be confirmed because it is obscured by slough or eschar. If slough or eschar is removed, a Stage 3 or Stage 4 pressure injury will be revealed. Stable eschar (i.e. dry, adherent, intact without erythema or fluctuance) on the heel or ischemic limb should not be softened or removed.

Pressure Injury Stages: (continued)

Deep Tissue Pressure Injury: Persistent non-blanchable deep red, maroon or purple discoloration

Intact or non-intact skin with localized area of persistent non-blanchable deep red, maroon, purple discoloration or epidermal separation revealing a dark wound bed or blood filled blister. Pain and temperature change often precede skin color changes. Discoloration may appear differently in darkly pigmented skin. This injury results from intense and/or prolonged pressure and shear forces at the bone-muscle interface. The wound may evolve rapidly to reveal the actual extent of tissue injury, or may resolve without tissue loss. If necrotic tissue, subcutaneous tissue, granulation tissue, fascia, muscle or other underlying structures are visible, this indicates a full thickness pressure injury (Unstageable, Stage 3 or Stage 4). Do not use DTPI to describe vascular, traumatic, neuropathic, or dermatologic conditions.

Additional pressure injury definitions.

Medical Device Related Pressure Injury: This describes an etiology.

Medical device related pressure injuries result from the use of devices designed and applied for diagnostic or therapeutic purposes. The resultant pressure injury generally conforms to the pattern or shape of the device. The injury should be staged using the staging system.

Mucosal Membrane Pressure Injury: Mucosal membrane pressure injury is found on mucous membranes with a history of a medical device in use at the location of the injury. Due to the anatomy of the tissue these ulcers cannot be staged.

Pus: Thick fluid composed of leukocytes, bacteria, and cellular debris

Shear: Sliding of skin over subcutaneous tissues and bones obstructing cutaneous capillaries, which may lead to ischemia

Sinus Tract: A course or pathway which can extend in any direction from the wound base; results in dead space with potential for abscess formation

Slough: Stringy, necrotic tissue; usually yellow, tan, gray, green, or brown

Strip: Removal of epidermis by mechanical means, usually tape

Ulcer: Loss of epidermis/dermis or mucous membrane with definite margins

Undermine: Skin edges of a wound that have lost supporting tissue under intact skin

Wound: A break in the integrity of the skin; an injury to the body which causes a disruption of the normal continuity of the body structures

Wound Margin: Rim or border of a wound

Skin Structure:

The **epidermis**, which is the outermost layer of the skin, is characterized as follows: avascular, varies in thickness (depending on body location), a dry structure which sheds cells and replaces itself every 4 to 6 weeks; approximately the thickness of a piece of plastic wrap

The **dermis** is located directly beneath the epidermis and is characterized as follows: provides strength and structural support through a vascular network of connective tissue blood vessels, nerves, hair, and nails. Sebaceous glands and sweat glands originate from this layer, which is thicker than the epidermis

Below the dermis is the **subcutaneous tissue** which is composed of major vessels, lymphatics, fat, and connective tissue. This area provides insulation and nutritional support for the skin. Located below the subcutaneous tissue are fascia, muscles, tendons, and bone. The thickness of the dermis and subcutaneous layers vary from person to person and on different parts of the body

Types of Ulcers

Pressure Ulcers: A pressure ulcer is localized injury to the skin and/or underlying tissue usually over a bony prominence as a result of pressure, or pressure in combination with shear and/or friction. A number of contributing or confounding factors are also associated with pressure ulcers; the significance of these factors is yet to be elucidated.

Venous Ulcers: These are the most common type of lower-leg ulcers in ambulatory people. The underlying cause involves vein damage or an incompetent calf muscle pump action which leads to venous hypertension. As a result, the blood pools in the lower extremities causing edema and leakage of fibrinogen and other blood products into the tissues. Trauma to the area or increased pressure within the tissues results in ulceration.

Arterial Ulcers: Arterial ulcers result from chronic or acute arterial insufficiency to the skin and subcutaneous tissue of the lower extremities. The most common cause is a progressive disease: atherosclerosis. The precipitating event leading to ulceration is usually trauma, such as a bumped toe or tight shoes. Arterial ulcers may occur alone or in combination with diabetes, venous stasis, and numerous other conditions. Multidisciplinary management of these patients with early intervention and close monitoring is key to prevention of more serious complications.

Neuropathic Ulcers: Neuropathic ulcers may occur in individuals with diabetes, spinal cord injury, Hansen's disease, or other conditions that result in loss of sensation in the legs and feet. Diabetic foot ulcers are most commonly caused by peripheral neuropathy and peripheral vascular disease. Multidisciplinary management of these patients with early intervention and close monitoring is key to prevention of more serious complications.

Reference: Miller-Keane Encyclopedia & Dictionary of Medicine, Nursing & Allied Health - Revised Reprint 7th Edition <http://www.us.elsevierhealth.com/dictionary/miller-keane-encyclopedia-dictionary-of-medicine-nursing-allied-health-revised-reprint-hardcover/9781416026044/>

The ISTAP Skin Tear Resource Kit has been designed to help implement a systematic approach to the prevention, management, and treatment of skin tears.

Hollister is proud to support the International Skin Tear Advisory Panel in the prevention, assessment, and treatment of skin tears. The support of Hollister has been provided through an unrestricted educational grant.

This information is provided for general medical educational purposes only.


More information is available at www.skintears.org



Skin Tears Defined


A wound caused by shear, friction, and/or blunt force resulting in separation of skin layers. A skin tear can be partial-thickness (separation of the epidermis from the dermis) or full-thickness (separation of both the epidermis and dermis from underlying structures.)¹

Skin Tear Risk Assessment Pathway²



Assess and Manage

Skin Tear Decision Algorithm³



References:

- 1 LeBlanc, K., Baranoski, S. et al. Skin Tears: State of the Science: Consensus Statements for the Prevention, Prediction, Assessment, and Treatment of Skin Tears. *Advances in Skin and Wound Care*. 2011;24(9Suppl):2-15.
- 2 Ellman, C. et al. (2013) Sustained Results: A learning management system platform which empowers and educates bedside staff.
- 3 LeBlanc K., Baranoski S., et al. International Skin Tear Advisory Panel: A Tool Kit to Aid in the Prevention, Assessment, and Treatment of Skin Tears Using a Simplified Classification System. *Advances in Skin and Wound Care*. 2013;26:459-476.



DID YOU KNOW...

Partnering for the Prevention and Management of Skin Tears.

We are proud to support the International Skin Tear Advisory Panel in the prevention, assessment, and treatment of skin tears. The support of Hollister has been provided through an unrestricted educational grant.

An ISTAP Skin Tear Resource Kit has been designed to help implement a systematic approach to the prevention, management, and treatment of skin tears.

This information is provided for general medical educational purposes only. More information is available at www.skintears.org.



Educational Materials

<i>FastCast OLS</i> Brochure	922780
Contact Layer Brochure	921818
<i>CalciCare</i> Calcium Alginate Brochure	922571
Skin Tears: State of the Science: Consensus Statements for the Prevention, Prediction, Assessment, and Treatment of Skin Tears	921822

Measuring Tools

<i>Restore</i> Wound Measuring Guide, Small	874366
<i>Restore</i> Wound Measuring Guide, Large	874367

Product Inservice Tools are available through your Hollister sales representative

Algorithms, Protocols, and Other Educational Tools are available on our website. For a complete description, visit the Educational Resources portion of our website, or contact your Hollister sales representative.

Ostomy Care
Ceramide Infused
Products

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Ostomy Accessories

Continence
Care

Bowel Care

Tube
Fasteners

Wound Care

Education

Clinical Poster Reprints

Literature No

Restore Contact Layer

Use of Negative Pressure Wound Therapy with a New Wound Contact Layer: Results of a Clinical Evaluation 907982

Restore Contact Layer, Silver

Use of a Contact Layer with Silver Under a Negative Pressure Dressing 910278

Specialty Organizations

The Association for the Advancement of Wound Care (AAWC)
www.aawconline.org
1.866.229.9999

Wound, Ostomy and Continence Nurses Society (WOCN)
www.wocn.org
1.888.224.9626

American Academy of Wound Management (AAWM)
abwmcertified.org
1.202.521.0368

American Burn Association
www.ameriburn.org
1.312.642.9260

Wound Healing Society
www.woundheal.org
1.763.765.2377

National Pressure Ulcer Advisory Panel
www.npuap.org
1.202.521.6789

American Diabetes Association
www.diabetes.org
1.800.DIABETES (342.2383)

National Spinal Cord Injury Association
www.spinalcord.org
1.800.962.9629

U.S. Department of Health and Human Services Centers for Medicare & Medicaid Services
www.cms.hhs.gov

National Guideline Clearinghouse
www.guideline.gov

European Wound Management Association
www.ewma.org
+45 70 20 03 05

Stock No	Stock No	Stock No	Stock No
1xxx	1837x	3777	520020
10xx	1838x	3778	520023
1120x	18382	3795	520025
1140x	1839x	3796	520027
1150x	18392	3797	529937
1340x	1840x	3799	529937R
1350x	18402	380x	529938
139xx	1841x	381x	529938R
1410x	18412	382x	529939
1420x	1842x	383x	529939R
143x	18422	384x	529940
1430x	1873x	387x	529940R
1440x	1875x	389x	529967
1450x	1890x	390x	529967R
1460x	1891x	506487	529968
1470x	1892x	506488	529968R
1480x	1900x	506489	529969
149xx	1905x	509341	529969R
1510x	1910x	509342	529974
1520x	1915x	509345	529975
1530x	1920x	509346	529976
1540x	1925x	509380	529977
1550x	1930x	509381	529978
1560x	1935x	509382	550762
1580x	1940x	509384	550763
1590x	1950x	509393	550764
1610x	31003	509394	550765
1640x	31005	509395	71xxx
1796	31009	509396	716x
1800x	318x	517210	7299
1801x	32004	517220	73xxx
1810x	32005	519921	7300
1811x	32006	519923	7331
1812x	32007	519925	74xxx
1813x	322x	519930	7700
1814x	323x	519932	771x
1815x	327x	519935	7721
1816x	33016	519953	7723
1817x	332x	519954	7724
1818x	334x	519955	7728
1819x	340x	519956	773x
1820x	360x	519957	7730
1828x	366x	519958	7736
18282	370x	519959	7737
1829x	372x	519963	7740
18292	373x	519965	7760
1832x	374x	520000	7766
1833x	376x	520001	7767
1835x	3761	520003	78xx
1836x	3774	520005	7800

Ostomy Care
Ceramide Infused
Products

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Ostomy Accessories

Continence
Care

Bowel Care

Tube
Fasteners

Wound Care

Index • Stock Number

Ostomy Care
Ceramide Infused
Products

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Ostomy Accessories

**Continence
Care**

Bowel Care

**Tube
Fasteners**

Wound Care

<u>Stock No</u>	<u>Stock No</u>	<u>Stock No</u>	<u>Stock No</u>
78500	84590	88700	9814
78501	846x	88800	9821
7905	8460	89xx	9822
7906	8469x	890x	9823
7910	84690	89001	9839
7917	847x	89002	9840
79300	8479xxx	89003	9873
79301	848x	89004	9880
79400	8489x	89005	Bxxx
79402	849x	8901	BxxxB
795xx	8499x	891x	FC0001
796xx	850xx	8931	FC0001-8
80xx	851xx	895xx	FC0002
80070	852xx	896x	FC0003
80110	853x	896xx	FC0006
8031	8531	898x	FC0006-8
8081	854xx	899x	FC0007
81xx	855x	90xxx	FC0008
81070	8551	91xxx	FC0011
81110	857x	910x	FC0012
8131	858xx	920x	FC0013
8181	859x	93xxx	FC0014
82xxx-30	861x	9342	FC0024
821xx	862xx	9343	
82100	863x	9345	
823xx	8631	9346	
82300	864x	9348	
82302	8641	9349	
824xx	8651	94xxx	
82400	866x	95xxx	
82402	867x	96xxx	
82500	872x	9645	
8261x	873x	9655	
8262xx	874x	97xxx	
827xx	875x	970x	
828x	877x	974xx	
829x	8770	975xx	
8331	88xx	976xx	
834x	8800	977x	
8371	8801	9779	
8381	881x	9780	
841x	883xx	9781	
8413xx	88300	9782	
843x	88302	9785	
8439xxx	884xx	9786	
844x	88400	9788	
8440	88401	9799	
8449x	88402	9800	
84490	88501	9805	
8450	887xx	9811	

Name

- A***adapt* Barrier Extenders
- Adapt CeraRing* Convex Barrier Rings
- Adapt CeraRing* Flat Barrier Rings
- Adapt CeraRing* Oval Convex Barrier Rings
- Adapt* Odor Eliminator and Lubricant
- Adapt* Ostomy Belts
- Adapt* Skin Barrier Paste
- Adapt* Skin Barrier Rings
- Adapt* Skin Barrier Strips
- Adapt* Skin Prepping Product
- Adapt* Stoma Powder
- Adapt* Adhesives and Removers
- Advance Plus Pocket* Touch Free Intermittent Catheter System
- Advance Plus* Touch Free Intermittent Catheter Kit
- Advance Plus* Touch Free Intermittent Catheter System
- AnchorFast Guard* Oral Endotracheal Tube Fastener
- AnchorFast* Oral Endotracheal Tube Fastener
- AnchorFast* Replacement Neck Band Straps
- Apogee* Family of Intermittent Catheters
- Apogee* HC Hydrophilic Intermittent Catheter
- Apogee* IC Intermittent Catheters, Firm Coudé Style
- Apogee* IC Intermittent Catheters, Firm Style
- Apogee* Intermittent Catheters, Other Styles
- Apogee Plus* System
- Apogee Plus* System Kit

Bedside Drainage Collection System with Anti-Reflux Valve

- C***alciCare* Calcium Alginate
- CalciCare* Calcium Alginate, Silver
- CenterPointLock* Closed Pouches
- CenterPointLock* Convex Skin Barriers with Floating Flange
- CenterPointLock* Drainable Pouches
- CenterPointLock* Flat Skin Barriers with Floating Flange
- CenterPointLock* Irrigator Sleeve
- CenterPointLock* Stoma Cap
- CenterPointLock* Urostomy Pouches
- CeraPlus* Convex Skin Barriers with Integrated Floating Flange
- CeraPlus* Flat Skin Barriers with Integrated Floating Flange

Drainable Pouch Clamp

- Drainable Pouches with *CeraPlus* Convex Barriers
- Drainable Pouches with *CeraPlus* Flat Barriers
- Drainable Pouches with *CeraPlus* Soft Convex Barriers
- Drain Tube Attachment Device (DTAD): Sterile

Ostomy Care
Ceramide Infused
Products

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Ostomy Accessories

Continence
Care

Bowel Care

Tube
Fasteners

Wound Care

Education

Index • Product Name

Name

E*veryday* Male External Catheter

Extended Wear Male External Catheter

Extension Tubing

Fabric Leg Bag Straps

FastCast OLS Total Contact Cast System Convenience Kit

Fecal Collectors

Feeding Tube Attachment Device (FTAD)

Female Urinary Pouch External Collection Device

Flexend (extended wear) Skin Barrier

Hollister Universal Catheter Access Port

Horizontal Tube Attachment Device (HTAD): Sterile

InView Male External Catheters

Irrigation Products

Irrigator Sleeve

Karaya 5 Closed Pouches

Karaya 5 Drainable Pouches

Karaya 5 Urostomy Pouches

Loop Ostomy Bridge

m9 Cleaner/Decrystallizer

m9 Odor Eliminators

Male Urinary Pouch External Collection Device

N*ew Image* Closed Pouches

New Image Convex Skin Barriers with Integrated Floating Flange

New Image Drainable Pouches

New Image Flat Skin Barriers with Integrated Floating Flange

New Image High Output Drainable Pouches

New Image Non Sterile Kits — Drainable

New Image Non Sterile Kits — Urostomy

New Image Sterile Kits — Drainable

New Image Sterile Kits — Urostomy

New Image Urostomy Pouches

O*nli* Ready to Use Hydrophilic Intermittent Catheter

P*ouchkins* Newborn Pouch

Pouchkins One-Piece Pouches

Pouchkins Pediatric Ostomy Belt

Pouchkins Premie Pouch

Pouchkins Two-Piece Pouches

Pouchkins Two-Piece Skin Barrier

Premier Closed Pouches with Flat Barriers

Premier Drainable Pouches with Convex Barriers

Premier Drainable Pouches with Flat Barriers

Premier Drainable Pouches with Soft Convex Barriers

Premier High Output Drainable Pouches

Name

Premier Non Sterile Kits — Drainable

Premier Non Sterile Kits — Urostomy

Premier Stoma Cap

Premier Urostomy Pouches with Convex Barriers

Premier Urostomy Pouches with Flat Barriers

Replacement Filters

Restore Calcium Alginate

Restore Calcium Alginate, Silver

Restore Cleanser & Moisturizer

Restore Contact Layer FLEX

Restore Contact Layer, Silver

Restore Extra Thin Hydrocolloid: Sterile

Restore Foam with Silicone Border

Restore Foam with Silicone with TRIACT Technology Non-Border

Restore Foam with TRIACT Technology Non-Border

Restore Hydrocolloid: Sterile

Restore Hydrocolloid with Foam Backing: Sterile

Restore Hydrogel

Restore LITE Foam with TRIACT Technology Non-Border

Restore Silver Foam with TRIACT Technology Non-Border

Restore Skin Cleanser

Restore Skin Conditioning Creme

Restore Wound Cleanser

Skin Barriers

Skin Barrier Pastes and Powders

Stoma Cap

Suction Tube Attachment Device (STAD)

TRIACT Foam with Silicone Border

Urinary Leg Bag

Urinary Leg Bag Combination Pack

Urostomy Drain Tube Adapter

Urostomy Pouches with CeraPlus Soft Convex Barriers

VaPro Hydrophilic Intermittent Catheter, Coudé Tip

VaPro Plus Pocket Touch Free Hydrophilic Intermittent Catheter

VaPro Plus Touch Free Hydrophilic Intermittent Catheter

Vented Urinary Leg Bag Combination Pack

Vertical Tube Attachment Device (VTAD): Sterile

Wound Drainage Collector

Ostomy Care
Ceramide Infused
Products

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Ostomy Accessories

Continence
Care

Bowel Care

Tube
Fasteners

Wound Care

Education

Index • HCPCS Codes

Education
Wound Care
Tube Fasteners
Bowel Care
Continence Care
Ostomy Care
 Ostomy Accessories
Ostomy Care
 Pouchkins Pediatric Pouching System
Ostomy Care
 Premier One-Piece Pouching System
Ostomy Care
 New Image Two-Piece Pouching System
Ostomy Care
 Ceramide Infused Products

HCPCS
 A4326
 A4328
 A4330
 A4331
 A4349
 A4351
 A4352
 A4353
 A4357
 A4358
 A4362
 A4363
 A4364
 A4367
 A4368
 A4371
 A4373
 A4385
 A4387
 A4388
 A4389
 A4390
 A4394
 A4397
 A4398
 A4399
 A4402
 A4405
 A4406
 A4407
 A4408
 A4409
 A4410
 A4411
 A4412
 A4414
 A4415
 A4416
 A4417
 A4419
 A4421
 A4424
 A4425

HCPCS
 A4428
 A4429
 A4430
 A4432
 A4435
 A4455
 A4456
 A5052
 A5054
 A5055
 A5056
 A5057
 A5061
 A5062
 A5063
 A5071
 A5114
 A5120
 A5122
 A5131
 A5200
 A6154
 A6196
 A6197
 A6199
 A6206
 A6207
 A6209
 A6210
 A6211
 A6212
 A6213
 A6234
 A6235
 A6236
 A6248
 A6250
 A6260
 A9270
 A9900
 A9999
 B9998

Ostomy Care
Ceramide Infused
Products

Ostomy Care
New Image Two-Piece
Pouching System

Ostomy Care
Premier One-Piece
Pouching System

Ostomy Care
Pouchkins Pediatric
Pouching System

Ostomy Care
Ostomy Accessories

Continence
Care

Bowel Care

Tube
Fasteners

Wound Care

Education

The Hollister logo, Adapt, Adapt CeraRing, Advance Plus, Advance Plus Pocket, AF300, AnchorFast, AnchorFast Guard, Apogee, Apogee Plus, Apogee IC, Apogee HG, CalciCare, CenterPointLock, CeraPlus, ComfortWear, Everyday, Extended Wear, FastCast OLS, Flextend, Flextend M, FlexWear, FormaFlex, Hollister, InView, Karaya 5, Lock 'n Roll, m9, New Image, Onli, Pouchkins, Premier, QuietWear, Restore, Secure Start, SoftFlex, TRIACT, VaPro, VaPro Plus and VaPro Plus Pocket are trademarks of Hollister Incorporated.

Endoform is a trademark of Aroa Biosurgery Limited.

Hydrofera Blue Ready is a trademark of Hydrofera LLC.

Remois is a technology of Alcare Co., Ltd.

All other trademarks and copyrights are the property of their respective owners.

Contact Us

Monday - Friday
8:00am - 5:00pm (CT)

For consumer or healthcare professional questions:

1.888.808.7456

For dealer, distributor, institution and
pharmacy questions or orders:

1.800.323.4060

FAX 1.847.680.1017



Hollister Incorporated
2000 Hollister Drive
Libertyville, Illinois 60048
1.800.323.4060

www.hollister.com

© 2018 Hollister Incorporated. Printed in USA.
923122-0918 US-00126